VOCO CediTEC® DT

CE 0482

EN Instructions for use MD EU Medical device

Product description:

CediTEC DT is a blank for the production of highly aesthetic permanent prosthetic teeth for removable dentures.

CediTEC DT contains 27% by weight inorganic fillers in a polymer matrix. Thanks to the composite technology, the material displays high abrasion resistance as well as high fracture resistance. The individual teeth or tooth sections are inserted into a prosthetic tooth luting agent.

CediTEC DT is available as a disc (dia. 98.4 x 20 mm).

Shades:

A1, A2, A3, BL (Bleach)

Indications:

Prosthetic teeth and tooth sections up to complete dental arches for removable dentures

Contraindication:

CediTEC DT contains methacrylates. **CediTEC DT** should therefore not be used on patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

CediTEC DT is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

CediTEC DT should only be used by a professionally trained dental practitioner.

Use:

Shade selection:

Select the shade using the VITA® shade system in daylight if possible.

Grinding/milling process:

Note: Designing the connector cross section height \geq width

CediTEC DT is suitable for both dry and wet processing. Select the corresponding blank size for the designed restoration and the grinding/milling parameters for **CediTEC DT**. When doing so, pay attention to the software settings of the respective CAD/CAM systems. "Diamond-coated tools" and composite parameters are recommended for the CAM processing. If the grinding/milling parameters are not already available in the CAD/CAM systems' software settings, they will need to be added before you proceed. Please contact the CAD/CAM system provider for assistance. The instructions for use from the manufacturers must be observed and followed. Check the finished restoration for defects such as cracks or chips. If the restoration displays any defects, it must be rejected.

Finishing:

Once the CAM process is complete, detach the restoration from the pin using a carbide tipped finishing bur or suitable cutting wheel. Sand the corresponding attachment points with fine-toothed carbide tipped finishing burs or fine diamond-coated grinders. Pay attention to any contact points.

Checking the accuracy of fit:

Check the prosthetic teeth's accuracy of fit before sticking them in position. Correct any occlusal interferences directly on the teeth.

Finishing/polishing:

The finishing/polishing can be performed in the same way as in analogue prosthetics.

Recommendation: **CediTEC DT** can be polished using a composite polishing paste in combination with goat's hair brushes and cotton/leather buffers.

Work using a low speed so as to avoid excessive abrasion.

The instructions for use from the manufacturers must be observed and followed.

Customisation:

For a highly aesthetic restoration, **CediTEC DT** can be customised, characterised or repaired at any time using a composite. Roughen the restoration surface via grinding or sandblasting (Al_2O_3 50-110 µm, 1-2 bar). Use a suction device to remove any dust. Remove abrasive material residues/dust carefully with a steam cleaner or ultrasound bath with water. Then dry the restoration with air. Apply a suitable adhesive system (e.g., **Futurabond U**) in accordance with the instructions for use. The instructions for use of the respective customisation systems must be observed.

Luting:

Preparation of adhesive surfaces:

For an optimal bond, blast the luting surface of the restoration with aluminium oxide (50 µm-110 µm, at 1 bar-2 bar). Use a suction device to prevent the accumulation of dust.

Remove abrasive material residues carefully with a steam cleaner or ultrasonic bath (70% ethanol). Then dry the restoration with air. Final cleaning with medical alcohol is possible. The instructions for use must be observed and followed.

CediTEC DT must be luted with a suitable luting agent system (e.g., **CediTEC Adhesive + CediTEC Primer**). The instructions for use must be observed and followed.

Note: – A transfer template or matrix can be used for the arrangement of the prosthetic teeth.
– **CediTEC DT** can be luted in a pressure pot at max. 55°C/2-6 bar.

Warnings/precautionary measures:

Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Storage instructions and application method:

Store at between 4°C and 23°C. If the material is stored in the refrigerator, bring it to room temperature before use. Do not store the product in a place where it is exposed to direct artificial light or sunlight. Do not use the product after its expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **CediTEC DT** must be reported to VOCO GmbH and the responsible authority.

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