VOCO Profibase

EN Instructions for use
MD EU Medical device

Product description:

Profibase is a light-curing acrylic available in the form of contoured plates for maxillae and mandibles. The material is ready for use and easy to handle. **Profibase** remains workable for a long period of time when exposed to artificial light and daylight, curing quickly under halogen and UV light. It retains its pleasant peppermint aroma even after curing.

Indications:

Fabrication of individual impression and functional trays for use in crown, bridge and total prosthetic technique to be used for the fabrication of bases of bite patterns as well as for dimensionally stable base plates for wax contouring in the total prosthetic technique.

Contraindications:

Profibase contains peppermint oil, methacrylates and BHT. **Profibase** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

Profibase is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

Profibase should only be used by a professionally trained dental practitioner.

Use:

1. Producing an individual tray

Allow **Profibase** to reach room temperature in good time before use. Mark the margins for the individual tray on the plaster model. If necessary, use spacers. If no spacers are required, the plaster model should be isolated. Remove the protective film and then place the plate on the model, mould it and trim the margins as indicated. The excess material can easily be used to shape the handle.

Light-curing:

Profibase is polymerised using a standard light unit (halogen or UV light). Depending on the thickness of the contoured tray, light-cure the reverse for 3 to 5 minutes. A longer polymerisation time will improve the physical properties of the material. As such, the curing time should be prolonged if necessary.

Finishing:

The tray produced in this way is ready for use immediately or can be further processed without delay. Oxygen in the air can cause a slight inhibition layer, which can be wiped away using a cloth moistened with alcohol.

2. Producing a base plate

Mark the margin for the planned denture base on the master model. Block out undercuts as usual. Isolate the model thinly with an alginate isolating material.

Carefully place a plate of **Profibase** onto the model, trim the margins with a penknife and polymerise the model with base plate for 5 minutes in a halogen or UVA light unit. One single exposure is sufficient.

Supplementary information:

Carbide burs are most suitable for processing the margins and the alveolar ridge. Introduce cross-shaped means of retention on the alveolar ridge before applying the wax. A thin layer of adhesive wax can then be applied.

Warnings, precautionary measures:

- Ensure sufficient and complete polymerisation! Exposure that is too short or not sufficiently intense will result in a loss in material strength.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Storage:

Store **Profibase** protected from light at temperatures between 4°C and 23°C. Do not use after the expiry date..

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Profibase** must be reported to VOCO GmbH and the responsible authority.

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