Instructions for use

**Product description:**

Fluoride varnish for the treatment of hypersensitivity, for deep fluoridation and caries prophylaxis.

**Presentation:**

Fluoride varnish made of synthetic resins for dental application.

**Composition:**

Sodium and calcium fluorides.

**Efficacy:**

**Bifluorid 12** is a fast-drying varnish which adheres well to dry enamel and dentine. This makes it possible to achieve long-term, intensive fluoridation with deep penetration of the treated surfaces, particularly at caries predilection sites.

**Bifluorid 12** forms a waterproof protective film insulating against thermal and chemical influences.

The objective of every application of **Bifluorid 12** must be to retain fluoride varnish as long as possible on the treated surfaces to allow the fluorides to work.

**Fields of application:**

For the treatment of hypersensitive necks of teeth and crown margins.

For fluoridation of dental enamel and desensitisation after grinding individual or several teeth.

For prophylaxis of secondary caries on cavity margins after fillings have been placed, especially when etching technique has been applied. Also after insertion of crowns, inlays and bridges.

**Method of use:**

Clean surfaces to be treated and dry thoroughly with air.

**Bifluorid 12** must be shaken well before each use. Solid particles must be well dispersed. The extra space provided in the bottle and two metal balls facilitate proper agitation.

Drop **Bifluorid 12** onto the brush, or preferably onto a **Pele Tim** foam pellet, and coat surfaces **thinly**.

Please note: Layers which are too thick will peel off too easily!

Dilute varnish if necessary.

**When using the SingleDose:**

Detach a *SingleDose* at the perforation and turn the printed side upwards. Pierce the foil with **Micro Tim**, enlarge the gap and moisten the brush with liquid through circular movement.

Allow the fluoride varnish to be absorbed for 10 - 20 s and then dry with air.

With correct application and proper mouth hygiene (**do not brush your teeth 12 - 24 h after application**) **Bifluorid 12** varnish will remain in place for several days. During this time the fluorides will act on the treated surface.

For hypersensitive dental necks, crown margins or similar indications, treatment should be repeated two or three times at intervals of 7 days.

Surfaces should be dried and coated again.

**Bifluorid 12** must be **sufficiently fluid** to run easily from the dropper after shaking. It may be necessary to dilute the varnish with a few drops of the solvent provided. This solvent should also be used for cleaning instruments etc.

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VOCO GmbH