VOCO Calcimol LC

EN Instructions for use MD EU Medical device

Product description:

Calcimol LC is a light-curing, radiopaque one-component material containing calcium hydroxide.

Indications:

- Indirect pulp capping
- As a cavity liner under all restorative materials

Contraindications:

Calcimol LC contains methacrylates, BHT and amines. **Calcimol LC** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

Calcimol LC is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

Calcimol LC should only be used by a professionally trained dental practitioner.

Use:

During application, ensure that the tooth surfaces are perfectly dry and avoid contamination. Instead of using **Calcimol LC**, use an aqueous calcium hydroxide product on dentine close to the pulp (< 0.5 mm) and for direct pulp capping. **Calcimol LC** can be applied directly from the syringe. To ensure functioning of the non-dripping syringe, the plunger must not be retracted during or after use. Dispense the required amount of **Calcimol LC** from the tube onto a mixing pad and apply using a suitable instrument. Cavity linings with a thickness in excess of 1 mm need to be built up in layers and each layer cured using blue light for at least 20 s.

Hold the output window of the light as close to the surface of the restoration as possible, as curing may otherwise be insufficient. Inadequate curing may lead to loss of strength and to discomfort.

Carefully remove any excess material prior to light-curing using suitable tools. Do not remove polymerised excess material with a probe but, instead, carefully use a rotating instrument and then place the restoration according to the manufacturer's instructions.

The attached application tip is intended for single use only. Remove it after use and seal the syringe again tightly with its original cap. The use of a protective barrier is recommended to protect the syringe against contamination with bodily fluids, dirt from hands or oral tissue.

If increased adhesion values are required (e.g. when using **Calcimol LC** for adhesive restorative treatments), a suitable adhesive can be used prior to the application of **Calcimol LC**. To this end, follow the adhesive manufacturer's instructions for use.

Warnings, precautionary measures:

- Phenolic substances, especially preparations containing eugenol or thymol, inhibit the curing of **Calcimol LC**. Therefore, do not use zinc oxide eugenol cements or other eugenol-containing substances in combination with **Calcimol LC**.
- Do not use **Calcimol LC** in the restoration margins. In the marginal area, the permanent restoration should always be placed on caries-free dental hard tissue which, when using the adhesive technique, must have been conditioned accordingly.
- **Calcimol LC** adheres to methacrylate-based products, e.g. crown and bridge materials, temporary sealing materials and certain impression materials (e.g. polyether), which can lead to removal of the cavity lining. In such cases, application of a separating agent is recommended (e.g. Vaseline, glycerine gel etc.). Remove the insulating material prior to placement of the final restoration.
- Once dispensed, **Calcimol LC** must not be exposed to direct light (operating light, sunlight) and must be applied as quickly as possible.
- Calcium hydroxide has a caustic effect on the skin and mucosal tissue, it must not be allowed to come into contact with the eyes!
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):

UDMA, Barium aluminium borosilicate glass, barium sulphate, fumed silica, aluminosilicate, TEGDMA, calcium hydroxide, DMAEMA, initiators, stabilisers, pigments

Storage instructions and application method:

Storage at 4°C-23°C protected from light and moisture. Store the syringe in the sealed aluminium pouch before first use. After opening the pouch, use the material within 6 months and before the shelf life has expired.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Calcimol LC** must be reported to VOCO GmbH and the responsible authority.

Note:

The Summary of Safety and Clinical Performance of **Calcimol LC** can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at www.voco.dental.

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