VOCO

Fissurit®/Fissurit® F

EN Instructions for use

MD EU Medical device

In accordance with DIN EN ISO 6874

Product description:

Fissurit is a transparent, light-curing fissure sealant. **Fissurit F** is a white, light-curing fissure sealant with fluoride.

Indications:

- Sealing/filling pits and fissures as well as small cavities
- Sealing/facing of damaged enamel surfaces
- Anchoring orthodontic appliances and covering caries predilection sites during orthodontic treatments
- Sealing composite or cement fillings (protective layer against humidity)
- Sealing carious lesions to stop them
- Restoring small defects in composite and amalgam fillings
- Sealing deciduous teeth

Contraindications:

Fissurit/Fissurit F contains methacrylates and BHT. **Fissurit/Fissurit F** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

Fissurit/Fissurit F is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

Fissurit/Fissurit F should only be used by a professionally trained dental practitioner.

Use:

Fissures should be sealed as early as possible, as soon as the eruption of the tooth allows treatment and moisture control is possible.

Preparation:

After isolating the operating site, thoroughly clean the surfaces to be sealed (rotating brush/cleaning paste), spray and dry. In case of incipient carious lesions or unclear diagnosis, extend the fissure with a fine instrument, if necessary.

Etching:

For example with **Vococid** gel or liquid. Please refer to the **Vococid** instructions for use.

Application from the bottle:

Apply and distribute **Fissurit/Fissurit F** carefully with a brush or other suitable instrument, remove any air bubbles (probe).

Syringe application:

Remove the protective cap from the syringe. Screw a new cannula onto the **Fissurit F** syringe. To avoid air bubbles, hold the syringe vertically and press the material up to the tip of the cannula.

Fissurit F flows out of the cannula evenly at light to medium pressure.

The plunger must not be retracted during or after use, otherwise the syringe may drip or run.

Apply **Fissurit/Fissurit F** into the fissure or cavity from the margin and then allow it to act for 15-20 s. Remove any air bubbles (e.g. using a fine probe).

The attached application tip is intended for single use only. Remove it after use and seal the syringe again tightly with its original cap. The use of a protective barrier is recommended to protect the syringe against contamination with bodily fluids, dirt from hands or oral tissue.

Light-curing:

Polymerise **Fissurit/Fissurit F** (20 s per surface) and check the sealing. Conventional polymerisation devices are suitable for light-curing the material. The light output should be at least 500 mW/cm² on halogen polymerisation devices and 300 mW/cm² on LED devices.

Finishing:

Re-seal ditchings, remove excess cautiously (e.g. finishing diamond), polish if necessary (rubber polisher). In the final step, apply a fluoride varnish to all the treated teeth or the entire set of teeth if possible.

Warnings, precautionary measures:

- Inadequate isolation from moisture is associated with the risk of a (partial) premature loss of the sealing. In such cases, the application of a fluoride varnish or sealing with a glass ionomer cement is recommended until perfect sealing is possible.

- If the etched surface is contaminated (e.g. by saliva), the adhesion may be reduced. In case of contamination, clean, etch again for approx. 10-20 s, spray and dry.

- Adhesion can be increased by using a bonding agent after etching (e.g. in case of difficult moisture control / extended sealing). Only expose **Fissurit/Fissurit F** to the operating light for a short time, as otherwise polymerisation will start (viscosity increases).

- Check the sulcus after application and remove any residues of material. Avoid contact with the mucosa as much as possible.

- Do not use cleaning or drying agents containing phenolic substances (e.g. eugenol, thymol) as they inhibit complete polymerisation.

- Carious lesions (e.g. around fillings) can be stopped by appropriate sealing if renewal of the filling is not desired or not possible. For this purpose, clean the areas affected, prepare as minimally as possible and etch in the enamel area. In such cases, sealing of the complete filling is recommended.

- In the case of orthodontic treatments in areas of reduced oral hygiene, sealing of these areas can provide good and permanent protection.

- Fluoridation should be continued in all cases to protect the unsealed tooth areas from caries.

- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):

Fissurit: BisGMA, UDMA, HEDMA, fumed silica, TEGDMA

Fissurit F: BisGMA, UDMA, HEDMA, fumed silica, sodium fluoride, TEGDMA, titanium dioxide

Storage instructions and application method:

Close **Fissurit/Fissurit F** containers again carefully immediately after use. Store **Fissurit/Fissurit F** at 4°C to 23°C and away from light. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Fissurit/Fissurit F** must be reported to VOCO GmbH and the responsible authority.

Note:

The Summary of Safety and Clinical Performance of **Fissurit/Fissurit F** can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at www.voco.dental.

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VC 60 AA1080 E1 0523 99 © by VOCO