### EN Instructions for use

**MD** EU Medical device

#### Product description:

**Grandio disc** is a nano-hybrid composite CAD/CAM disc for the fabrication of permanent single-tooth restorations.

**Grandio disc** contains 86 % w/w inorganic fillers in a polymer matrix. The restorations produced must be placed with an adhesive composite luting agent in combination with a dentine-enamel bonding agent.

**Grandio disc** is available in one size (ø 98,4 × 15 mm).

#### Shades (two degrees of translucency):

High Translucent (HT): A1, A2, A3, A3.5  
 Low Translucent (LT): A1, A2, A3, A3.5, B1, C2, BL  
 Multicolor: A1, A2, A3, A3.5, B1, C2

#### Indications:

- Crowns, inlays, onlays, veneers
- Implant-supported crowns

#### Contraindications:

**Grandio disc** contains methacrylates. **Grandio disc** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

#### Patient target group:

**Grandio disc** is suitable for use in all patients without any age or gender restrictions.

#### Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

#### User:

**Grandio disc** should only be used by a professionally trained dental practitioner.

#### Use:

##### Shade selection:

Use the VITA® shade system (**GrandioSO**) to determine the shade against the cleaned but unprepared tooth prior to anaesthesia and preferably in daylight conditions.

##### Core and cavity preparation:

As a general rule, the core/cavity preparation should be performed in accordance with the conventional rules for all-ceramic restorations. That means rounding off interior corners and edges and using a shoulder preparation with rounded interior edges or a chamfer preparation. The temporary restoration and any contamination should be mechanically removed (e.g. using pumice and a polishing brush) and the area subsequently rinsed with water.

Observe the minimum thicknesses stated below.

**Note:** Do not create any undercuts.

##### Crowns, inlays and onlays (see Fig. 1 - 3):

- Preparation angle of 4° - 6°
- Minimum occlusal wall thickness 1.5 mm
- Minimum buccal wall thickness 0.8 mm
- Minimum wall thickness below load-bearing cusp 1.5 mm
- Minimum cervical wall thickness 0.8 mm

##### Veneers (see Fig. 4):

- Minimum cervical wall thickness 0.4 mm
- Minimum labial wall thickness 0.5 mm
- Minimum incisal wall thickness 0.5 mm

##### Cement space:

70 µm (± 10 µm)

##### Grinding/milling process:

Select the correct disc size for the planned restoration as well as the milling parameters for **Grandio disc** or for composites. When doing so, pay attention to the software settings of the respective CAD/CAM systems.

**Processing requires diamond instruments.** First of all, the settings for disc sizes and the milling parameters need to be created, if they are not already included in the software settings of the CAD/CAM systems. Please contact the CAD/CAM system provider for assistance. The operating instructions from the manufacturers must be observed and followed.

Please note the software settings of the respective CAD/CAM systems and the direction of the multicolour restoration.

Check the processed restoration for defects such as fissures or material chipping. If the restoration displays any defects, it must be rejected.

**Note: Grandio disc is an already completely cured composite and must not be fired.**

##### Polishing:

The restoration can be finished and polished with standard composite polishers (e.g. medium, fine or extra fine diamond grinding tools, polishers), this can be done either intraorally or extraorally.

##### Adhesive luting:

###### Preparing the restoration:

For an optimal bond, blast the luting surface of the restoration with aluminium oxide (25 µm - 50 µm) at 1.5 bar - 2 bar. Use a suction device to remove the dust produced.

**Note:** Do not etch with hydrofluoric acid or phosphoric acid.

Clean the pretreated restoration using an ultrasonic bath or steam cleaner and then dry with oil-free air. Final cleaning with medical alcohol is possible. Treat the luting surface with a silane coupling agent such as **Ceramic Bond**. Apply **Ceramic Bond** to the inside of the work and allow it to work for 60 seconds. Then dry for 5 seconds with oil-free air. The instructions for use must be observed and followed.

###### Adhesive luting of the restoration:

The processed restoration must be luted with a composite-based, adhesive luting agent. The instructions for use must be observed and followed. **VOCO** recommends the dual-curing luting system **Bifix QM** in combination with **Ceramic Bond** for an optimal bond to the restoration and a suitable dentine-enamel bond (**Futurabond DC/Futurabond U**). Alternatively, it is also possible to use a purely light-cured, flowable composite/ORMOCER® (e.g., **GrandioSO Flow**) if the restoration has thin walls (≤ 2 mm).

##### Modifications:

**Grandio disc** is a nano-hybrid composite CAD/CAM restorative material which can be modified, characterised or repaired at any time using a composite/ORMOCER®. Roughen the surface of the restoration (for example with a diamond bur) and apply a suitable adhesive system (e.g., **Futurabond U**) in accordance with the instructions for use. Then apply the composite/ORMOCER® (e.g. **GrandioSO**, **GrandioSO Flow**, **FinalTouch**) in accordance with the respective instructions for use.

##### Warnings, precautionary measures:

– Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

#### Constituents (in descending order according to content):

Barium aluminium borosilicate glass, silicon dioxide, polymethacrylate, stabilisers, pigments

#### Storage:

Storage at 4°C-23°C. Do not store the product in a place where it is exposed to direct artificial light or sunlight. Do not use after the expiry date.

#### Disposal:

Dispose of the product in accordance with local regulations.

#### Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Grandio disc** must be reported to VOCO GmbH and the responsible authority.

#### Note:

The Summary of Safety and Clinical Performance of **Grandio disc** can be found in the European database on medical devices (EUDAMED) - <https://ec.europa.eu/tools/eudamed>. Detailed information can also be found at [www.voco.dental](http://www.voco.dental).

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VC 60 AA6051 E1 1122 99 © by VOCO
ORMOCER® is a registered trademark of Fraunhofer Gesellschaft zur Förderung der angewandten Forschung e.V., München, Germany.
VITA® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co.