Application capsule

Product description:
Ionolux AC is a light-curing, radiopaque glass ionomer restorative material in user-friendly capsules. The material has outstanding aesthetics, good handling properties and high chemical adhesion.
Ionolux AC is available in shades A1, A2, A3, A3.5 and B1. The usable quantity of cement per capsule is at least 0.17 ml.

Indications:

- Class III and V restorations, primarily cervical fillings and root caries
- Restoration of deciduous teeth
- Small class I fillings
- Temporary fillings
- Core build-ups
- Linings

Contraindications:

Ionolux AC contains polyacrylic acid, BHT and methacrylates. Ionolux AC should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

Ionolux AC is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

Ionolux AC should only be used by a professionally trained dental practitioner.

Use:

Preparation:

In principle, cavity preparation should be minimally invasive to conserve healthy tooth substance. There is no need for undercuts.

Class V cavities should have a minimum depth of 1 mm. Remove preparation residues with water spray, then dry with air so that the surface remains slightly glossy (do not overdry). Place a matrix, if necessary.

If the residual dentine thickness is insufficient, the area close to the pulp should be covered with a calcium hydroxide product.

Activating the application capsule and mixing the material:

Only extract the application capsule from the tubular bag immediately before use. During activation, the liquid from the plunger must be pressed into the inside of the capsule as fully as possible. To achieve this, the end of application capsule where the plunger of liquid protrudes is slowly pressed down onto a hard surface (e.g. a table) as far as the plunger will move (see fig. 1). The capsule must only be activated directly before use and loaded immediately into a high-frequency mixer with approx. 4,000 oscillations per minute. Mix the activated application capsule for 10 s, remove it from the mixer without delay and insert it in the holder of an application device (e.g. AC Applicator, see corresponding instructions for use).

Application:

Carefully open the outlet by lifting the cannula until you feel a stop (see fig. 2), and apply the material directly into the cavity within 10 s of the end of mixing by compressing the lever on the applicator slowly and evenly. The cannula can be aligned as required by turning the whole application capsule inside the holder. Working time and consistency may be influenced by the respective type of capsule mixer and mixing time.

Fig. 1: Activating the capsule

Fig. 2: Opening the capsule

Avoid unnecessary waiting times between application and finishing since the setting process of the material may make it more difficult to work with. Apply the selected shade of Ionolux AC in layers with a maximum thickness of 2 mm, adapt with a suitable instrument and light-cure afterwards. At room temperature, Ionolux AC has a working time of 2.5 min. Higher temperatures shorten the working time, lower temperatures prolong it.

Light-curing:

Conventional polymerisation devices are suitable for light-curing the material. The curing time when using an LED-/halogenlight with a minimum light-output of 500 mW/cm² is 20 s per layer.

Hold the light emission tip of the device as close as possible to the surface of the filling. If the distance is more than 5 mm, the curing depth may be reduced.

Quantity per capsule [ml]	min. 0.17
Net weight per capsule [g]	min. 0.35
Mixing time [s]	10
Working time [min]	min 2.5
Light-curing [s]	20

Finishing:

After light-curing, remove excess material and finish and polish the filling as usual.

Warnings, precautionary measures:

- Ensure that the application capsule is mixed immediately after activation.
- In case of contact with oral mucosa or skin, remove immediately with a cotton swab soaked in alcohol. Rinse with water.

- In case of eye contact, rinse well with plenty of water and consult an ophthalmologist.
- Phenolic substances, especially preparations containing eugenol or thymol, inhibit the curing of Ionolux AC. Therefore, do not use zinc oxide eugenol cements or other eugenol-containing substances in combination with Ionolux AC.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):
Fluorosilicate glass, polyacrylic acid, HEMA, water, UDMA, GlyDMA, fumed silica, initiators, stabilisers, pigments

Storage instructions and application method:
Store Ionolux AC in a dry place at temperatures of 4 °C - 23 °C. Only remove it from the foil immediately before use. Do not use after the expiry date.

Disposal:
Dispose of the product in accordance with local regulations.

Reporting obligation:
Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of Ionolux AC must be reported to VOCO GmbH and the responsible authority.

Note:
The Summary of Safety and Clinical Performance of Ionolux AC can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at www.voco.dental.

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VC 60 BB2115 E1 0822 99 © by VOCO