VOCO**Rebilda® Post GT**

Clustered glass fibre-reinforced composite post

Instructions for use

Product description:

The radiopaque, translucent Rebilda Post GT is made of glass fibre-reinforced composite and exhibits elasticity behaviour similar to that of dentine. The adhesive technique makes it possible to obtain a metal-free restoration which is durable and of a high aesthetic standard. Rebilda Post GT consists of several fine individual posts (ø 0.3 mm) in different quantities, held together by a sleeve. Rebilda Post GT is available in four sizes.

Size	blue	red	green	black
Number of individual posts	4	6	9	12
Length (without sleeve)	20 mm	20 mm	20 mm	20 mm
idealised (round) diameter	0.8 mm	1.0 mm	1.2 mm	1.4 mm

Indication:

Core build-up restorations with root posts for endodontically treated teeth.

Application instructions:

A) Pre-treatment of root canal and Rebilda Post GT

- Successful completion of endodontic treatment
- Utilisation of a rubber dam is recommended
- Selection of the appropriate post size based on X-ray
- Correct determination of the required working length (apically, a minimum of 5 mm endodontic filling must be retained, ½ of the post should be located in the root canal and ½ in the core build-up)
- Removal of the root canal filling material using, for example, a Gates-Glidden bur, a Peeso reamer and/or other hot instrument down to the proposed depth. To ensure that the adhesive luting process is not disrupted, there must be no residues of gutta-percha and/or sealers left on the walls of the root canal.
- If necessary, widen the root canal using a drill (1000 - 2000 min⁻¹ instrument speed; green angle piece). Avoid preparation beyond the apex and/or lateral perforation of the root.
- To ensure optimal drilling performance, intermittently rinse the root canal and the drill to remove dentine residues
- Clean Rebilda Post GT with alcohol before the try-in
- Check the position in the mouth
- Check by X-ray
- Clean Rebilda Post GT with alcohol and dry
- Carry out silanisation of post (e.g. Ceramic Bond – allow to act for 60 s and then dry with oil-free air)
- Finally, rinse the root canal out with alcohol (> 75 %). [Do not use hydrogen peroxide (H_2O_2), other disinfectants or EDTA solutions, as any residues may compromise the adhesive strength and the curing reaction of the adhesive luting material]
- Then dry with paper points

B) Adhesive luting with Rebilda DC (in combination with a dual-curing adhesive)

1. Bonding

For luting Rebilda Post GT, use a dual-cure adhesive (e.g. Futurabond U) in accordance with the respective instructions for use.

2. Using the cartridge/QuickMix syringe

It is recommended to apply Rebilda DC directly into the root canal, either from the cartridge (using the mixing tip type 19, which can be pre-bent by hand if necessary) or from the QuickMix syringe, using the intraoral tip type 1. Always apply the material starting from the apical end, and leave the tip immersed in the luting composite during application. Do not use a lentulo! Additionally, moisten Rebilda Post GT with Rebilda DC. Insert Rebilda Post GT ensuring to leave small amounts of excess cement. Remove any excess cement whenever not using Rebilda DC for core build-up.

3. Removing the sleeve/shortening the posts

a) Shortening the posts before polymerisation

Shorten Rebilda Post GT below the sleeve with a fine diamond bur at high speed, using suction but no water cooling. Remember at this point that the core build-up material will later need to cover Rebilda Post GT completely. Adjust the individual posts of Rebilda Post GT to fit the geometry of the root canal, using a suitable instrument (e.g. a spreader).

To fix the post and cure the dual-cure adhesive applied to the occlusal surface, light-cure for 40 s.

b) Shortening the posts after polymerisation

Hold Rebilda Post GT in place using tweezers and remove the sleeve with a second pair of tweezers. Adjust the individual posts of Rebilda Post GT to fit the geometry of the root canal, using a suitable instrument (e.g. a spreader). To fix the post and cure the dual-cure adhesive applied to the occlusal surface, light-cure for 40 s. Shorten Rebilda Post GT with a fine diamond bur at high speed, using suction but no water cooling. Remember at this point that the core build-up material will later need to cover Rebilda Post GT completely.

If using a luting composite from another manufacturer, proceed according to the respective information provided by the manufacturer.

C) Core build-up

- Follow instructions for use for the selected core build-up material (e.g. Rebilda DC).

Cleaning, disinfection and sterilisation (Endodontic posts are not sterile when delivered):

- Rebilda Post GT can be disinfected with medical alcohol

Information, precautionary measures:

- Rebilda Post GT contains methacrylates and BHT. Rebilda Post GT should not be used in case of known hypersensitivities (allergies) to any of these ingredients.
- Do not use Rebilda Post GT in cases of less than 2 - 3 mm of supragingival dental tissue.
- Rebilda Post GT must not be used for bruxism patients or those with a deep overbite (cover bite).
- The posts are intended for single use only and must not be reused.
- Avoid using materials containing eugenol, which also includes some materials for root canal fillings. They may interfere with the curing of luting composites.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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