voco

Twinky Star® Flow

EN Instructions for use
MD EU Medical device
 In accordance with DIN EN ISO 4049

Product description:
Twinky Star Flow is a flowable, coloured compomer with glitter effect. This radiopaque restorative material for cavities in deciduous teeth contains anorganic fillers at a level of 65 % by weight inside a methacrylate matrix and is cured with halogen light /LED light (blue light).

Twinky Star Flow is available in the non-dripping syringe and is used in combination with a dentine enamel bond.

Shades:

Blue, Pink

Indications:

- Fillings in deciduous teeth

Contraindications:

Twinky Star Flow contains methacrylate and BHT. **Twinky Star Flow** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

Twinky Star Flow is suitable for use in all patients with first dentition of any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

Twinky Star Flow should only be used by a professionally trained dental practitioner.

Use:

Preparation:

Clean the teeth being treated. If applicable, mark occlusal contact points. Bring material to room temperature before use.

Shade selection:

Choose one of the 2 colours indicated in the shade guide

Cavity preparation:

In principle, the cavities should be prepared according to the rules for adhesive restorative treatment; preparation should be minimally invasive in order to protect healthy dental hard tissue.

Follow special preparation forms for deciduous teeth!

Cleaning:

- Remove any residues in the cavity with a water jet.
- Avoid contamination with blood or saliva after cleaning.
- Remove any excess water with a faint airjet. Do not overdry the dentine! The dentine surface should be moist but not too wet (wet bonding).
- Secure a dry working field.

Pulp protection:

It is recommended to apply a calcium hydroxide lining in areas close to the pulp.

Matrix placement:

Take appropriate measures for areas where the applied **Twinky Star Flow** could flow out. Translucent matrices that are wedged in the approximal area are advantageous. Minimal separation facilitates the shaping of the approximal contact and insertion of the matrix.

Bonding materials:

Twinky Star Flow is applied with a dentine enamel bond, using the adhesive technique. Any light-curing bonding materials can be used. Observe the relevant instructions for use during preparation (etching technique) and processing.

Application:

Lock the enclosed application cannula on the syringe with a clockwise turn and directly apply **Twinky Star Flow**. In order to ensure the functionality of the non-dripping syringe do not pull back the plunger during or after use.

For fillings with a thickness of over **2 mm**, apply and cure in layers.

The attached application tip is intended for single use only. Remove it after use and seal the syringe again tightly with its original cap. The use of a protective barrier is recommended to protect the syringe against contamination with bodily fluids, dirt from hands or oral tissue.

Light-curing:

Conventional polymerisation devices are suitable for light-curing the material. Depending on light output, polymerisation times are as follows:

LED-/halogen lamp	all colours
$\geq 1000 \text{ mW/cm}^2$	20 s
$500 - 1000 \text{ mW/cm}^2$	40 s

Bring the light aperture as close as possible to the light-curing unit as otherwise the curing will be less effective. Inadequate curing can cause discolouration and complaints.

Finishing:

After removing the matrices, the filling can be trimmed and polished immediately while being cooled (e.g. fine or extra-fine diamond burs, polisher). The tooth should be fluorinated as a final step.

Warnings, precautionary measures:

- Phenolic substances, especially preparations containing eugenol or thymol, inhibit the curing of **Twinky Star Flow**. Therefore, do not use zinc oxide eugenol cements or other eugenol-containing substances in combination with **Twinky Star Flow**.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):

Barium aluminium borosilicate glass, strontium aluminium borosilicate glass, TEGDMA, BisGMA, UDMA, fumed silica, fluorosilicate glass, initiators, stabilisers, pigments, glimmer

Storage instructions and application method:

Storage at 4°C-23°C. Seal the syringes immediately after removing the material in order to avoid exposure to light and subsequent polymerisation.

Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Twinky Star Flow** must be reported to VOCO GmbH and the responsible authority.

Note:

The Summary of Safety and Clinical Performance of **Twinky Star Flow** can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at www.voco.dental.

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