# VOCO

## Fit Test® C&B

**EN** Instructions for use  
**MD** EU Medical device

### Product description:

**Fit Test C & B** is an addition-curing silicone for testing the accuracy of fit of crowns, bridges, inlays/onlays, secondary components and attachments as well as temporary blocking during the permanent cementing of dentures. For dentists, **Fit Test C & B** represents a diagnostic aid for optimal adaptation of the final restoration. The blue staining makes it possible to identify problem areas on both metal and all-ceramic restorations reliably.

**Fit Test C & B** impresses with simple handling and particularly high tear resistance.

### Indications:

Control of precise fit, marginal tightness and marginal gap of

- full or partial crowns
- bridges
- inlays/onlays
- cast post abutments/post crowns
- metal bases
- implant-supported crowns and bridges
- attachments and other secondary parts of dentures

### Contraindications:

- The use of **Fit Test C & B** is also contraindicated when non-hardening, silicone-based base materials (e.g., rebasings) are present, as a chemical bond with such materials is possible.
- The following materials can prevent or delay the setting reaction of **Fit Test C & B**: condensation-curing silicone, polysulphide impression materials, materials containing eugenol, incompletely cured acrylics and latex gloves. These materials must not be used when working with **Fit Test C & B**.
- **Fit Test C & B** contains platinum catalysts. **Fit Test C & B** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

### Patient target group:

**Fit Test C & B** is suitable for use in all patients without any age or gender restrictions.

### Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

### User:

**Fit Test C & B** should only be used by a professionally trained dental practitioner.

### Use:

#### Preparation:

Clean and dry the restoration, preparation and abutment as usual. For technical reasons, there may be slight differences in the fill levels of both chambers prior to initial use. You should therefore remove the cap and continue to press out the material until an equal amount has been ejected from both outlets. Apply the mixing tip correspondingly (observe the slot!) and fix in position by turning 90° clockwise. The material is mixed as it is ejected from the tip. All QM syringes come with a detachable piston. The piston can easily be pulled backwards out of the syringe and then reinserted. The seals in the syringe do not change their position whilst this is being done. This prevents cross-contamination in the syringe. The warranted product characteristics can only be achieved when using the supplied original mixing tips, which are also available as accessories. The attached mixing tip is intended for single use only. After use, store the QuickMix syringe tightly sealed. Check that material flows freely through the openings before reusing.

#### Application:

**Fit Test C & B** can be applied into the cleaned restoration or cavity either directly from the mixing tip or, in hard-to-reach areas, using a suitable instrument. When doing so, ensure that all surfaces are covered with **Fit Test C & B** and no air bubbles are trapped. Position and secure the restoration (do not let the patient bite down!).

#### Curing:

The working time from the beginning of the mixing is at least 1 min; the intraoral setting time is at most 2 min – use excess material to monitor the setting. Once set, carefully remove the restoration with the silicone film adhering. Small, single-surface inlays can be removed with adhesive wax if necessary.

#### Assessment:

Even layers of silicone indicate a good accuracy of fit. Uneven layers, in contrast, require correction of the restoration/preparation at the identified points using standard techniques.

#### Warnings, precautionary measures:

- Avoid getting uncured silicone on objects such as clothing as it is difficult to remove. The working times specified are based on a room temperature of approx. 23°C and relative humidity of approx. 50%. Lower temperatures extend and higher temperatures shorten the times required.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

#### Storage:

Store at 4°C-23°C. Do not use after the expiry date.

#### Disposal:

Dispose of the product according to local regulations.

#### Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Fit Test C & B** must be reported to VOCO GmbH and the responsible authority.

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VC 60 BB2095 E1 0422 99 © by VOCO