# VOCO V-Print® splint comfort

**EN** Instructions for use  
**MD** EU Medical device

## Product description:

**V-Print splint comfort** is a light-curing resin for the generative manufacturing of thermoflexible splints (including therapeutic splints) as well as thermoflexible auxiliary parts and functional parts for dental diagnostics using CAD/CAM technology.

## Indications:

- Therapeutic splints
- Auxiliary parts and functional parts for diagnostics
- Bleaching splints (home bleaching)
- Palatal plates

## Contraindications:

**V-Print splint comfort** contains (meth)acrylates and phosphine oxide. **V-Print splint comfort** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

- Palatal plate with velar extension (spur)

## Patient target group:

**V-Print splint comfort** is suitable for use in all patients without any age or gender restrictions.

## Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

## User:

**V-Print splint comfort** should only be used by a professionally trained dental practitioner.

## Hardware and software requirements:

|   |   |
|---|---|
| CAD software <sup>1</sup><br>dental scanner | Software for the planning and design of therapeutic splints as well as splints, auxiliary parts and functional parts for use in dental diagnostics.<br>The software and dental scanner must satisfy local and current medical device specifications and allow for issuance of the patient-specific design as an STL data set. |
| CAM software                                | Software for preparation of the print order. The part will not be modified during this process. Structures that facilitate the 3D printing are simply created.<br>For example:<br>– Autodesk Netfabb version 2020 or later for SolFlex 3D printing.   |

<sup>1</sup>The designation **Software as Medical Device SaMD** includes standalone (autonomous) software that is a medical device (MD) and not part of one.

|                         |   |
|-------------------------|---|
| Manufacturing equipment | For example:<br>VOCO SolFlex 170<br>VOCO SolFlex 350<br>VOCO SolFlex 650<br>VOCO SolFlex 170 HD |
| Post-curing devices     | For example:<br>Otoflash G171   |

See also: accompanying list of resources or [www.voco.dental/3dprintingpartners](http://www.voco.dental/3dprintingpartners)

All manuals and/or operating instructions for the respective programmes, and for device, materials and/or parts manufacturers, which are required for the manufacturing process, must be observed.

Clarify ahead of time whether the programmes, devices and/or objects that you intend to use have been designed and approved for the corresponding applications.

**CAUTION:** Non-authorised changes to the process equipment, parameters, or software could result in the **V-Print splint comfort** end object not satisfying specifications.

## Use:

### Preparation:

For an indication-appropriate CAD construction, the following design conditions must be observed:

- Minimum wall thickness 1 mm / Minimum wall thickness in case of an occlusal construction of adjusted splints: 2-3 mm
- Recommended gap: 70 µm
- Round off internal edges on the outside.

Prepare a print job using CAM-software.

### Processing:

**Note:** Use separate material containers and cleaning baths for each printing material, in order to prevent cross contamination.

The materials container should be filled immediately before the start of the printing process. It is important to ensure that the material is free of bubbles to the extent possible, and filled to the fill level mark.

Start the print job observing the parameters that you previously selected.

Once the printing process has ended, a dripping time of approximately 10 minutes is recommended. Next, carefully detach the printed objects from the build platform.

In the following steps, the printed objects will need to be cleaned, dried and post-exposed, in order to guarantee the required product characteristics. A detailed explanation of the steps outlined above can be found under **Post-processing**.

**Recommendation:** After completing your work, transfer the remaining material from the material tray to the original container (use a stainless steel sieve if necessary). This serves to check the material tray and also enables optimum storage of the print material.

### Post-processing:

#### Cleaning

For cleaning purposes, use isopropanol (purity ≥ 98%) as a cleaning solution, in a cleaning device. An unheated ultrasonic bath or an unheated stirring bath may be used as a cleaning device.

The printed objects must be cleaned in two, or optionally in three steps. Position the unclean printed objects inside the cleaning bath so that any openings point downwards. Use tweezers or appropriate lowering baskets for transferring the printed objects into the detergent solution.

Please ensure that the printed objects do not come into contact with one another during cleaning.

|                         |  |   |
|-------------------------|--|---|
|                         | Ultrasonic bath  | Stirring bath                             |
| Pre-cleaning (optional) | Carefully pre-clean the printed objects by submerging them several times in a beaker with isopropanol. |   |
| Preliminary cleaning*   | 5 minutes<br>– may be used multiple times  | 5 minutes<br>– may be used multiple times |
| Final cleaning          | 3 minutes<br>– fresh cleaning bath   | 3 minutes<br>– fresh cleaning bath        |

\*Note: The bath's cleaning efficacy decreases with increased use. Resin residue on the surface could indicate that the cleaning efficacy of the bath has diminished, or that the parts came into contact with one another. When the cleaning efficacy decreases, the respective bath must be replaced.

Next, the printed objects must be dried carefully using compressed air. If there is any resin residue on the printed object after the final cleaning, or if residue escapes from the undercut when drying, the printed object can be briefly immersed once again in the final cleaning bath. Next, repeat the drying process.

#### Preparation for post-exposure:

Obstructive support structures can be removed before the post-exposure process by using a rotary instrument as close to the printed object as possible, carefully and without exerting pressure. Use a suction device. Remove any remaining resin dust carefully using compressed air. Then rinse printed objects for a few seconds with fresh isopropanol and remove any dust using a fine, clean brush if necessary. Then dry the printed objects again carefully using compressed air.

#### Post-exposure:

Conduct the post-exposure a minimum of 15 minutes after the most recent contact with isopropanol. A protective gas atmosphere is not mandatory, but is recommended for easier post-processing.

It is important to ensure that the printed objects do not overlap or contact each other, as post-exposure would be negatively affected by the shadows that are cast.

Post-exposure can be conducted using the following devices:

| Post-exposure device   | Programme        |  |
|--|------------------|--|
| For example:<br>Xenon photoflash unit<br>Otoflash G171<br>(optional with N2) | 2 x 2000 flashes | After 2000 flashes, observe a cooling phase of at least 2 minutes with open lid. Next, turn over and light-cure with another 2000 flashes. |

See also: accompanying list of resources

#### Finishing:

Generally work with low contact pressure and reduced speed. This guarantees consistent results and reduces processing marks.

In order to sand the support stubs use a rubber polisher or a fine-toothed carbide bur, for example. The bur can also be used for additional finishing of special structures.

In order to achieve form-fit grinding, e.g. between the support base and the print object, it is recommended to process the surface in the corresponding area with coarser or finer rubber polishers. A corresponding result can also be achieved with sandpaper, possibly with different grit sizes.

In order to obtain a high-gloss finish, polish the object first with a pumice stone. Then, thoroughly remove any pumice dust under running water using a brush. Finally, polish the object to a high-gloss using a buffer and high-gloss polishing paste, without applying excessive pressure.

#### Final cleaning:

Clean the object thoroughly. Coarse residue can be removed with the steam jet cleaner. The arising deformation is reversible. The final cleaning can be conducted by briefly placing the object in an unheated ultrasonic water bath. In order to remove oily or fatty contaminants, a surfactant solution may be used in place of water.

#### Disinfection:

Objects manufactured from **V-Print splint comfort** may be disinfected using alcohol- or aldehyde-based disinfectants (e.g. ethanol (≥ 70 %), MD 520 by Dürr, Cavex ImpreSafe by Cavex). Observe the manufacturer's Instructions for use.

#### Warnings, precautionary measures:

- Only use **V-Print splint comfort** intraorally in a fully cured state. Pay attention to the finishing process.
- Contact between uncured **V-Print splint comfort** and the skin/mucous membranes and eyes can cause irritation and should be avoided.
- The wearing of protective clothing is recommended. Furthermore, it is important to ensure that no vapours and/or dusts are inhaled. The wearing of a suitable mask and/or the use of suction devices is recommended.
- For a more comfortable intraoral insertion or to check the fitting on a model it is recommended to heat the splint in warm water (40°C-60°C).
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for the intended purpose.

#### Storage:

Storage at 15°C-28°C. Reseal bottle immediately after use. The material will cure if exposed to light. Do not use after the expiry date.

#### Disposal:

Dispose of the product in accordance with local regulations.

#### Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition, and a serious risk to public health that arises or could have arisen in association with the use of **V-Print splint comfort** must be reported to VOCO GmbH and to the responsible authority.

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VC 60 AA6126 E1 0725 99 © by VOCO