VOCO

V-Print splint

Instructions for use

Product description:

V-Print splint is a light-curing acrylic for the generative production of therapeutic splints as well as splints, auxiliary parts and functional parts for use in dental diagnostics.

Indications:

- Therapeutic splints
- Auxiliary parts and functional parts for diagnostics
- Bleaching splints (home bleaching)

Safety and precautions:

Avoid direct contact between the skin and liquid or uncured **V-Print splint**. The wearing of protective clothing is recommended. This should encompass gloves (nitrile), safety goggles and a long-sleeved lab coat. Further information on the handling of the product can be found on the safety data sheet.

V-Print splint can trigger sensitisation or irritation of the skin and/or respiratory organs. Repeated and/or prolonged exposure through contact and/or inhalation, including of dusts, can trigger irritations or even result in serious illness.

Note: Use separate material containers and cleaning baths for each printing material so as to rule out the possibility of cross-contamination.

Preparation:

V-Print splint is a light-curing acrylic designed for the generative production of articles using a 385 nm LED DLP printer. Information on approved printers is available from the corresponding device manufacturers and from **VOCO**.

The following construction conditions must be observed for indication-appropriate CAD constructions:

- Minimum wall thickness 1.5 mm / recommended wall thickness: 3 mm
- Round off internal edges on the outside.

Prepare a print job using slicing software. In addition to the material-dependent construction specifications in these instructions for use, please also observe the dependencies on the positioning, type of support and adaptation from our other documents for your construction.

V-Print splint can be printed in the following layer thicknesses:

- 50 µm, 75 µm and 100 µm

Processing:

The material container should be filled directly before printing is started. When doing so, ensure the material is as free from bubbles as possible and pay attention to the fill level.

Note: Do not shake the material before starting printing.

Start the print job observing the previously selected parameters.

We recommend allowing a dripping time of approx. 10 minutes once printing is complete. Then detach the printed objects carefully from the build platform.

Following this, the printed objects must be cleaned, dried and post-exposed in order to guarantee the requisite product characteristics. A detailed description of the two steps outlined above can be found under **Post-processing**. Post-process the printed object without delay so as to avoid distortion of its shape and ensure the required precision.

Remove the excess material from the material container following printing, so as to allow you to clean and check the container. The print material can be returned to its original packaging. If necessary, a separate container can also be used (HDPE, not light-transmissive, airtight).

Please ensure that the excess material intended for further use is free from contaminants and polymerised components. To this end, a stainless steel sieve or similar instrument should always be used when transferring material.

Post-processing:

Cleaning:

We recommend a solvent-resistant cleaning mat for cleaning.

Cleaning should be performed with isopropanol (≥ 98% purity) in an unheated ultrasound bath. The respective operating manuals and/or instructions for use provided by the manufacturers must be observed.

The printed objects must be cleaned in two or optionally in three stages.

Pre-cleaning:

- **Optional:** Pre-clean printed objects carefully for approx. 15 seconds by repeatedly immersing them in a beaker of isopropanol. When doing so, hold the component by a support using tweezers.
- Pre-clean the printed objects for 3 minutes in a reusable isopropanol ultrasound bath. Position uncleaned printed objects in the ultrasound bath in such a way that any openings face downwards.

Note: The bath's cleaning efficacy decreases the more it is used. Replace the bath when the cleaning efficacy decreases.

Main cleaning:

Clean the printed objects for a further 2 minutes in a fresh isopropanol ultrasound bath.

Then dry them carefully with compressed air.

If there are any resin residues on the printed object after the main cleaning or any exits from the undercuts during drying, the printed object can be immersed in the main cleaning bath again briefly before being dried again.

Preparation for post-exposure:

The support structures do not need to be removed before post-exposure.

Inconvenient support structures can be removed prior to the post-exposure carefully and without exerting pressure using a rotary instrument (e.g. cutting wheel) as close to the printed object as possible. Use a suction device. Remove any remaining acrylic dust carefully using compressed air. Then rinse printed objects for a few seconds with fresh isopropanol and remove any dust using a fine, clean brush if necessary. Dry the printed objects carefully again using compressed air.

We recommend waiting 15 minutes after the last contact with isopropanol before initiating the post-exposure. **Note:** Ensure that the printed objects are not exposed to any direct light during this time.

Post-exposure:

No protective gas atmosphere is required. However, ensure that the printed objects do not cast shadows over or touch each other, as this could result in incomplete post-curing.

The post-exposure can be performed with the following devices:

- **Otoflash G171 xenon photoflash unit** 2 × 2,000 flashes
- **LC-3DPrint Box¹ UV light box** 30 minutes

In printed objects post-exposed with the Otoflash G171, a cooling-off period of at least 2 minutes with the lid open must be included after the first 2,000 flashes.

Once the cooling-off period has elapsed, the printed objects should be turned and light-cured with a further 2,000 flashes.

Insufficient cooling can distort the shape of the printed object.

Finishing:

As a general rule, work with a low contact pressure and speed. This guarantees constant results and also reduces the risk of undesirable signs of processing.

Use a fine-toothed carbide bur or similar to sand the support stubs. The same bur can also be used for subsequent finishing of special structures.

For precise sanding, e.g. between the support strut and the printed object, it is recommended that the surface be processed in the corresponding areas with sandpaper, if necessary with different grain sizes. A corresponding result can also be achieved with coarser or finer silicone polishers.

To create a high-lustre surface, the object should first be pre-polished with pumice stone. Then remove any pumice dust and remnants under running water using a brush. Finally, use a buffer and high-lustre polishing compound without exerting excessive pressure to polish the object to a high-lustre finish.

Final cleaning:

Clean the object thoroughly. Begin by removing large residues with the steam jet. The subsequent final cleaning can be performed via brief storage in an unheated water ultrasound bath. A surfactant solution can be used instead of the water to remove oily and fatty substances.

Disinfection:

Objects produced from **V-Print splint** can be disinfected with alcohol- or aldehyde-based disinfectant solutions (e.g. ethanol, MD 520 from Dürr, Cavex Impre Safe from Cavex). Observe the manufacturer's instructions for use.

Information, precautionary measures:

V-Print splint may only be used intraorally in its completely cured state. Pay attention to the finishing process.

V-Print splint contains acrylates, BHT and TPO. **V-Print splint** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.

Contact between uncured **V-Print splint** and the skin/mucous membranes and eyes can cause mild irritation and should be avoided. The wearing of protective clothing is recommended. In addition, care must be taken not to inhale any vapours and/or dusts. The wearing of a suitable mask and/or the use of extraction systems is recommended. Further information on the handling of the product can be found on the safety data sheet.

Splints made from **V-Print splint** are unsuitable for the production of temporary restorations.

Storage:

Store at between **15 °C and 28 °C**. Seal bottle again immediately after use. The material will cure if exposed to light. Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. VOCO recognizes its responsibility to replace products if proven to be defective. VOCO does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith. Descriptions and data constitute no warranty of attributes and are not binding.

CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

For questions or comments, please call 1-888-658-2584.

Keep this material out of reach of children.

For dental use only.

An explanation of the symbols used in labeling can be found at www.voco.com/us/symbols

¹ Or a similar device. Please consult your device's manufacturer for further information.

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VC 60 006044 E1 0918 99 © by VOCO
VOCO GmbH Phone +49 (4721) 719-0