Grandio®SO x-tra

EN Instructions for use

[MD] EU Medical device

In accordance with DIN EN ISO 4049

Product description:

GrandioSO x-tra is a light-curing, radiopaque nano-hybrid restorative material. It is characterised by an increased depth of cure while at the same time showing reduced shrinkage stress, which means that it can be cured in layers of 4 mm thickness.

GrandioSO x-tra contains 86 % w/w inorganic fillers and is used with a dentine enamel bond.

Shades:
Universal (U), A1, A2, A3

Indications:

- Class I and II posterior restorations
- Base in class I and II cavities
- Class V restorations
- Locking, splinting of loose teeth
- Repairing veneers, enamel defects and temporary C&B-materials
- Extended fissure sealing
- Restoration of deciduous teeth
- Core build-up

Contraindications:
GrandioSO x-tra contains methacrylate and BHT. GrandioSO x-tra should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:
GrandioSO x-tra is suitable for use in all patients without any age or gender restrictions.

Performance features:
The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:
GrandioSO x-tra should only be used by a professionally trained dental practitioner.

Use:

Preparation:
Clean the teeth being treated. If applicable, mark occlusal contact points. Bring material to room temperature before use.

Shade selection:
Clean the teeth prior to shade selection. The shade is selected in comparison with the tooth while it is still moist.

Cavity preparation:
In principle, the cavities should be prepared according to the rules for adhesive restorative treatment; preparation should be minimally invasive in order to protect healthy dental hard tissue. Clean and dry the cavity before treatment. Non-carious cervical lesions do not have to be prepared; thorough cleaning is sufficient here.

Creating a dry working field:
Ensure that the work area is sufficiently dry. Use of a rubber dam is recommended.

Matrix placement:

Use of matrices is recommended for cavities with an approximal portion. Translucent matrices that are wedged in the approximal area are advantageous. Minimal separation facilitates the shaping of the approximal contact and insertion of the matrix.

Pulp protection:
A suitable pulp-protective liner that can be covered with a stable cement, if necessary, should be applied in close proximity to the pulp.

Bonding material:

GrandioSO x-tra is applied with a dentine enamel bond, using the adhesive technique. Any light-curing bonding materials can be used. Observe the relevant instructions for use during preparation (etching technique) and processing.

Application:
Apply GrandioSO x-tra in layers that are a maximum of 4 mm thick, adapt with a suitable instrument and light-cure afterwards.

Light-curing:

Conventional polymerisation devices with a wavelength in the range of 400 - 500 nm are suitable for light-curing the material.

Depending on the light output and shade the curing time for increments of 2 - 4 mm is:

LED-/halogen light	Universal (U)	A1, A2, A3
≥ 1000 mW/cm²	10 s	20 s
≥ 500 mW/cm²	20 s	40 s

The curing time for increments of **up to 2 mm** is 10 s for all shades when using an LED-/halogen light with a minimum light output of 500 mW/cm².

Bring the light aperture as close as possible to the light-curing unit as otherwise the curing will be less effective. Inadequate curing can cause discolouration and complaints.

Finishing:
The restoration can be finished and polished immediately after removing the matrices (e.g., fine or extra-fine diamond burs, polisher). The tooth should be fluoridated as a final step.

Warnings, precautionary measures:

- Phenolic substances, especially products containing eugenol and thymol, lead to curing failures in restorative composites. Consequently, the use of zinc oxide/eugenol cements or other materials containing eugenol should be avoided in combination with restorative composites.
- Each GrandioSO x-tra Caps unit should be used for one patient only and should be applied exerting uniform but not excessive pressure.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):

Barium aluminium borosilicate glass, silicon dioxide, ytterbium trifluoride, BisGMA, BisEMA, methacrylated silica, TCDDMA, initiators, stabilisers, colour pigments

Storage instructions and application method:

Storage at 4°C-23°C. Seal the syringes immediately after removing the material in order to avoid exposure to light and subsequent polymerisation. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of GrandioSO x-tra must be reported to VOCO GmbH and the responsible authority.

Note:

The Summary of Safety and Clinical Performance of GrandioSO x-tra can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at www.voco.dental.

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