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VOCO


Perfect Bleach®

EN Directions for fabrication of the vacuum formed tray


The correct fit of the tray is important for successful treatment. Rough, sharp or overlapping edges and cracks of the tray can cause injury or irritation of the oral mucosa. In the event of subsequent breakage, there is a risk of aspiration for the user. In case of known allergic reaction to the plastic EVA (ethylene vinyl acetate), avoid skin contact. For single use only (thermoforming). Hollow spaces in the occlusal or incisal areas, and margins that do not fit tightly may cause the whitening gel to leak and irritate the gingiva if the tray is overfilled.

Fabricating the tray:

1. Make a plaster model from the impression with the model being trimmed to approx. 4 mm below the gingival line. The underside of the model should be parallel to the occlusal plane.
2. Allow the model to dry for at least 2 h. The drying process can be accelerated with heat, e.g. in a drying cabinet or using warm airflow. The **Block Out Gel LC** adheres better to dry plaster, which is advantageous when vacuum forming the tray.
3. Apply the **Block Out Gel LC** on the model to the labial surfaces of the teeth to be treated, maintaining a minimum distance of 0.5 to 1.0 mm from the gingival line (Fig. 1). Do not cover the occlusal surfaces or incisal edges with **Block Out Gel LC**. Apply **Block Out Gel LC** in layers that are approx. 0.5 mm thick. Carefully remove any excess **Block Out Gel LC**!

4. Polymerise **Block Out Gel LC** with a curing light for composites or in a light-box. Cure **Block Out Gel LC** on each tooth for 20 s (light-box: follow manufacturer's instructions for curing time). Use a cloth moistened with ethanol to remove the inhibition layer after curing. Do not use acetone or other agents that could partially dissolve the **Block Out Gel LC**!
5. Place the model thus prepared in the vacuum forming unit and vacuum form the enclosed foil over it.
6. Heating time
The vacuum forming foil (ERKOFLEX 95) is an ethyl vinyl acetate foil that is 1.50 mm thick. This is vacuum formed with the spacer foil that points in the direction of the model. Heat the sheet until it sags slightly and becomes transparent. Set the vacuum forming unit accordingly. The spacer foil should only be removed after vacuum forming or finishing.
7. After the tray has cooled, remove the tray whilst being careful not to damage the model. Remove the rough excess first and then cut along the gingival line (Fig. 2).

8. Place the tray thus prepared on the model, check the edges and make any necessary adjustments. To avoid injuries to the mucosa, polish any protruding or sharp edges using suitable grinding and polishing discs. To adjust areas that protrude or are too short, heat the areas (e.g. small burner) and correct with a moistened finger.
9. Clean and disinfect the tray after completion.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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VOCO GmbH Phone +49 (4721) 719-0