VOCO V-Print® c&b temp

EN Instructions for use
MD EU Medical device

Product description:

V-Print c&b temp is a light-curing resin for the generative production of highly aesthetic temporaries / long-term temporaries in the CAD/CAM technology. The maximum clinical wearing time is 12 month.

V-Print c&b temp is a thixotropic material with an anorganic filler content of 26 % w/w.

Thanks to the Viscosity Change Technology, multiple swirling is sufficient to fill the tray with **V-Print c&b temp** efficiently.

Thanks to the composite technology, the material displays good abrasion resistance as well as high fracture resistance and is thus also ideal for long-term temporaries. **V-Print c&b temp** is fluorescent. The restorations can be inserted with a temporary cement or an adhesive technique.

Indications:

Long-term temporaries like crowns, bridges and mock-ups

Contraindications:

V-Print c&b temp contains (meth)acrylates and phosphine oxide. **V-Print c&b temp** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

V-Print c&b temp is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

V-Print c&b temp should only be used by a professionally trained dental practitioner.

Shade selection:

Use the VITA® shade system to determine the shade against the cleaned but unprepared tooth prior to anaesthesia and preferably in daylight conditions.

Core and cavity preparation:

As a general rule, the core/cavity preparation should be performed in accordance with the conventional rules for all-ceramic restorations. That means rounding off interior corners and edges and using a shoulder preparation with rounded interior edges or a chamfer preparation.

Hardware and software requirements:

CAD software ¹ dental scanner	Software for the planning and design of crowns and bridges. The software and dental scanner must satisfy local and current medical device specifications and allow for issuance of the patient-specific design as an STL data set.
CAM software	Software for preparation of the print order. The part will not be modified during this process. Structures that facilitate the 3D printing are simply created. For example: - Autodesk Netfabb version 2020 or later for SolFlex 3D printing.

¹The designation **Software as Medical Device SaMD** includes standalone (autonomous) software that is a medical device (MD) and not part of one.

Manufacturing equipment	For example: VOCO SolFlex 170 VOCO SolFlex 350 VOCO SolFlex 650 VOCO SolFlex 170 HD
Post-curing devices	For example: Otoflash G171

See also: accompanying list of resources or www.voco.dental/3dprintingpartners

All manuals and/or operating instructions for the respective programmes, and for device, materials and/or parts manufacturers, which are required for the manufacturing process, must be observed.

Clarify ahead of time whether the programmes, devices and/or objects that you intend to use have been designed and approved for the corresponding applications.

CAUTION: Non-authorised changes to the process equipment, parameters, or software could result in the **V-Print c&b temp** end object not satisfying specifications.

Use:

Preparation:

For an indication-appropriate CAD construction, the following design conditions must be observed:

Minimum wall thickness:

- occlusal 1.5 mm
- circular 1.0 mm

Bridges with a maximum of one pontic can be produced.

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- Note: Design of the cross-sectional areas: height ≥ width

Prepare a print job using CAM-software.

Processing:

Note: Use separate material containers and cleaning baths for each printing material, in order to prevent cross contamination.

Multiple swirling of the bottle improves the flow properties of the material and should be done immediately before printing.

It is important to ensure that the material is free of bubbles to the extent possible, and filled to the fill level mark.

Start the print job observing the parameters that you previously selected.

Once the printing process has ended, a dripping time of approximately 10 minutes is recommended. Next, carefully detach the printed objects from the build platform.

In the following steps, the printed objects will need to be cleaned, dried and post-exposed, in order to guarantee the required product characteristics. A detailed explanation of the steps outlined above can be found under **Post-processing**.

Recommendation: After completing your work, transfer the remaining material from the material tray to the original container (use a stainless steel sieve if necessary). This serves to check the material tray and also enables optimum storage of the print material.

Post-processing:

Cleaning

Remove unpolymerised resin residues on the print objects using an isopropanol (purity ≥ 98 %) soaked brush.

Next, the printed objects must be dried carefully using compressed air. If there is any resin residue on the printed object after the final cleaning, or if residue escapes from the undercut when drying, the printed object should again be cleaned with an isopropanol (purity ≥ 98 %) soaked brush.

Preparation for post-exposure:

Obstructive support structures should be removed before the post-exposure process by using a rotary instrument as close to the printed object as possible, carefully and without exerting pressure. Use a suction device. Remove remaining resin dust carefully with compressed air and, if necessary, with a brush and isopropanol (purity ≥ 98 %). Carefully dry the printed objects once again with compressed air.

Post-exposure:

Conduct the post-exposure a minimum of 15 minutes after the most recent contact with isopropanol. It is important to ensure that the printed objects do not overlap or contact each other, as post-exposure would be negatively affected by the shadows that are cast.

Post-exposure can be conducted using the following devices:

Post-exposure device	Programme	
For example: Xenon photoflash unit Otoflash G171	2 x 2000 flashes	After 2000 flashes, observe a cooling phase of at least 2 minutes with open lid. Next, turn over and light-cure with another 2000 flashes.

See also: accompanying list of resources

Further print objects should only be post-exposed after the device has cooled down for 10 minutes in order to avoid discolouration of the material.

Finishing/polishing:

Use, for example, a fine-toothed carbide cutter to grind down the support attachments. This can also be used for subsequent elaboration of special structures

In order to achieve form-fit grinding, e.g. between the support base and the print object, it is recommended to process the surface in the corresponding area with coarser or finer rubber polishers. A corresponding result can also be achieved with sandpaper, possibly with different grit sizes.

Polish restoration chairside with conventional composite polishers. A multiphase polishing system is advantageous.

A polishing paste can be used latside in combination with goat's hair brushes and cotton/leather buffers.

Use only light pressure and wiping motions to control the removal. The instructions for use from the manufacturers must be observed and followed.

Customisation:

For a highly aesthetic temporary restoration, the restorations can be customised, characterised or repaired at any time using a composite/ORMOCER®. Roughen the restoration surface via grinding or sandblasting (Al₂O₃ 50-100 µm, 1-2 bar). Remove abrasive material residues/dust carefully with an ultrasonic bath (70% ethanol) or steam cleaner. Then dry the restoration with air. Apply a suitable adhesive system (e.g., **Futurabond U**) in accordance with the instructions for use. Using **GrandioSO**, **Flow** or **Heavy Flow**, for example, in combination with **FinalTouch**, you can customise the restorations quickly and simply with purely light-curing techniques.

The instructions for use from the manufacturers must be observed and followed.

Luting:

Preparing the restoration

For an optimal bond, roughen the luting surface of the restoration using aluminium oxide (50-100 µm) at 1-2 bar or roughen with a fine HM milling device. Use a suction device to remove the dust produced. Remove abrasive material residues carefully with an ultrasonic bath (70% ethanol) or steam cleaner. Then dry the restoration with air. Final cleaning with medical alcohol is possible. The instructions for use must be observed and followed.

Temporary luting for crowns and bridges

Lute **V-Print c&b temp** with a temporary cement (e.g., **Provicol QM**). **Note:** Use a temporary, eugenol-free cement if the permanent restoration is to be subsequently luted with an adhesive.

Adhesive luting for crowns and bridges

For a wearing period of > 30 days, the restoration can be luted using a composite-based adhesive luting material (e.g., **Bifix QM**). The instructions for use must be observed and followed.

- Note:**
- Adhesive luting causes time-consuming removal of the temporary restoration.
 - Regular check-ups and follow-up examinations are required for long-term use.

Warnings, precautionary measures:

- Only use **V-Print c&b temp** intraorally in a fully cured state. Pay attention to the finishing process.
- Contact between uncured **V-Print c&b temp** and the skin/mucous membranes and eyes can cause irritation and should be avoided.
- The wearing of protective clothing is recommended. Furthermore, it is important to ensure that no vapours and/or dusts are inhaled. The wearing of a suitable mask and/or the use of suction devices is recommended.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for the intended purpose.

Storage:

Storage at 15°C-28°C. Reseal bottle immediately after use. The material will cure if exposed to light. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition, and a serious risk to public health that arises or could have arisen in association with the use of **V-Print c&b temp** must be reported to VOCO GmbH and to the responsible authority.

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VC 60 AA6897 E1 0725 99 © by VOCC