VOCO

VOCO Ionofil Molar AC/Quick

GBIn accordance with
DIN EN ISO 9917-1 : 2007**CE** 0482

Instructions for use

Product description:

VOCO Ionofil Molar AC/Quick is a radiopaque glass ionomer cement for fillings in application capsules (AC) for direct application.

VOCO Ionofil Molar AC/Quick adheres chemically to enamel and dentine and continuously releases fluoride ions.

VOCO Ionofil Molar AC is available in shades A1, A2 and A3.

VOCO Ionofil Molar AC Quick is available in shades A1, A2 and A3.

From every application capsule 0.15 ml cement can be extracted.

Indications:


- Restoration of V-shaped defects and cervical enamel or root erosions including class V cavities
- Restoration of deciduous teeth (especially class I)
- Fillings of class II cavities (not occlusion-bearing) by means of tunnel preparation
- Build-up fillings and linings, core build-up
- Temporary long-term treatment of class I and II cavities

Method of application:

Preparation:

The adhesion of **VOCO Ionofil Molar AC/Quick** to enamel and dentine prevents the formation of marginal gaps and creates retentions. Undercuts give additional hold to the filling. The layer thickness of fillings should not be less than 1 mm. Avoid feather-shaped margins. In case of doubt choose conventional preparation. Clean the cavity (e.g. 3 % hydrogen peroxide), thoroughly rinse with water and dry (warm air). Do not overdry dentine. If the remaining dentine is thicker than 1.5 mm additional lining is not necessary. Cover areas in close proximity to the pulp with calcium hydroxide to avoid pulp irritations which may occur in cases of insufficient dentine thickness.

The capsule must be activated, mixed and opened before application.

Activation:

For activation, the liquid from the sachet underneath the clamp has to be emptied as completely as possible into the inner chamber of the capsule. This is done with a capsule activator. Place the capsule in the activator so that the application cannula protrudes from it and the clamp faces the movable part of the activator.

Hold the lever for 2 s in completely compressed position.

Mixing:

Fix the activated capsule in a capsule mixer (4,000 - 4,500 oscillations per min) and mix for approximately 10 s. Immediately remove the capsule from the mixer. The use of capsule mixers of other manufacturers may result in shorter or longer mixing times due to different mixing frequencies.

Opening:

Open the twist cap with application cannula with a half turn to the right until you notice a stop (see stamp on protective cap).

Application:

Insert the opened capsule into the holding bracket of an applicator. By turning the capsule in the holding bracket the cannula can be aligned in the desired direction. Compress the lever slowly and evenly to apply the material directly into the cavity within 30 s after completion of mixing. Do not apply excessive pressure.

VOCO Ionofil Molar AC:

Mixing time	working time	setting time
10 s	2 min	4 min

VOCO Ionofil Molar AC Quick:

Mixing time	working time	setting time
10 s	1.5 min	2.5 min

After mixing, **VOCO Ionofil Molar AC/Quick** has a working time of approximately 1.5 min (room temperature). During the placing of a filling the working area must be kept dry. The use of a rubber dam is recommended.

Glass ionomer cements are susceptible to humidity (breath, saliva) or, respectively, dehydration during the setting time, especially directly after application. If possible, cover immediately with a strip or matrix or apply a protective varnish (see instructions for use **Final Varnish LC**).

Finishing and polishing (flexible polishing discs) of fillings should be conducted no earlier than 7 min (**VOCO Ionofil Molar AC**) or, respectively, 2.5 min (**VOCO Ionofil Molar AC Quick**), after placement of the filling. Recommendation: Re-apply protective varnish afterwards.

The best results of fine finishing will be achieved approximately 24 h after placement of the filling.

Additional advice/precautions:

VOCO Ionofil Molar AC/Quick contains polyacrylic acid and fluoroalicates. **VOCO Ionofil Molar AC/Quick** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.

VOCO Ionofil Molar AC/Quick should be mixed only immediately before application in order to have the full working time available for inserting and shaping and to avoid cracks from working the material too late.

Ensure that the application capsule is mixed immediately after activation.

Storage instructions and application method:

Store **VOCO Ionofil Molar AC/Quick** in a dry place at temperatures of 15 °C - 23 °C. Only remove it from the foil immediately before use. Do not use after expiry date.

Tightly close **Final Varnish LC** container immediately after use. Store at room temperature and protected from light.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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