VOCO

Fissurit/Fissurit F

GB

CE 0482

In accordance with DIN EN ISO 6874

Instructions for use

Product description:

Fissurit is a transparent light-curing fissure sealant.
Fissurit F is a white, light-curing fissure sealant with fluorides.

Indications:

- Sealing/filling pits and fissures as well as small cavities
- Sealing/facing of damaged enamel surfaces
- Anchoring orthodontic appliances and covering caries predilection sites during orthodontic treatments
- Sealing composite or cement fillings (protective layer against humidity)
- Sealing carious lesions to stop them
- Restoring small defects in composite and amalgam fillings
- Sealing deciduous teeth

Directions:

Fissures should be sealed as soon as the eruption of the tooth allows treatment and a dry work field can be secured. After drying, thoroughly clean the surfaces to be sealed (rotating brush/cleaning paste), spray and dry. In case of incipient carious lesions or unclear diagnosis extend the fissure with a fine instrument.

Etching:

Apply **Vococid** etching gel (blue colour for control), work in with a short-haired disposable brush and let act for 15 - 30 s according to the condition of the enamel.
Thoroughly rinse and dry. The etched enamel should have an opaque and matt-white appearance. Avoid subsequent contamination. Please refer also to the **Vococid** instructions for use.

Application:

Carefully apply and disperse **Fissurit/Fissurit F** with a brush or any other suitable instrument. Remove any air bubbles (e.g. with a fine probe).

Syringe application:

Remove the protective cap from the syringe. Place a new cannula on the **Fissurit F** syringe. To avoid air bubbles, hold the syringe vertically and press the material up to the tip of the cannula.
Fissurit F homogeneously flows from the cannula upon light to medium pressure.

In order to ensure the functioning of the non-dripping syringe do not pull back the plunger during or after use.

Apply **Fissurit/Fissurit F** from the margin into the fissure or cavity and let penetrate for 15 - 20 s. Remove any air bubble (e.g. with a fine probe).

Light-curing:

Polymerise **Fissurit/Fissurit F** for 20 s per surface and check the sealing. To light-cure this material conventional polymerisation devices are suited. The light output should be a minimum of 500 mW/cm² on halogen polymerisation devices and 300 mW/cm² on LED devices.
Re-seal deficit areas, remove excess cautiously (e.g. finishing diamond) and polish if necessary (rubber polisher). Finally apply a fluoride varnish on the treated and neighbouring, preferably on all teeth.

Advice and precautionary measures:

A work field which cannot be dried sufficiently involves the risk of a (partial) premature loss of the sealing. In these cases, the application of a fluoride varnish or sealing with a glass ionomer cement is recommended until a perfect sealing is possible.

If the etched surface is contaminated (e.g. by saliva), the adhesion of **Fissurit/Fissurit F** may be reduced. In case of contamination, clean, etch again for approx. 10 - 20 s, spray and dry.

Adhesion may be improved after etching by using a bonding (difficult drying or extended sealing).

Expose **Fissurit/Fissurit F** to operating light only a short time as otherwise polymerisation will start (viscosity increases).

Check sulcus after application and remove all material residue.

Fissurit/Fissurit F contains Bis-GMA, diurethane dimethacrylate, BHT, benzotriazolderivate. **Fissurit F** also contains sodium fluoride. **Fissurit/Fissurit F** should not be used in case of known hypersensitivities (allergies) to any of these ingredients. Avoid contact with the mucosa as far as possible. Do not use cleaning and drying agents containing phenolic substances (e.g. eugenol, thymol) as they inhibit a complete polymerisation.

Carious lesions (e.g. around fillings) can be stopped by corresponding sealing if renewal of the filling is not desired or possible. For this purpose, clean the corresponding areas, prepare as minimally as possible and etch in the enamel area. Sealing the complete filling is recommended in these cases. In cases of orthodontic treatments in areas of reduced oral hygiene, sealing those areas can provide good and permanent protection. Fluoridation should be continued in all cases to protect the unsealed tooth areas from caries.

Storage:

Immediately and carefully close the **Fissurit/Fissurit F** syringe after use. Remove the cannula of the syringe and attach syringe cap. Before next use of the syringe replace the direct application tip with a new one. Store **Fissurit/Fissurit F** at temperatures between 4 °C - 23 °C and protect from light. Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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VC 60 AA1080 E1 0415 99 © by VOCO
VOCO GmbH