VOCO**Ufi Gel® hard**

Instructions for use

**Product description:**

Ufi Gel hard is a cold-curing permanently hard relining material for dentures and is polymethacrylate-based. It is simple and quick to use for direct as well as indirect relinings.

The **Ufi Gel hard** liquid contains no methylmethacrylates, thus minimising the risk of allergies and irritations of the mucosa. The set additionally contains a conditioner to achieve a permanent and stable bond.

Fields of application:

- hard, permanent, total or partial relining to restore the function of partial and total dentures
- lengthening of denture margins

Application:**1. Preparation of the denture**

Check occlusion and carry out corrections, if necessary. Clean denture thoroughly with a brush and dry. Roughen all areas to be relined incl. buccal and labial margins with a suitable bur. Then clean and dry. Insulate areas not to be brought into contact with **Ufi Gel hard** (e.g. artificial teeth, labial and buccal areas below roughened surfaces), e.g. with vaseline. **Clasps, anchors and attachments of partial dentures have to be blocked out with thin-flowing silicone or wax with a low melting point.**

2. Application of the conditioner

Apply conditioner with enclosed brush on all surfaces to be relined and let dry in the air (approx. 30 s). First coat the labial and buccal surfaces, then the base of the denture. The conditioner remains effective for about 10 min after application.

Attention: Close the bottle tightly immediately after use because of high volatility.

The brush can be cleaned with alcohol, for example.

3. Dosing and mixing

Prior to first use, exchange the transportation cap against the dropper. **Ufi Gel hard** is mixed in a ratio of 1 ml liquid to 3 ml (= 1.8 g) powder. This corresponds to 2 graduation marks of the dropper to 1 graduation mark of the glass cylinder.

2 graduation marks liquid : 1 graduation mark powder

A thicker consistency, e.g. for lengthening the margins of the denture, is achieved by increasing the amount of powder.

Extract liquid with the dropper and drop it into a PP mixing cup. Shake the powder briefly and dispense it into the glass cylinder. For exact dosing smooth the surface of the powder by lightly tapping the side of the glass cylinder. Place powder into the liquid and mix carefully with a plastic spatula until homogeneous (approx. 10 s). Prevent formation of bubbles by stirring slowly and along the side of the mixing cup. Make any potential bubbles rise to the surface by tapping the mixing cup. Let the material soak until a workable consistency is achieved (approx. 1.5 min after commencement of mixing).

4. Application of Ufi Gel hard

Apply the material evenly onto the prepared margin and/or base of the denture with a plastic spatula. Prevent build-up of excess material, remove with a suitable instrument if necessary. Re-insert denture and have the patient exert slight pressure in terminal occlusion for 1 min. Then carry out functional, chewing and swallowing movements for 2 min. Ensure that no material flows into the throat when relining an upper denture, especially at the transition from hard to soft palate (A-line).

4.1 Partial and total denture with undercuts

Remove denture 5 min after commencement of mixing and remove excess material immediately with scissors or a scalpel. Insert again into the mouth for a final occlusion check and allow to cure for about 2 - 3 min.

Instead of re-inserting the denture, curing can also be completed in warm water, e.g. in a pressure pot at approx. 40 °C. **Do not let the material cure in contact with air** since atmospheric oxygen will cause an uncured inhibition layer on the surface.

4.2 Total dentures without undercuts

Excess material can be removed intra- or extraorally before final polymerisation.

Intraorally: remove excess material after 5 min (from beginning of mixing) with a suitable instrument. Let the denture cure for further 2 - 3 min in the mouth until **Ufi Gel hard** is completely cured. Extraorally: see section 4.1

5. Finishing the relining

The relined denture can be finished and polished with the usual instruments (tunsten carbide bur, rubber polisher, polishing disc).

beginning of mixing	2 min	mix and apply
insert denture into the mouth	1 min	exert pressure
	2 min	functional movements
remove denture if necessary	1.5 min	remove excess material
insert denture into the mouth	2 min	curing
remove denture	1 min	finishing, polishing
	9.5 min	

extraorally intraorally intra- or extraorally

Indications, precautionary measures:

- Store bottles of liquid and conditioner carefully closed and in upright position.
- Avoid contact of liquid or conditioner with skin. Rinse sites of contact thoroughly with water and soap.
- In case of contact with the eyes rinse thoroughly with water and consult an ophthalmologist.
- **Ufi Gel hard** contains hydroxyethyl methacrylate, benzoyl peroxide, acetone. **Ufi Gel hard** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.
- The indicated time periods refer to a room temperature of 25 °C and a temperature of 35 °C in the mouth. Polymerisation will be slower at lower temperatures, higher temperatures will accelerate polymerisation.
- Excessively long or violent mixing may lead to air bubbles and an un-homogeneous consistency, which produces a rough surface.
- Use the mixing cup several times.

Storage:

Store at temperatures of 4 °C to 23 °C. Do not use after expiry date.

Cleaning indications:

Dentures relined with **Ufi Gel hard** can be cleaned with the usual cleansers and procedures. This applies to domestic as well as professional cleaning.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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VC 60 AA2211 E1 0717 99 © by VOCO
VOCO GmbH