
VOCO

CE 0482

Bifix® Hybrid Abutment**EN Instructions for use**
MD EU Medical device**Product description:**

Bifix Hybrid Abutment is a radiopaque, composite-based, dual-curing luting system for the permanent luting of abutment crowns and mesostructures made of metal, ceramics or polymers (composite or PMMA) on prefabricated adhesive or titanium bases (extraoral use), customised abutments made of titanium or zirconium dioxide or dental hard tissue (combination solutions).

Bifix Hybrid Abutment is available in two highly opaque, highly masking shades (white/opaque and universal/opaque) and in a translucent shade. Colour-matched try-in pastes are available to allow shade selection in advance.

Indications:

Permanent extraoral and intraoral luting of superstructures (incl. mesostructures) made of metal, ceramics, zirconium dioxide and polymers (composite, hybrid ceramics or PMMA) on prefabricated adhesive/titanium bases (extraorally only), on customised titanium and zirconium dioxide bases (both extraorally and intraorally) as well as on dental hard tissue.

Contraindications:

Bifix Hybrid Abutment contains methacrylates, benzoyl peroxide, BHT and amines. **Bifix Hybrid Abutment** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

Bifix Hybrid Abutment is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

Bifix Hybrid Abutment should only be used by a professionally trained dental practitioner.

Use:**1. Shade selection (optional)**

For intraoral use, remove temporary restoration if present. Clean the abutment or tooth, rinse thoroughly and gently dry with air.

The shade effect of **Bifix Hybrid Abutment** on the restoration can be simulated by inserting the definitive restoration with **Bifix Hybrid Abutment Try-In paste** in the intended shade. The shade of **Bifix Hybrid Abutment Try-In paste** corresponds to the shade of the respective cured **Bifix Hybrid Abutment**. To check the shade, apply **Bifix Hybrid Abutment Try-In paste** to the inner surface of the workpiece (type 41 application tip). Then, insert the workpiece carefully, using slight pressure. Do not perform an occlusion check (risk of fracture of the restoration).

After checking the shade effect, remove the workpiece, rinse **Bifix Hybrid Abutment Try-In paste** off the adhesive base, abutment or tooth and off the restoration thoroughly with water spray and then dry. Take care not to overdry the dentine.

2. Preparation for extraoral use

For detailed information on the correct use, please also consult the instructions for use of the respective bonding agent, as well as the manufacturer's specifications on preparation of the respective restoration materials.

2.1 Preparation of adhesive or titanium base or of customised titanium or zirconium dioxide abutment

Protect the implant-titanium base connecting surface and the screw channel with wax or silicone, or seal with Teflon tape or a foam pellet. Sandblast the joining surface with a suitable abrasive material (observe the manufacturer's instructions accordingly). Remove the wax and abrasive material residue thoroughly using an instrument, steam jet cleaner and/or ultrasonic bath. Then dry with air.

Apply a suitable bonding agent to the thus prepared surface (e.g. **Ceramic Bond**: apply, allow to react for 60 seconds and dry thoroughly with air). Do not touch or otherwise contaminate the conditioned surface.

To protect the implant-titanium base joining surface and the screw channel prior to the cementation (see 5), apply wax or silicone again or seal with Teflon tape or a foam pellet.

2.2 Preparation of workpieces made of silicate ceramic

Clean the workpiece in an ultrasonic bath and using a steam jet cleaner, then dry with air. Wax can be applied to protect the outer surfaces. Treat the joining surfaces with a hydrofluoric acid etching agent (observe the manufacturer's instructions), then rinse thoroughly with running water and dry with air.

Apply a suitable bonding agent to the inside of the workpiece (e.g. **Ceramic Bond**: apply, allow to react for 60 seconds and dry carefully with air). Do not touch or otherwise contaminate the conditioned surface.

2.3 Preparation of metal, composite and zirconium dioxide workpieces

Wax can be applied to protect the outer surfaces. Sandblast the joining surface with a suitable abrasive material (observe the manufacturer's instructions). In order to control the blasting process, the joining surface can be colour-marked before sandblasting. Remove any abrasive material residue carefully in the ultrasonic bath or using a steam jet cleaner and then dry the joining surface with air.

Apply a suitable bonding agent to the inside of the workpiece (e.g. **Ceramic Bond**: apply, allow to react for 60 seconds and dry carefully with air). Do not touch or otherwise contaminate the conditioned surface.

2.4 Preparation of PMMA workpieces

Observe the manufacturer's instructions if sandblasting is required.

Carefully remove any abrasive material residues, e.g. using the ultrasonic bath, and dry the joining surface with air. Next, condition the joining surface with a suitable PMMA bonding agent, following the instructions for use. Do not touch or otherwise contaminate the conditioned surface.

3. Preparation for intraoral use

For detailed information on the correct use, please consult the instructions for use of the respective bonding agent, as well as the manufacturer's specifications on preparation of the respective restoration materials. The working area must be kept dry. The application of a rubber dam is recommended.

3.1 Preparation of customised titanium/zirconium dioxide abutment

If not already prepared at the dental laboratory, sandblast the joining surface with a suitable abrasive material. Follow the instructions for use of the respective manufacturer. Carefully remove any abrasive material residues, e.g. in an ultrasonic bath, and dry the adhesive surface.

Apply a suitable bonding agent to the thus prepared surface (e.g. **Ceramic Bond**: apply, allow to react for 60 seconds and dry carefully with air). Do not touch or otherwise contaminate the conditioned surface.

3.2 Preparation of workpieces made of silicate ceramic

Prior to the preparation, the inner joining surface of the work to be inserted must be conditioned with a hydrofluoric acid etchant. This is usually done in advance in the dental laboratory. If the restoration has not yet been conditioned, prepare it using a suitable etchant (observe the manufacturer's instructions), then rinse thoroughly with running water and dry with air.

Apply a suitable bonding agent to the inside of the workpiece (e.g. **Ceramic Bond**: apply, allow to react for 60 seconds and dry carefully with air). Do not touch or otherwise contaminate the conditioned surface.

3.3 Preparation of metal, composite and zirconium dioxide workpieces

If not already done at the dental laboratory, sandblast the joining surface with a suitable abrasive material (observe the manufacturer's instructions). Carefully remove any abrasive material residues, e.g. in an ultrasonic bath, and dry the adhesive surface. Apply a suitable bonding agent to the inside of the workpiece (e.g. **Ceramic Bond**: apply, allow to react for 60 seconds and dry carefully with air). Do not touch or otherwise contaminate the conditioned surface.

3.4 Preparation of PMMA workpieces

If not already done at the dental laboratory, observe the manufacturer's instructions if sandblasting is required. Carefully remove any abrasive material residues, e.g. in an ultrasonic bath, and dry the joining surface. Then condition the joining surface with a suitable PMMA bonding agent, following the manufacturer's instructions for use. Do not touch or otherwise contaminate the conditioned surface.

3.5 Dental hard tissue

Remove temporary restoration if present. Clean the tooth, rinse thoroughly and lightly dry with air. **Bifix Hybrid Abutment** must not be applied to the exposed pulp or dentine close to the pulp. Use a suitable base fill (e.g. calcium hydroxide product and/or glass ionomer cement) to protect the pulp.

Note:

The working area must be kept clean and dry. The use of a rubber dam is recommended. The shade appearance can be affected by the drying of teeth under a rubber dam.

Bifix Hybrid Abutment must be used with a bonding agent suitable for dual-curing materials (e.g. **Futurabond U**). Observe the adhering instructions for use.

4. Disinfection/sterilisation of extraorally adhered superstructures

Prior to its insertion the adhered superstructure should be disinfected with an agent containing ethanol or can be sterilised (observe the manufacturer's instructions for individual components).

5. Luting of workpiece with Bifix Hybrid Abutment

If stored in a refrigerator, bring the material to room temperature before use. Before luting the workpiece, isolate any attachments or other connectors, using Vaseline or similar.

If necessary, block out the workpiece's screw channel.

Mixing tips:

Type 9 or 11. The type 11 mixing tip can be equipped with a type 4 or 1 intraoral tip. Remove the protective cap from the QuickMix syringe. Prior to every use, express material until it flows evenly from both openings. Then attach a mixing tip and lock (with a 90° turn clockwise).

The warranted product characteristics can only be achieved when using the supplied original mixing tips, which are also available as accessories. The attached mixing tip is intended for single use only. After use, store the QuickMix syringe tightly sealed. Check that material flows freely through the openings before reusing.

Do not mix the material until immediately before luting. The working time at room temperature (23°C) is approximately 2 minutes from the start of mixing. The material is automatically mixed without bubbles or errors in the tip as it is dispensed and can be applied directly to the prepared contact surfaces. Insert the restoration and affix uniformly by applying slight pressure. **The extraoral setting time is approximately 7 minutes at room temperature (23°C); the intraoral setting time is approximately 4 minutes (37°C).** Larger quantities of excess material can be removed with a foam pellet, disposable brush or floss. Dental floss is recommended in interdental spaces. However, small quantities of excess material should initially be left on the luting joint (see 6.) so as to avoid deficits.

6. Excess removal

As is the case for all resin-based materials, surfaces in contact with atmospheric oxygen do not cure completely, like **Bifix Hybrid Abutment** surfaces. The following methods are available to avoid deficits in the attachment joint:

6.1 Self-curing mode

Two variants are available:

a. Leave small amounts of excess material in place and remove after complete curing, during subsequent finishing procedure.

b. Mask restoration margins with glycerine gel in the gel phase, after removal of excess material, and wait for complete curing. Once polymerisation is complete, rinse off the glycerine gel with water.

6.2 Light-curing

With small amounts of excess material, removal can be performed after light-curing. To do so, use an LED/halogen lamp with a light output of at least 1,000 mW/cm² to polymerise the areas of the luting joint for 1-2 seconds for the translucent shade and at least 10 seconds per segment for opaque shades. Hold the tip of the lamp as close as possible to the surface to be cured. The restoration must be fixed in the desired position during the light polymerisation. Mask restoration margins with a minimum of 10 sec per segment, and polymerise once again for a minimum of 10 sec per segment. Once polymerisation is complete, rinse off the glycerine gel with water.

7. Finishing

Removal of the excess material is possible immediately after polymerisation. Carefully remove cured excess material using a finishing diamond or suitable silicone polisher. Finish and polish the approximal areas with finishing and polishing strips.

Warnings, precautionary measures:

– Only use suitable bonding agents.

– Some enamel-dentine bonding agents are not compatible with dual-curing composites. In case of doubt, observe the manufacturer's instructions for use for the respective bonding agent. **Futurabond U** is ideal for an optimal adhesive bond.

– In case of eye contact, rinse well with plenty of water and consult an ophthalmologist.

– Phenolic substances, especially products containing eugenol and thymol, interfere with the curing of filling composites. For this reason, the use of zinc oxide eugenol cements or other materials containing eugenol in combination with **Bifix Hybrid Abutment** should be avoided.

– Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):

Barium aluminium borosilicate glass, titanium dioxide, HEDMA, BisGMA, fluorosilicate glass, fumed silica, initiators, stabilisers, pigments

Storage:

Storage at 4°C-23°C. Do not use after the expiry date. Storage in a refrigerator is recommended.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Bifix Hybrid Abutment** must be reported to VOCO GmbH and the responsible authority.

Note:

The Summary of Safety and Clinical Performance of **Bifix Hybrid Abutment** can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at www.voco.dental.

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VC 60 AA210 E1 0724 99 © by VOCO