VOCO Futurabond® DC

EN Instructions for use
MD EU Medical device

Product description:

Futurabond DC is a dual-curing, self-etching bond reinforced with nano-particles for creating a durable bond between tooth substance and light-, self- and dual-curing restoratives without marginal leakage. This all-in-one adhesive consists of 2 components, which eliminates the need for a separate etching step yet still achieves the same adhesive properties as total-etch bonds.

Futurabond DC tolerates residual moisture and can be used universally with all light-curing, self-curing and dual-curing composites, compomers and ORMOCER[®]s.

Indications:

- Direct self-curing or dual-curing composite restorations and core build-ups
- Direct light-curing composite/compomer/ORMOCER[®]-based restorations
- Fixing of root pins/posts with dual-curing or self-curing composite cements
- Indirect restorations; when using dual-curing/self-curing composite cements for bonding of inlays, onlays, crowns and bridges

Contraindications:

Futurabond DC contains organic acids, methacrylates and BHT. **Futurabond DC** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

Futurabond DC is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

Futurabond DC should only be used by a professionally trained dental practitioner.

Use:

1. Creating a dry working field

A dry working field is essential due to the subsequent application of the restorative. Avoid contamination with blood or saliva. Use of a rubber dam is recommended.

2. Preparation

Prepare the cavity according to the principles of the adhesive restoration technique. Rinse cavity thoroughly with water. Remove excess moisture with a faint air jet. Do not overdry dentine. The surface of the cavity should remain slightly moist. If the tooth has not been prepared (e.g. cervical defect) the surface must be thoroughly cleaned and the enamel margin roughened with a finishing diamond. Areas in close proximity of the pulp should be protected with a suitable calcium hydroxide-based liner.

3. Mixing

Bring material to room temperature before use.

Futurabond DC is dual-curing and therefore must be applied immediately after mixing.

When using the bottled version:

Mix 1 drop of Liquid 1 and 1 drop of Liquid 2 on a mixing palette with a **Single Tim** for approximately 2 s (produces a self-etching adhesive).

When using SingleDose:

Detach a **SingleDose** at the perforation and turn the printed side upwards. To activate, hold the **SingleDose** blister between thumb and forefinger and press on the part of the blister marked "press here" in blue. This injects Liquid 1 into Liquid 2 in the dispensing bay in a 1:1 ratio. Pierce the foil of the dispensing bay with the **Single Tim**, enlarge the hole with a circular movement and stir well to create a homogeneous mixture.

4. Futurabond DC application

Apply the adhesive in a layer of medium thickness to the enamel/dentine and rub into the tooth surface for 20 s. Please note: The mixed bond is light-curing, thus should not be exposed to intensive ambient light. The operating light should be dimmed during application.

5. Drying

Dry the adhesive layer for at least 5 s with an air syringe.

6. Polymerisation

6.1. Direct restoration

Light-curing, self-curing and dual-curing materials

The bonding layer must be polymerised with blue light (halogen or LED lights) for 10 s when using light-curing, self-curing and dual-curing restoratives.

6.2. Indirect restoration

Light-curing luting composites

The bonding layer must be polymerised with blue light (halogen or LED lights) for 10 s when using light-curing luting materials.

Self-curing or dual-curing luting materials

When using dual-curing or self-curing luting materials, **Futurabond DC** also provides excellent adhesion without light-curing. **Futurabond DC** can also be light-cured for 10 s with blue light (halogen or LED lights) as an optional step.

Luting of root pins

1. Massage the mixed **Futurabond DC** (liquid 1 and 2) into the occlusal surface of the core and into the root canal for 20 s with an applicator. Disperse solvent with an oil-free air-jet for 5 s. Remove excess material by dabbing it off with an absorbent strip of paper. Do not light-cure the bonding layer.
2. Having prepared the root post in accordance with the manufacturer's instructions, insert it into the root canal using a suitable luting material such as **Rebilda DC** or **Bifix QM** (see respective instructions for use). The aim is to obtain a slight excess of cement.
3. Light-cure for at least 40 s to fix the root post and cure the **Futurabond DC** applied to the occlusal surface.

7. Advice on luting Maryland Bridges

For luting Maryland bridges, the enamel must either be roughened or conditioned with a phosphoric acid etching gel before applying **Futurabond DC**.

Self-etching bonding materials are less effective on unprepared enamel. Overapplication of bonding material to enamel that has not been ground or bevelled can lead to margin discoloration. For optimal adhesion results with **Futurabond DC**, the enamel should be prepared in advance in class III, IV, V cavities and for diastema closures.

Warnings, precautionary measures:

- Do not fold blister before activation. Blisters with obvious damage such as creases or kinks should not be used.
- Do not damage either the blister foil or tab on either the top or the bottom side with fingernails.

- Do not use any type of instrument for activation. Only the thumb and index finger should be used for activation with the thumb placed exactly over the mark on the label and the index finger placed exactly on the opposite side of the package.
- Activation and subsequent puncturing of the foil with a **Single Tim** should not take place within close proximity to the eyes.
- Never simultaneously press the reservoir from both sides or press the liquid back and forth between the reservoirs.
- Phenolic substances, especially preparations containing eugenol or thymol, inhibit the curing of **Futurabond DC**. Therefore, do not use zinc oxide eugenol cements or other eugenol-containing substances in combination with **Futurabond DC**.
- Avoid contact with the oral mucosa. Contact of the oral mucosa with **Futurabond DC** may result in a whitish discolouration caused by the coagulation of proteins. This is a reversible occurrence that subsides within a few days.
- In case of eye contact, rinse well with plenty of water and consult an ophthalmologist.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):

Methacrylate phosphoric acid ester, ethanol, water, BisGMA, HEMA, TEGDMA, TMPTMA, initiators, stabilisers

Storage:

Tightly close **Futurabond DC** bottles after use. Store in an upright position to guarantee flow-back of the liquid. Storage at 4°C-23°C. Do not expose **Futurabond DC** to direct sunlight or the chairlight. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Futurabond DC** must be reported to VOCO GmbH and the responsible authority.

Note:

The Summary of Safety and Clinical Performance of **Futurabond DC** can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at www.voco.dental.

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VC 60 AA1163 E1 0226 99 © by VOCO