EN	DE	FR	ES	PT	IT
EL	NL	DA	FI	NO	SV

VOCO
Grandio®SO

EN Instructions for use
MD EU Medical device
In accordance with DIN EN ISO 4049

Product description:

GrandioSO is a light-curing, radiopaque nano-hybrid restorative material. **GrandioSO** contains 89% by weight inorganic fillers and is applied with a dentine enamel bond. It is available in screw syringes and caps for direct application.

Shades:

A1, A2, A3, A3.5, A4, B1, B2, B3, C2, D3, Incisal, OA1, OA2, OA3.5, BL, GA3.25, GA5

Indications:

- Class I to V restorations
- Reconstruction of traumatically damaged anteriors
- Facetting of discoloured anteriors
- Correction of shape and shade for improved aesthetic appearance
- Locking, splinting of loose teeth
- Repairing veneers
- Restoration of deciduous teeth
- Core build-up under crowns
- Composite inlays

Contraindications:

GrandioSO contains methacrylate and BHT. **GrandioSO** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

GrandioSO is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

GrandioSO should only be used by a professionally trained dental practitioner.

Use:

Preparation:

Clean the teeth being treated. If applicable, mark occlusal contact points. Bring material to room temperature before use.

Shade selection:

Clean the teeth prior to shade selection. The shade is selected in comparison with the tooth while it is still moist.

Cavity preparation:

In principle, the cavities should be prepared according to the rules for adhesive restorative treatment, preparation should be minimally invasive in order to protect healthy dental hard tissue. For anterior restorations, angle off the enamel margins. Then clean and dry the cavity. Caries-free lesions in the cervical area do not need to be prepared; thorough cleaning is sufficient in this case.

Ensure that the working field is sufficiently dry. The use of a rubber dam is recommended. It is recommended that shaping aids be used in cases where a portion of the cavity is in the proximal region. It helps to use translucent female parts, wedged in the proximal region. Keeping separation to the minimum makes it easier to shape the proximal contacts and attach a female part. Suitable protection for the pulp should be applied close to the pulp and covered with a layer of stable cement if necessary.

Bonding materials:

GrandioSO is applied with a dentine enamel bond, using the adhesive technique. Any light-curing bonding materials can be used. Observe the relevant instructions for use during preparation (etching technique) and processing.

Application:

Apply the selected shade of **GrandioSO** (see shade selection) in layers of no more than 2 mm thickness; adapt using a suitable instrument and then cure.

Light-curing:

Conventional polymerisation devices are suitable for light-curing the material. Curing time if an LED/halogen lamp with a light output of at least 500 mW/cm² is used is 20 s for the shades A1, A2, A3, A3.5, A4, B1, B2, B3, C2, D3, Incisal, BL, GA3.25 40 s for the shades OA1, OA2, OA3.5, GA5

Curing time if an LED/halogen lamp with a light output of at least 800 mW/cm² is used is

- 10 s for the shades A1, A2, A3, B1, Incisal, BL
- 20 s for the shades A3.5, A4, B2, B3, C2, D3, GA3.25, GA5
- 40 s for the shades OA1, OA2, OA3.5

Bring the light aperture as close as possible to the light-curing unit as otherwise the curing will be less effective. Inadequate curing can cause discolouration and complaints.

Finishing:

Restorations can be finished and polished immediately after removing the female parts (e.g. fine and extra-fine diamond burs, abrasive discs). The tooth should then be fluoridized.

GrandioSO can be used for direct or indirect inlays using the known procedures. The physical stabilisation values are improved by the usual external (additional) curing.

Warnings, precautionary measures:

- Phenolic substances, especially products containing eugenol and thymol, lead to curing failures in restorative composites. Consequently, the use of zinc oxide/eugenol cements or other materials containing eugenol should be avoided in combination with restorative composites.
- Each **GrandioSO** Caps unit should be used for one patient only and should be applied exerting uniform but not excessive pressure.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Components (in descending order according to content):

Barium aluminium borosilicate glass, silicon dioxide, BisGMA, TEGDMA, BisEMA, initiators, stabilisers, pigments

Storage instructions and application method:

Storage at 4°C-23°C. Seal the syringes immediately after removing the material in order to avoid exposure to light and subsequent polymerisation. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **GrandioSO** must be reported to VOCO GmbH and the responsible authority.

Note:

The Summary of Safety and Clinical Performance of **GrandioSO** can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at www.voco.dental.

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VC 60 AA2600 E1 0123 99 © by VOCO