

## Direct, complex reinforcement of endodontically treated, pericervical dentine with supercomposite for molar preservation

Author: Dr Clarence Tam

A 65-year-old Canadian dentist presented in my practice with fourth-degree root caries in tooth 26 (upper left first molar). A distal-occlusal restoration was placed, also extending into the root region. The pulp was not opened during the excavation. Ten weeks later, after a sudden progression of the symptoms, the decision was taken to perform a root canal treatment, and the patient was referred to an endodontist.

Once the endodontic treatment and temporary restoration were complete, the patient presented in the practice again for further treatment (Fig. 1).

The patient requested a restoration made of a composite with optimal mechanical properties so that the coronal structure is stabilised, the cusps covered and thus reinforced, and as much of the pericervical dentine as possible is preserved, as this factor ultimately determines the fracture resistance and prognosis of a tooth.

In order to achieve the layer thickness required for this type of restoration, an occlusal reduction of 2 mm is required to ensure the longevity of the composite restoration. This was done by firstly producing grooves as depth markings (Fig. 2), then removing the temporary restoration and the old amalgam restorations and finally reducing the occlusal surface (Fig. 3+4).

To seal the pulp chamber, a layer of glass ionomer restorative material was applied as a fluoride-releasing, protective liner.

Then just a shoulder measuring 1 mm in depth and 1.5 mm across was prepared all around the tooth (Fig. 5).

This method allowed preservation of all the remaining pericervical dentine, with the result that the tooth still possesses high fracture resistance even after endodontic treatment.

Following micro air abrasion with 27 µm aluminium oxide, the tooth was etched for 15 seconds with 33% orthophosphoric acid so as to create an even larger bonding surface. Futurabond U (VOCO) was then applied and massaged in for 20 seconds so as to achieve maximum wetting and consequently adhesion of the filling.

### Technology

The buccal and lingual walls of the cusps were sculpted using the nanohybrid composite GrandioSO (VOCO) in A3 shade and placed by hand first of all, with the height of the neighbouring cusps being used as a reference (Fig. 6). Figure 7 shows the occlusal view of the initial buccal and lingual cusps and the preliminary formation of the ridge between the distobuccal and mesiopalatal cusps, which will later be angled. The marginal ridges are formed with the Triodent/Dentsply V4 sectional matrix system only when the buccal and lingual walls of the cusps are finished. A class I restoration is produced (Fig. 8). The mesiobuccal cusp is finished and the fissures are given a contrasting shade (Fig. 9).

The restoration is completed with successive cusp formation on the central mesial, mesiopalatal and distopalatal sides (Fig. 10). Finally, occlusal adaptations are performed (Fig. 11) before the restoration with its five cusps is polished and finished (Fig. 12).

## Material discussion

### Selection of bonding agent: Futurabond U (VOCO).

This eighth-generation bonding agent is characterised by outstanding shear strength and is suitable for a range of applications: total etch, selective etch and self-etch. It has dual-curing properties, meaning that no separate activator needs to be added for dark curing (self-cure). Even after thermocycling, the bond to enamel is significantly better than with comparable products and achieves a remarkable 33.8 MPa in the thermocyclic self-etch technique. The bond to dentine is also notable, and even in the thermocyclic self-etch technique is a full 8 MPa higher than with a conventional alternative.

**Selection of composite: GrandioSO (VOCO), shade A3.** GrandioSO was selected for its excellent flexural strength of 187 MPa, which is considerably higher than that of dentine (165.6 MPa). Even after thermocycling, it is very biomimetic, at 158 MPa. As the material is employed in a region subject to very high loads, its compressive strength is decisive for avoiding cohesive compression shear fractures.

At 439 MPa, it is far higher than that of enamel (384 MPa) and also dentine (297 MPa). The microhardness of the surface and abrasion resistance are two further factors which ensure that the surface remains intact for as long as possible. The surface hardness of GrandioSO is almost twice as high as that of the other composites tested (*Behrend, D., University of Rostock, 2010*), and at 211 MHV is very close to natural enamel, which has a surface hardness of 350-450 MHV. It proved possible to use just a single shade (A3) for the entire restoration, as its properties and the high filler content permit an outstanding chameleon effect.

## Conclusion

The end result is an extraordinarily stable, direct acrylic restoration, which protects the underlying, damaged tooth structure and, thanks to the strong bond and its physical properties, possesses considerable durability. In addition, it enables great aesthetics.

## IMAGES



Fig. 1: Preoperative situation following completed endodontic treatment



Fig. 2: Production of the grooves



Fig. 3: Removal of the old restorations and cusp reduction



Fig. 4: Application of a layer of glass ionomer restorative material as a fluoride-releasing cavity liner



Fig. 5: Creation of a broad shoulder as a measure for the retention shape



Fig. 6: Sculpting of the buccal and lingual walls of the cusps



Fig. 7: Occlusal view



Fig. 8: Marginal ridge formation after finishing of the buccal and lingual walls of the cusps



Fig. 9: Production of the mesiobuccal cusp and contrasting shade for the fissure system



Fig. 10: Finishing of the restoration



Fig. 11: Occlusal adaptations



Fig. 12: Polished restoration

### About the author



Dr Clarence Tam heads a practice in Auckland, New Zealand, which specialises in cosmetic and restorative dentistry. She is originally from Canada, where she completed her Doctor of Dental Surgery and General Practice Residency at the University of Western Ontario and the University of Toronto, respectively. Clarence is the Chairperson and Director of the New Zealand Academy of Cosmetic Dentistry. She is currently the only person in Australasia to hold Board-Certified Accredited Member Status with the American Academy of Cosmetic Dentistry. She frequently and continually lectures internationally.

### Contact details

Dr Clarence Tam, HBSc, DDS  
Cosmetic and General Dentistry  
Auckland (New Zealand)  
E-mail: [clarence.tam@gmail.com](mailto:clarence.tam@gmail.com)  
[www.clarencetam.co.nz](http://www.clarencetam.co.nz)