**Material-related characteristics**

Chemical composition [wt. %]

Polymethylmethacrylate (PMMA)	> 99
Coloring pigments	< 1

**Physical specifications**

Flexural strength	[20795-1]	[MPa]	65
Flexural modulus	[20795-1]	[MPa]	2030
Fracture toughness	[20795-1]	[MPam <sup>1/2</sup> ]	2,6
Fracture work	[ISO 20795-1]	[J/m <sup>2</sup> ]	> 2200
Charpy impact strength, +23°C	[ISO 179/1eU]	[kJ/m <sup>2</sup> ]	45
Density	[ISO 1183]	[g/cm <sup>3</sup> ]	1.17
Water sorption		[µg/mm <sup>2</sup> ]	< 24
Solubility	[ISO 20795-1]	[µg/mm <sup>2</sup> ]	< 0,3
Residual monomer content	[ISO 20795-1]	[%]	< 0,5

**1. Intended Purpose**

**CediTEC DB** polymer milling blanks are intended for the fabrication of restorations.

**2. Indications**

**CediTEC DB** is a pre-colored dental milling blank made of impact resistant PMMA for the manufacture of denture bases for removable dentures. **CediTEC DB** is suitable for long-term use in the oral cavity for up to 10 years.

**3. Contraindications**

- Intolerance to the contained components
- Crowns, bridges, model casting with clasps, bars, narrow sublingual bars and transversal connectors

**4. Safety information**

Avoid inhalation of milling dust during processing. In case of Inhalation: Move subject to fresh air. Give oxygen or artificial respiration as required.

Avoid all contact with fluids. Wear gloves, protective goggles and a face mask to avoid skin irritation. In case of eye contact: Flush eyes promptly with copious amounts of water for 15 minutes and consult a physician. Wash skin with soap and water.

In case of ingestion: Contact your regional Poison Control Center immediately.

**5. Handling and storage**

Only use clean, defect-free material from undamaged packaging and make sure that the blanks are stored at room temperature, in their original carton, dry and protected from light. Avoid contact with fluids.

**6. Construction**

This medical device may only be processed by trained personnel. The denture base should be positioned centrally in the milling blank and should not be less than a minimum wall thickness of 1 mm.

**7. Milling and further processing**

Milling must be carried out with a milling strategy adapted to PMMA or milling parameters for **CediTEC DB**.

Only use one-edged milling tools in the milling machine or wet grinding procedures! **Further processing:**

Carbide cutters with FSQ toothing and black-ring cutters for titanium are suitable for roughing. The separation of the denture bases from the milling blank and the fine finishing can be performed with simple, cross-toothed carbide millers. The connectors are cut out carefully and without high pressure. The polishing is carried out in the usual manner with polishing agents customary for denture base materials. When polishing the material, ensure that the contact pressure is kept as low as possible, otherwise the heat generated can lead to smearing / melting.

**8. Preparation of bonding surfaces**

Check the fit of the prosthetic teeth before bonding. In case of interfering contacts, make the correction on the tooth. For an optimum bond, roughen the bonding surfaces of the restoration with aluminum oxide (1 - 2 bar / 50 to 125 µm). Carefully remove blasting agent residues mechanically or by means of a water-ultrasonic bath (free of solvents). The restoration must then be dried thoroughly. A suitable luting material (e.g. **CediTEC, VOCO**) must be used to secure denture teeth in the workpiece. Details can be found in the corresponding instructions for use.

**9. Individualisation, relining and repair**

Individualisations, relinings and repairs of **CediTEC DB** can be carried out with conventional auto-polymerizing denture base resin / light-curing composites. Roughen the restoration surface by blasting (Al<sub>2</sub>O<sub>3</sub> / 1 - 2 bar / 50 to 125 µm) or with emery paper. Carefully remove dust residues mechanically or by means of a water-ultrasonic bath (free of solvents). Then dry the restoration with air. Apply a suitable adhesive system. Details can be found in the corresponding instructions for use.

**10. Cleaning**

- Ultrasonic bath (max. 40 °C / 104 °F) with mild cleaning agent

Do not use the following methods for cleaning:

- Steam blasting
- Autoclave method (e.g. sterilization / disinfection)
- Organic solvents and acids (e.g. acetone, hydrochloric acid)
- Highly concentrated alcohols
- Chemical denture cleaners
- Abrasive cleaning agents

**11. Disposal**

Remains of PMMA blanks can be disposed of in the regular household garbage. **Note:** Any serious incident that has occurred in relation to the devices should be reported to the manufacturer and the competent authority to the Member State in which the user and/or patient is established.

**Explanation of symbols**

Manufacturer	Date of manufacture	Use-by date
Batch Code	Catalog number	Medical Device
Keep dry	Keep away from sunlight	Consult instructions for use

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VC 60 AA6195 E1 1023 99 © by VOCO
**VOCO GmbH**