

Quick Up®

SELF-CURING LUTING MATERIAL FOR ATTACHMENTS
AND SECONDARY COMPONENTS IN DENTURES



Quick Up®

ONE SET FOR ALL LUTING NEEDS

In many cases the hold of full dentures is unsatisfactory. The main reasons for this are bone atrophy, resorption processes and changes to the denture base, which is made up of both hard and soft tissue.

Implants are very frequently used for fixing full dentures. Small implant diameters and simplified less time-consuming surgical procedures, make this improvement in the hold of the denture very attractive both to patients and dentists. In most cases the existing denture can continue to be worn. This saves patients time, inconvenience and expense. The dentist is able to concentrate on the implantation and can make a removable connection with the implants simply by fitting the corresponding attachments or secondary components (female parts) to the existing denture.

The Quick Up Set contains all the components required for attaching female parts, attachments and other attachment elements to composite-based dentures. The reattachment of secondary components in dentures is therefore a further indication for Quick Up.

In developing Quick Up we attached great importance to making the product as easy to use as possible, so that the bonding in of secondary components would be a stress-free procedure for the dentist and their patient delivering predictable results.

The steps for bonding in secondary components always follow the same procedure: the denture base around the implants must be relief-ground immediately after the placement of the implants in order to facilitate irritation free healing. In the next step the secondary components (female parts, attachments) must be attached to the denture base after the implants have healed or immediately after placement in the case of immediate loading.

This procedure must be carried out on the patient to ensure a more accurate fit and to enable a check to be made before the final intraoral bonding in of the denture. Any corrections necessary can then still be made without any difficulty. Moreover, the patient is spared the inconvenience of having to manage without a denture. This procedure saves the dentist a lot of time.

Clinical procedure in three steps



The marked areas are hollowed out to a sufficiently large extent as a preparation for bonding in the secondary components

into the denture base. This is done so that the denture can be inserted without any difficulty over the implants / root post caps with the secondary components attached. The dentist checks that the secondary components are securely seated into the denture base as well as carrying out the clinical check with Fit Test C & B, a special control and covering silicone.

The purpose of these checks is to enable the dentist to make any corrections required without being pressed for time, before using Quick Up and the final bonding in stage.

The additional application steps comprise merely removing Fit Test C & B from the denture base and if necessary, blocking out or covering the implants or primary components / male parts and the gingiva with Fit Test C & B in the patient's mouth.

The special adhesive is then applied to the denture base. The previously roughened secondary components / attachments are fitted in the patient's mouth and the holes in the denture base provided for the attachments are only filled approx. 2 / 3 full with Quick Up.

The denture is inserted and brought into occlusion under the supervision of the dentist.

Quick Up cures in just 3.5 minutes from the start of mixing, the denture can then be removed from the patients mouth.

The attachments or secondary components / female parts have now already been bonded into the denture base. Ditchings are now quickly and systematically filled with the light-curing Quick Up LC. This procedure comprising the application of both self-curing Quick Up Material followed by correction with light-curing Quick Up LC is the best way of preventing a disastrous gumming-up of the primary and secondary components in the mouth. Finishing, polishing and the removal of the control silicone Fit Test C & B used for blocking out conclude the clinical procedure.

Quick Up®

LUTING OF ATTACHMENTS IN JUST THREE STEPS

First step:

- Mark the position of the secondary components on the denture base
- Determine the space required and check the secondary components



Quick Up®

QUICK UP® ADHESIVE PROVIDES A SECURE HOLD IN THE DENTURE

Second step:

- Wet prepared holes in the denture base with adhesive
- Place attachments / secondary components on the primary components in the mouth
- Cover areas to be protected such as the gingiva and joints required with the enclosed control and covering silicone Fit Test C & B
- Bond in the prepared secondary components in the denture base with Quick Up




Quick Up®

CORRECTIONS EASILY MADE WITH QUICK UP® LC

Third step:

- Fill ditchings with Quick Up LC, if necessary
- Finishing and polishing



Quick Up® – the better alternative

Quick Up enables attachments and secondary components to be bonded into the patient's denture quickly and easily. Alternative systems that use soft relining materials instead of (metal-)attachments for anchoring don't reach a comparable stability and suffer from the fact that they have to be revised more frequently. Difficulties in the mechanical or adhesive retention, as well as a generally shorter life-time of soft materials are reasons for this.

In contrast, the definite luting of (metal-)attachments with Quick Up offers a definite solution with higher stability and less need for control and revisions.

What users are saying

Dentists who already use Quick Up for Locator® systems and attachments were very positive about the composition of the set and the individual components. This is what they had to say:

- very good handling thanks to the small application syringe (especially intraorally)
- shade almost identical with acrylic resin
- low material consumption
- very simply designed product

A dentist and lecturer from the USA wrote in similar vein: "Many thanks for the opportunity of using Quick Up. My experiences with this material have so far been extremely positive. I am impressed by the adhesive and the form of application, as well as by the light-curing correction material which is ideal for filling small ditchings."

Quick Up®

SIMPLE, FAST, STRESS-FREE

Indications

Gingiva-coloured composite for luting attachments, such as ball, Locator and telescope retention elements in prostheses with acrylic bases

Reluting secondary elements in prostheses with acrylic bases, e.g. connecting bars

Advantages

- Very simple chairside use
- High strength values
- Complete set for immediate use, including checking and block out material for primary structures that require protection, e.g. implants
- Ideal combination of self-curing luting material and light-curing correction material

Presentation

- REF 1625 Set QuickMix syringe 7.5 g, Quick Up LC syringe 2 g, bottle 4 ml adhesive, Fit Test C & B QuickMix syringe 5 ml, accessories
- REF 1626 QuickMix syringe 7.5 g, accessories
- REF 1627 Quick Up LC syringe 2 g, accessories
- REF 1628 Adhesive bottle 4 ml, accessories
- REF 2261 Dispenser – QuickMix
- REF 2189 Mixing tips type 10, 50 pcs.
- REF 2144 Application cannulae type 44, 100 pcs.
- REF 2247 Single Tim, application brushes, 100 pcs.
- REF 2315 Mixing palettes, 20 pcs.



VOCO GmbH
 Anton-Flettner-Straße 1-3
 27472 Cuxhaven
 Germany

Freecall: 00 800 44 444 555
 Fax: +49 (0) 4721-719-140

info@voco.com
 www.voco.dental

Available from: