# Admira® Fusion

## EN Instructions for use

MD EU Medical device

In accordance with DIN EN ISO 4049

### Product description:

Admira Fusion is a light-curing, radiopaque nano-hybrid ORMOCEP® restorative material. Admira Fusion contains 84% by weight inorganic fillers and is applied with a dentine enamel bond. It is available in screw syringes and caps for direct application.

### Shades:

A1, A2, A3, A3.5, A4, B1, B2, B3, C2, D3, OA1, OA2, OA3, OA3.5, BL, BXL, Incisal, GA3.25, GA5

### Indications:

- Class I to V restorations
- Base in class I and II cavities
- Reconstruction of traumatically damaged anterior teeth
- Facetting of discoloured anterior teeth
- Correction of shape and shade for improved aesthetic appearance
- Locking, splinting of loose teeth
- Repairing veneers, small enamel defects and temporary C&B-materials
- Extended fissure sealing
- Restoration of deciduous teeth
- Core build-up
- Composite inlays

### Contraindications:

Admira Fusion contains methacrylate and BHT. Admira Fusion should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

### Patient target group:

Admira Fusion is suitable for use in all patients without any age or gender restrictions.

### Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

### User:

Admira Fusion should only be used by a professionally trained dental practitioner.

### Preparation:

Clean the teeth being treated. If applicable, mark occlusal contact points. Bring material to room temperature before use.

### Shade selection:

To select the shade, compare the cleaned tooth with the enclosed shade guide prior to anaesthesia. To do so, moisten the Admira Fusion sample shades and the tooth and determine the shade, preferably in daylight conditions.

### Cavity preparation:

In principle, the cavities should be prepared according to the rules for adhesive restorative treatment; preparation should be minimally invasive in order to protect healthy dental hard tissue. For anterior restorations, angle off the enamel margins. Then clean and dry the cavity. Caries-free lesions in the cervical area do not need to be prepared; thorough cleaning is sufficient in this case.

Ensure that the working field is sufficiently dry. The use of a rubber dam is recommended. It is recommended that shaping aids be used in cases where a portion of the cavity is in the proximal region. It helps to use translucent matrices, wedged in the proximal region. Keeping separation to the minimum makes it easier to shape the proximal contacts and attach a matrix. Suitable protection for the pulp should be applied close to the pulp and covered with a layer of stable cement if necessary.

### Bonding materials:

Admira Fusion is applied with a dentine enamel bond, using the adhesive technique. Any light-curing bonding materials can be used. Observe the relevant instructions for use during preparation (etching technique) and processing.

### Application:

Apply the selected shade of Admira Fusion (see shade selection) in layers of no more than 2 mm thickness; adapt using a suitable instrument and then cure.

### Light-curing:

Conventional polymerisation devices are suitable for light-curing the material. Curing time if an LED/halogen lamp with a light output of at least 500 mW/cm² is used is 20 s for the shades A1, A2, A3, A3.5, A4, B1, B2, B3, C2, D3, BL, BXL, Incisal, GA3.25

40 s for the shades OA1, OA2, OA3, OA3.5, GA5

Bring the light aperture as close as possible to the light-curing unit as otherwise the curing will be less effective. Inadequate curing can cause discolouration and complaints.

### Finishing:

Restorations can be finished and polished immediately after removing the matrices (e.g. fine and extra-fine diamond burs, abrasive discs). The tooth should then be fluoridized.

Admira Fusion can be used for direct or indirect inlays using the known procedures. The physical stabilisation values are improved by the usual external (additional) curing.

### Warnings, precautionary measures:

- Phenolic substances, especially products containing eugenol and thymol, lead to curing failures in restorative composites. Consequently, the use of zinc oxide/eugenol cements or other materials containing eugenol should be avoided in combination with restorative composites.

- Each Admira Fusion Caps unit should be used for one patient only and should be applied exerting uniform but not excessive pressure.

- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

### Constituents (in descending order according to content):

Barium aluminium borosilicate glass, ORMOCEP® resin, silicon dioxide, initiators, stabilisers, pigments

### Storage instructions and application method:

Storage at 4°C-23°C. Seal the syringes immediately after removing the material in order to avoid exposure to light and subsequent polymerisation. Do not use after the expiry date.

### Disposal:

Dispose of the product in accordance with local regulations.

### Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of Admira Fusion must be reported to VOCO GmbH and the responsible authority.

### Note:

The Summary of Safety and Clinical Performance of Admira Fusion can be found in the European database on medical devices (EUDAMED - <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at [www.voco.dental](http://www.voco.dental).

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ORMOCER® is a registered trademark of Fraunhofer Gesellschaft zur Förderung der angewandten Forschung e.V., München, Germany.