VOCO

Perfect Bleach® Office+

EN Instructions for use
MD EU Medical device

CE 0482

In-office tooth whitening system

Product description:

Perfect Bleach Office+ is a gel for teeth whitening carried out by the dentist in his dental office. The 4:1-double syringe contains 5 g (4 ml) of gel. During dispensing from the self-activating syringe, base and activator are mixed to a gel with a hydrogen peroxide content of 35 %.

Indications:

Perfect Bleach Office+ is suited for whitening individual or multiple discoloured, vital or devitalised teeth.

The whitening success depends on the natural tooth shade as well as the intensity and cause of the discolouration. It is possible to whiten discolourations of organic origin, which generally have the following causes:

- discolouration as a result of trauma
- discolouration caused by medication (e.g. tetracycline)

Discolourations of inorganic origin, such as components of amalgam (silver), are not whitened, nor are restorations.

Contraindications:

- During pregnancy and lactation
- In patients with allergies or intolerances to one or several of the ingredients
- In patients who have periodontal disease
- In patients with exposed dentine, defective filling margins or untreated caries
- In hypersensitive patients

Application must be carried out by a dentist or under direct supervision of a dentist.

Discontinue treatment if the patient experiences excessive pain!

Method of Use:

Preparation for treatment:

- Discuss the cause(s) of the tooth discolouration with the patient (especially external causes).
- Professionally clean to remove calculus, surface contaminations and anomalies. Determine the type and intensity of the discolouration and evaluate if whitening treatment is indicated.
- Inform the patient in detail about the options, risks and side effects. Also discuss potentially required subsequent treatments, such as replacement of fillings or other tooth-coloured restorations, since these will not be whitened during the whitening treatment.
- Insufficient fillings and restorations as well as untreated carious lesions must be replaced or receive temporary treatment before whitening.
- Do not use anaesthesia!
- The patient, practitioner and assistants must all wear protective eyewear during the application of **Perfect Bleach Office+**. The practitioner and dental assistants must also wear gloves and a mask. Do not allow **Perfect Bleach Office+** to come in contact with skin, eyes or mucosa.

A) External whitening:

1. Insert cheek and lip retractor or similar.
2. Apply **LC Dam** to protect the gingiva (see the instructions for use for **LC Dam**). A rubber dam can additionally be used, which should be inserted before the application of **LC Dam**.
3. Remove safety clip, carefully remove syringe cap from the syringe and extrude material until it is delivered uniformly from both chambers. Then attach the pointed mixing cannula and extract the gel from the tip until an evenly coloured mass with a gel-like consistency emerges from the cannula. Discard the initial inhomogeneous mass.
4. Apply **Perfect Bleach Office+** directly from the mixing tip onto the labial surfaces of the teeth in a 1 to 2 mm thick layer. The material applied to the teeth can be evenly distributed and moved around every 5 min with an instrument to improve effectiveness.
5. Leave whitening gel to act on the teeth for approximately 10 to 15 min. The gingiva should be checked for irritation during the entire duration of the treatment.
6. Thoroughly aspirate the whitening gel from the teeth with a surgical aspirator.
7. Steps 4 through 6 can be repeated twice if the teeth have not reached the desired shade. A fresh mixing cannula should be used for each new application and the gingiva protection should be checked.
8. When the teeth have reached the desired shade, remove any gel remaining after aspiration of **Perfect Bleach Office+** with gauze or a cotton pellet. Rinse with plenty of water afterwards.
9. The attached mixing tip is intended for single use only. Detach it after usage and seal the syringe again with the original lid (cap).

B) Internal whitening of endodontically treated teeth:

1. Before the whitening treatment, an x-ray has to be taken to verify that the root canal filling is closely sealed.
2. First, expose the entrance to the pulp chamber and root canal and carefully remove the existing restorative. Ensure that sufficient space is created in the area of the coronal pulp to accommodate the whitening gel.
3. The root canal filling should be removed to a maximum depth of ca. 2 mm below the cervical gingiva and should afterwards be tightly covered with a lining material.
4. Create a dry working field around the affected tooth. Using a rubber dam is recommended to provide adequate protection against inhalation or swallowing of the highly reactive whitening gel. For additional protection of the gingiva apply **LC Dam** according to instructions for use.
5. Then, using the blunt mixing cannula with the intraoral tip, fill the cavity with **Perfect Bleach Office+** (see step A3), which should be left in the cavity for approximately 10 to 15 minutes.
6. During the reaction time of the preparation, the cavity can be sealed with a suitable, secure temporary obturation material and the patient can for example take a seat in the waiting room. The patient must be informed that he or she may neither eat nor drink during the reaction time of the whitening gel, since this could result in undesired leaking of the whitening gel through loss of the temporary.
7. After the reaction time, carefully remove the whitening gel from the cavity using a surgical aspirator.
8. To increase the whitening effect, steps 5 through 7 can be repeated twice.
9. When the teeth have reached the desired shade, the cavity has to be rinsed thoroughly with plenty of water after removing the whitening gel.

10. Seal the cavity with a suitable provisional restorative (e.g. a glass ionomer cement). This filling should remain in situ for a minimum of 14 days before the definitive treatment with a composite can be carried out.

After the treatment:

- The most common undesired effects are temporary hypersensitivity of the treated teeth and periodontal irritation. These symptoms disappear within a short period of time after ending or interrupting the treatment. Carry out remineralisation of the teeth with a fluoride preparation, if desired.
- The teeth may slightly darken again during the first few weeks after whitening. Tooth-shaded restorations should therefore be carried out two weeks after whitening at the earliest.
- Strongly staining or acidic food and drink as well as tobacco should be avoided within the first days after the application of **Perfect Bleach Office+**, since increased staining may result otherwise.

Information, precautionary measures:

Perfect Bleach Office+ contains hydrogen peroxide, potassium nitrate, sodium fluoride, sodium hydroxide, glycol derivatives, thickening agent and colourant. **Perfect Bleach Office+** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.

The pure activator (red component) is a corrosive agent. The mixed whitening gel acts as an irritant when in contact with skin, eyes or oral mucosa. If contact occurs, immediately rinse the affected area with plenty of water. Also seek medical attention in case of contact with the eyes. Do not swallow (causes nausea and heartburn). Keep out of reach of children!

The warranted product characteristics can only be achieved when using the supplied original mixing tips, which are also available as accessories.

Storage:

Refrigerate **Perfect Bleach Office+** at a temperature between 2 °C and 8 °C. Do not store in the freezer. Protect from direct sunlight and heat. Extended storage at higher temperatures may occasionally cause bubbles to form in the gel. This may lead to inhomogeneous, more liquid mixtures, especially in the last quarter of the double syringe. Do not use the material in this case. Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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VC 60 CC1672 E1 1123 99 © by VOCO
VOCO GmbH