VOCO


VOCO Structur® CAD

Minimum wall thicknesses

	Single crown	Abutment crown	Anterior bridge	Posterior bridge
Cervical wall thickness	0.6 - 0.8 mm	—	0.8 mm	1.0 mm
Occlusal wall thickness	1.2 - 1.5 mm	1.5 mm	1.2 - 1.5 mm	1.2 - 1.5 mm
Circular wall thickness	—	0.8 mm	—	—
Shoulder with Ti base	—	0.4 mm	—	—
Connector cross-section 1 pontic	—	—	10 - 12 mm ²	12 - 15 mm ²
Connector cross-section 2 pontics	—	—	12 mm ²	16 - 20 mm ²

Abutment

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MD EU Medical device

Product description:

Structur CAD is a blank for the production of highly aesthetic temporary work and long-term temporaries. The maximum clinical period of wear is 3 years. **Structur CAD** contains 27% by weight inorganic fillers in a polymer matrix. Thanks to the composite technology, the material displays good abrasion resistance as well as high fracture resistance and is thus also ideal for long-term temporaries. **Structur CAD** is fluorescent. The restorations can be inserted with a temporary cement or an adhesive technique. **Structur CAD** is available as a block (15.5 x 19 x 39 mm) and a disc (Ø 98.4 x 20 mm).

Shades:

A1, A2, A3

Indications:

- Temporary bridges spanning up to two pontics
- Temporary abutment crowns
- Temporary crowns

Contraindications:

Structur CAD contains methacrylates. **Structur CAD** should not be used for patients with a known hypersensitivity (allergy) to these constituents.

Patient target group:

Structur CAD is suitable for application on all patients without any age or gender restrictions.

Product performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

Application

Structur CAD should only be applied by a professionally trained dental practitioner.

Shade selection:

Use the VITA® shade system to determine the shade against the cleaned but unprepared tooth prior to anaesthesia and preferably in daylight conditions.

Core and cavity preparation:

As a general rule, the core/cavity preparation should be performed in accordance with the conventional rules for all-ceramic restorations. That means rounding off interior corners and edges and using a shoulder preparation with rounded interior edges or a chamfer preparation.

Observe the minimum thicknesses given at the start of this IFU.

Note: Do not create any undercuts.

Minimum wall thicknesses (see Table):

Note: Pontic design : height ≥ width

Grinding/milling process:

Structur CAD is suitable for both dry and wet processing. Select the corresponding block/disc size for the designed restoration and the grinding/milling parameters for the **Structur CAD**. When doing so, pay attention to the software settings of the respective CAD/CAM systems. **"Diamond-coated tools" and composite parameters are recommended** for the CAM processing. If the grinding/milling parameters are not already available in the CAD/CAM systems' software settings, they will need to be added before you proceed. Please contact the CAD/CAM system provider for assistance.

The instructions for use from the manufacturers must be observed and followed. Check the finished restoration for defects such as cracks or chips. If the restoration displays any defects, it must be rejected.

Finishing/polishing:

Once the CAM process is complete, detach the restoration from the lugs using a carbide tipped finishing bur or suitable cutting wheel. Sand the corresponding areas with fine-toothed carbide tipped finishing burs or fine diamond-coated grinders. Pay attention to any contact points.

Note: No post-processing of the shoulder to the Ti base on abutment crowns. Polish restoration chairside with conventional composite polishers. A multiphase polishing system is advantageous.

A polishing paste can be used labside in combination with goat's hair brushes and cotton/leather buffers.

Use only light pressure and wiping motions to control the removal. The instructions for use from the manufacturers must be observed and followed.

Luting:

Preparing the restoration

For an optimal bond, roughen the luting surface of the restoration using aluminium oxide (50-100 µm) at 1-2 bar or roughen with a fine HM milling device. Use a suction device to remove the dust produced. Remove abrasive material residues carefully with an ultrasonic bath (70% ethanol) or steam cleaner. Then dry the restoration with air. Final cleaning with medical alcohol is possible. The instructions for use must be observed and followed.

Temporary luting for crowns and bridges

Lute **Structur CAD** with a temporary cement (e.g., **Provilcol QM**). **Note:** Use a temporary, eugenol-free cement if the permanent restoration is to be subsequently luted with an adhesive.

Adhesive luting for crowns and bridges

For a wearing period of > 30 days, the restoration can be luted using a composite-based adhesive luting material (e.g., **Bifix QM**). The instructions for use must be observed and followed.

- Note:**
- Adhesive luting causes time-consuming removal of the temporary restoration.
 - Regular check-ups and follow-up examinations are required for long-term use.

Extraoral luting of the abutment crown on the titanium base

- Observe the manufacturer's specifications for the preconditioning of the Ti base (e.g., sandblasting, cleaning).
 - Precondition the adhesive surface of the Ti base with a suitable bonding agent.
 - Prepare the adhesive surface of the crown with a suitable bonding agent (e.g., **Ceramic Bond**).
- When using Ceramic Bond:**
Apply **Ceramic Bond** to the adhesive surface with the **Single Tim** or a disposable brush, allow it to act for 60 seconds and dry carefully with air.
- Use a suitable opaque luting composite to stick the restoration and Ti base together. Apply the luting composite to the abutment crown and the Ti base alike.
- The instructions for use from the manufacturers must be observed and followed.

Customisation:

For a highly aesthetic temporary restoration, the restorations can be customised, characterised or repaired at any time using a composite/ORMOCER®. Roughen the restoration surface via grinding or sandblasting (Al₂O₃, 50-100 µm, 1-2 bar). Remove abrasive material residues/dust carefully with an ultrasonic bath (70% ethanol) or steam cleaner. Then dry the restoration with air. Apply a suitable adhesive system (e.g., **Futurabond U**) in accordance with the instructions for use. Using **GrandioSO, Flow** or **Heavy Flow**, for example, in combination with **FinalTouch**, you can customise the restorations quickly and simply with purely light-curing techniques.

The instructions for use from the manufacturers must be observed and followed.

Information, precautionary measures:

Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application.

Storage instructions and application method:

Store at between 4°C and 23°C. If the material is stored in the refrigerator, bring it to room temperature before use. Do not store the product in a place where it is exposed to direct artificial light or sunlight. Do not use the product after its expiry date.

Disposal:

Dispose of the product according to local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's condition and a serious risk to public health that arise or could have arisen in association with the use of **Structur CAD** must be reported to VOCO GmbH and the responsible authority.

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VC 60 BB6071 E1 0723 99 © by VOCO
ORMOCER® is a registered trademark of Fraunhofer Gesellschaft zur Förderung der angewandten Forschung e.V., München, Germany.
VITA® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co.