**Product description:**

**V-Print model fast** is a light-curing resin for the generative production of fast printable models, especially for the dental thermoforming technique (e.g. for aligner or retainer splints).

**Colour:**

Blue

**Indication:**

– Dental models

**Contraindications:**

**V-Print model fast** contains (meth)acrylate and phosphine oxide. **V-Print model fast** should not be used on patients with known hypersensitivities (allergies) to these ingredients.

**Performance features:**

The product's performance features satisfy the requirements of its intended use and the relevant product standards.

**Use / User:**

**V-Print model fast** should only be used by a professionally trained dental practitioner.

**Hardware and software requirements**

CAD software¹ dental scanner	Software for planning and designing dental models. The software and dental scanner must satisfy current medical device specifications and allow for issuance of the patient-specific design as an STL data set.
CAM-software	Software for preparation of the print order. The part will not be modified during this process. Structures that facilitate the 3D printing are simply created. For example: - Autodesk Netfabb version 2018 or later for SolFlex 3D printing.

¹The designation **Software as Medical Device SaMD** includes standalone (autonomous) software that is a medical device (MD) and not part of one.

Manufacturing equipment	For example: VOCO SolFlex 170 (FlexVat) VOCO SolFlex 350 (FlexVat & PowerVat) VOCO SolFlex 650 (PowerVat)
Post-curing devices	For example: Otoflash G171

See also: accompanying list of resources or [www.voco.dental/3dprintingpartners](http://www.voco.dental/3dprintingpartners)

All manuals and/or operating instructions for the respective programmes, and for device, materials and/or parts manufacturers, which are required for the manufacturing process, must be observed.

Clarify ahead of time whether the programmes, devices and/or objects that you intend to use have been designed and approved for the corresponding applications.

**CAUTION:** Non-authorized changes to the process equipment, parameters, or software could result in the **V-Print model fast** end object not satisfying specifications.

**Preparation:**

Information on printer approvals can be requested from the respective equipment manufacturers or from **VOCO** ([www.voco.dental/3dprintingpartners](http://www.voco.dental/3dprintingpartners)).

For an indication-appropriate CAD construction, the following design conditions must be observed:

**Recommendation:** For a material-saving printing process models that are not filled a wall thickness of 2 mm is recommended.

Prepare a print job using slicing software. In addition to the material-dependent construction specifications in these Instructions for use, please also observe the dependencies of the positioning, support type, and fit found in our other documents, for your construction. The pertinent documents can be downloaded from the **VOCO** website.

**V-Print model fast** is designed for prints with higher layer thickness and thus for a faster printing process.

**Processing:**

**Note:** Use separate material containers and cleaning baths for each printing material, in order to prevent cross contamination.

**Note:** Do **not** shake the material before starting the printing process. The materials container should be filled immediately before the start of the printing process. It is important to ensure that the material is free of bubbles to the extent possible, and filled to the fill level mark.

Start the print job observing the parameters that you previously selected. Once the printing process has ended, a dripping time of approximately 10 minutes is recommended. Next, carefully detach the printed objects from the build platform.

In the following steps, the printed objects will need to be cleaned, dried and post-exposed, in order to guarantee the required product characteristics.

A detailed explanation of the steps outlined above can be found under **Post-processing**. After use, **V-Print model fast** can be returned to its original or a similar container (HDPE, not light-transmissive, air-tight).

**V-Print model fast** may also be stored in material containers under dust-proof conditions and away from light. The manufacturer's specifications on the materials container – storage of remaining material – provide information on whether the printing material can be stored in the materials container that you used.

In all cases, including storage in the materials container, make sure that the remaining material is free of contaminants and polymerised residue before further use. Thus, when transferring the printing material, use a stainless steel sieve or clean the material using the 3D printer, as the case may be.

**Recommendation:** Once your work is completed, transfer the remaining material from the materials container into the original container. This allows for the materials container to be inspected and facilitates optimal storage of the printing material.

**Post-processing:**

**Cleaning**

For cleaning purposes, use isopropanol (purity ≥ 98%) as a cleaning solution in a cleaning device. An unheated ultrasonic bath or an unheated stirring bath may be used as a cleaning device.

The printed objects must be cleaned in two, or optionally in three steps. Position the unclean printed objects inside the cleaning bath so that any openings point downwards.

Use tweezers or suitable submersible baskets to fill the baths. Please ensure that the printed objects do not come into contact with one another during cleaning.

	Ultrasonic bath	Stirring bath
Pre-cleaning (optional)	Carefully pre-clean the printed objects by submerging them several times in a beaker with isopropanol.	
Preliminary cleaning*	3 minutes – may be used multiple times	3 minutes – may be used multiple times
Final cleaning	2 minutes – fresh cleaning bath	2 minutes – fresh cleaning bath

\*Note: The bath's cleaning efficacy decreases with increased use. When the cleaning efficacy decreases, the respective bath must be replaced. Resin residue on the surface could indicate that the cleaning efficacy of the bath has diminished, or that the parts came into contact with one another.

Next, the printed objects must be dried carefully using compressed air. If there is any resin residue on the printed object after the final cleaning, or if residue escapes from the undercuts when drying, the printed object can be briefly immersed once again in the final cleaning bath. Next, repeat the drying process.

**Preparation for post-exposure:**

Obstructive support structures can be removed before the post-exposure process by using a rotary instrument as close to the printed object as possible, carefully and without exerting pressure. Use a suction device. Carefully remove any remaining plastic dust using compressed air. Then, rinse the printed objects with fresh isopropanol for a few seconds. Carefully dry the printed objects once again with compressed air.

**Post-exposure:**

Conduct the post-exposure a minimum of 15 minutes after the most recent contact with isopropanol. A protective gas atmosphere is not required. It is important to ensure that the printed objects do not overlap or contact each other, as post-exposure would be negatively affected by the shadows that are cast.

Post-exposure can be conducted using the following devices:

Post-exposure device	Programme	
Xenon photoflash unit Otoflash G171	2x 2000 flashes	After 2000 flashes, observe a cooling phase of at least 2 minutes with open lid. Next, turn over and light-cure with another 2000 flashes.
UV light box LC-3DPrint Box²	30 minutes	Position the objects uniformly on the support disc. Avoid the formation of shadows.

²Or a device identical in construction. If needed, please consult your device manufacturer.

**Finishing:**

In general, please work with a low contact pressure and reduced speed. This guarantees consistent results and fewer processing marks.

In order to sand the support stubs use a fine-toothed carbide bur, for example. The bur can also be used for additional finishing of special structures.

In order to achieve precise sanding, e.g. between the support stub and printed object, it is recommended that the surface be sanded in the corresponding area with sandpaper, if necessary of different grain sizes. A similar result can also be obtained using coarser or finer silicone polishers.

**Cleaning of printing objects:**

Foreign materials (e.g. dental wax) can be removed by steaming the printing objects. Avoid excessive overheating. Don't put printing objects into the extraction unit as this could cause deformation of the printing object.

**Warnings, precautionary measures:**

Only use **V-Print model fast** in a fully polymerized state.

Pay attention to the finishing process.

Contact between uncured **V-Print model fast** and the skin/mucous membranes and eyes can cause mild irritation and should be avoided. The wearing of protective clothing is recommended. Furthermore, it is important to ensure that no vapours and/or dusts are inhaled. The wearing of a suitable mask and/or the use of suction devices is recommended. Further information on handling can be found in the safety data sheet.

Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for the intended purpose.

**Storage:**

Storage at **15°C–28°C**. Reseal bottle immediately after use. The material will cure if exposed to light. Do not use after the expiry date.

**Disposal:**

Dispose of the product in accordance with local regulations.

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VC 60 AA6125 EI 0721 99 © by VOCO
VOCO GmbH                      Phone +49 (4721) 719-0