VOCO

Die silicone

USA | EN Instrucciones for use

Carefully read instructions prior to use

Product description:

This **die silicone** is an addition-curing silicone for the extraoral fabrication of dental models. Indirect composite restorations can be fabricated chairside on this silicone model. The **die silicone's** special curing formula allows you ample time to apply the material. The material is fully cured after 4 minutes, at which time you can then start modelling the restorations.

Indications:

Extraoral fabrication of cavity models and models of individual, or groups of teeth.

Contraindications:

This **die silicone** contains a platinum catalyst. The **die silicone** should therefore not be used in case of known hypersensitivities (allergies) to this ingredient.

Method of application:

Preparation:

Insert the cartridge into a type 2 **VOCO** dispenser or comparable dispenser. Due to technical reasons there may be slight differences in the fill levels of the two cartridge chambers prior to initial use. In such a case, remove the cartridge cap and carefully extract material until it is delivered evenly from both openings. Then attach a type 20 mixing tip and lock with a 90° turn clockwise. The material is extruded by pressing the lever and is automatically mixed in the correct ratio. After use, do not remove the mixing tip from the cartridge. Replace the used mixing tip with a new one only immediately before the next application. When changing the mixing tips, check the openings of the cartridge to ensure that the material can flow freely.

Processing:

To fabricate a cavity model, or a model of an individual tooth or several teeth, first take an impression of the prepared area with an alginate. Silicones or polyethers cannot be used for this purpose. Clean and dry the impression in the usual way. Ensure all drops of moisture have been removed as these will lead to artefacts in the model and thus in the restoration. Then fill the impression with the **die silicone** within the working time. **Before use, discard the material (roughly the size of a bean) that is initially extruded from the mixing tip.**

Starting at the deepest point, place the material directly into the impression and fill the cavity from the bottom upwards leaving the tip immersed. When applying the material, ensure that it flows evenly to avoid air bubbles. **Working time at room temperature (73 °F (23 °C)) is at least 30 seconds.** Higher temperatures reduce and lower temperatures extend the working and setting times. The completed silicone model can be removed from the impression 4 minutes after commencement of mixing. After the material has set it is extremely hard, enabling you to start modelling the composite restoration immediately. Follow the instructions for use of the composite regarding preparation, processing and curing. Indirect composite restorations can then be adhesively attached in the usual way, following the respective manufacturer's instructions for use.

Warnings/precautionary measures:

- Do not use in combination with condensation-curing silicones, polyethers or materials containing eugenol.
- Setting may be impaired by the use of latex gloves, hand creams, cleaning agents, mouth rinses, retraction solutions, materials containing methacrylates, or residues of other materials. These should therefore be removed before applying the material.
- Avoid incompletely cured silicone material coming into contact with clothing.
- Please note that any residues of material left in the mixing tips will cure and set.
- The warranted product characteristics can only be achieved when using the (original Sulzer Mixpac) mixing tips as approved by **VOCO** and available as accessories.

Storage:

Store at 39 °F - 73 °F (4 °C - 23 °C). Do not use the **die silicone** after expiry date.

Product Ordering Information:

Cartridge 2 x 50ml, mixing tips type 20 REF 1186

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. **VOCO** recognizes its responsibility to replace products if proven to be defective. **VOCO** does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith. Descriptions and data constitute no warranty of attributes and are not binding.

CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

For questions or comments, please call 1-888-658-2584.

Keep this material out of reach of children.

For dental use only.

An explanation of the symbols used in labeling can be found at www.voco.dental/symbols

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VOCO GmbH Phone +49 (4721) 719-0