Description:

GrandTEC is a glass strand pre-impregnated with light-curing methacrylate resin that remains pliable until its polymerization.

GrandTEC-glass fibers are cured with conventional polymerization lights with blue light (Halogen, LED or Pac Lights). The impregnated fibers can be used on a model or placed in an impression as well as on natural teeth. The cured **GrandTEC**-glass strands increase the stability of loosened teeth. They also provide a foundation for the replacement of missing teeth and splinting existing teeth.

When incorporated into temporary crowns or bridges, they provide reinforcement to the temporary treatment. The **GrandTEC**-glass fiber must always be covered with at least one layer of a (flowable) composite for all applications!

GrandTEC can be used chairside or labside.

A combination of both application methods is recommended for some indications.

Indications:

- Stabilization of teeth:
 - After orthodontic treatment
 - After periodontal treatment
- For support in treating trauma:
 - For immobilizing and splinting avulsed and loosened teeth
 - Fragment anchoring after fracture
- For the reinforcement/fabrication of provisional bridges made from composite and provisional C & B material
- For temporary treatment after extraction and replacement of a missing tooth with the use of a plastic tooth
- For the temporary or semi-permanent treatment of a tooth space with the use of an extracted natural tooth
- For temporary treatment after placing an implant during the osseointegration of the implant

Contraindications:

GrandTEC may not be used:

- In acute periodontitis, as long as the necessary treatment has not been carried out
- In an untreated dentition with progressive, multiple carious lesions
- For replacing missing teeth in the molar range

Side effects:

There are no side effects known at our present level of knowledge.

Precautions:

- Please carefully read these instructions for use!
- Familiarize yourself with the product and its application before the clinical application. Models fabricated from impressions taken in actual clinical situations are well-suited for practicing. Practicing on a model is an excellent way to familiarize oneself with the application of **GrandTEC**-glass fiber without time constraints or the pressure to succeed.
- Ensuring a dry working field is essential. A rubber dam should be employed, if possible. The application of **GrandTEC** should be re-evaluated, if a dry working field cannot be guaranteed.
- Avoid contact between the oral mucosa and unpolymerized **GrandTEC**-glass strands.
- Contact between unpolymerized composite and skin or the eyes can lead to irritation.
- Do not use in patients with a proven allergy to methacrylate.
- Wearing gloves when working with **GrandTEC** is recommended to avoid contact with unprotected skin.
- Pay particular attention to the instructions for polymerizing the selected (flowable) composite.
- The glass fiber may require shortening with sharp scissors to the correct length before application in the mouth or on the model. Shortening the glass fiber should be carried out with the protective film intact.
- Always employ water cooling and suction when finishing a restoration that contains **GrandTEC**-glass fibers.
- Work quickly after removing the glass fiber from the packaging, since any type of light triggers polymerization. As a result, modelling the glass fiber may become difficult or entirely impossible.

Instruments for the application:

These are, for example: modelling instruments employed in the application of composite, especially the Heidemann's spatula, blunt modelling and plugging instruments, tweezers with smooth, not diamond, tips, sharp scissors for shortening the glass fiber, a small calliper, periodontal probe with wax wire, and Wedjets or similar for measuring and assigning the necessary length to the glass fiber. Additional (flowable) composites and compatible bonds are needed for processing **GrandTEC**. A fine diamond drill may be necessary for finishing (e. g., for correcting the occlusion). Rotating instruments and polishers are always needed for smoothing and creating a high gloss on the surface of the completed work.

Handling advice:

1. Splinting of teeth, direct method (chairside-use):

1. Splinting teeth is carried out on the oral or vestibular surfaces of the teeth to avoid disturbing the occlusion. The teeth must be caries-free.
2. Determine the length of the needed glass fiber using one of the instruments described above in the section Instruments for the application (e. g., a small calliper, periodontal probe, wax wire or similar).
3. Carefully clean the surfaces of the teeth to be included in the splinting. Use a rubber dam to ensure a dry working area.
4. Remove a **GrandTEC**-glass strand from the blister, leaving it in its protective film, and shorten it to the previously determined length. Lay the prepared glass strand aside, protected from light.
5. Subsequently carry out the conditioning of the enamel according to the acid-etch technique.
6. Apply an enamel-dentin bond and also wet the approximal tooth surfaces. Polymerize the applied bond according to the manufacturer's instructions.
7. Apply a flowable composite to the prepared tooth surfaces, but do not polymerize yet!
8. Remove the prepared **GrandTEC**-glass strand from the protective film with tweezers and position it on surfaces wetted with the flowable composite. **GrandTEC** can be modelled with modelling instruments, slender spatulas or similar instruments, following the tooth contour.
9. Subsequently polymerize in increments and carefully protect areas of the glass strand that have not yet been polymerized with a suitable instrument, such as a mirror or wide spatula. The wattage and polymerization times are provided in a table at the end of these instructions for use.
10. In the next step, completely cover the glass strand with the flowable composite and carefully polymerize. Pay particular attention to the approximal areas here: These must be completely filled with the flowable composite, so that every area of the incorporated **GrandTEC**-glass fiber strip is covered with composite.
11. Remove the rubber dam. Check the splinting for possible occlusal interferences, smooth and polish the flowable composite.

2. Splinting teeth, indirect method (labside-use) using a deep-drawing template:

Please note: you will need to ensure that the material being used to fabricate the deep-drawing template will not form a bond with the dental composite. Polyethylene is well suited for this purpose.

1. Take an impression of the jaw containing the teeth that should be splinted. Fabricate a model from the impression.
2. Mark the boundaries of the intended splinting on the model, while taking the patient's occlusion into consideration.
3. Create the desired shape of the splint from wax on the model. Keep in mind that the **GrandTEC**-glass strand must be completely covered with composite and prepare the model accordingly.
4. Subsequently fabricate a deep-drawing tray on the prepared model with a thickness of approximately 1.5 mm.
5. When the tray is deep-drawn, remove any residual wax and carefully clean the tray.
6. Clearly identify the outline of the planned splinting inside the deep-drawing tray during inspection. Insert the tray into the patient's mouth and make any necessary corrections.
7. Prepare the tooth surfaces as described in sections 1.3, 1.5, 1.6 and 1.7.
8. Shorten the **GrandTEC**-glass fiber to the required length. Wax wire or other similar tools can also be utilized for measuring here. Leave the glass fiber in its protective sleeve during this process.
9. Apply a small amount of a flowable composite in the deep-drawing tray where recesses were created on the model for the splinting.
10. Position the shortened **GrandTEC**-glass fiber in this composite layer and subsequently insert the deep-drawing tray into the patient's mouth.
11. Polymerize the composite through the deep-drawing tray in individual subsections.
12. Subsequently remove the deep-drawing tray and apply additional flowable composite onto the **GrandTEC**-glass strands where needed. Carefully polymerize afterwards.
13. Complete the splinting as described in section 1.11.

3. Splinting teeth with the replacement of an extracted tooth, direct method (chairside-use):

Important preliminary remark:

The **GrandTEC**-glass fiber permits an adhesive bond between two natural teeth, even over a space. With the assistance of **GrandTEC**-glass fiber, a natural, appropriately prepared tooth crown of a patient's extracted tooth or a prosthetic (replacement) tooth can be attached: Even direct modelling of a tooth with composite in the space is possible.

It is important to extend the adhesive surfaces as far apart as possible on the involved teeth, including the approximal and lingual and/or palatal surfaces of the anchoring teeth, to achieve the highest stability.

1. Inspect the clinical situation of the region where **GrandTEC** will potentially be used. Measure the length between the anchor teeth to determine the necessary length of the glass fiber, as described above.
2. Prepare the tooth that is being inserted into the space: Create the attachment sites (preparation) for the **GrandTEC**-glass fiber, making indentations and reducing the lingual or palatal surfaces, if necessary. The preparations are the same for both natural and artificial teeth.

Use of an extracted tooth:

If using an extracted tooth from the patient, first remove the root and carefully seal the root canal with composite according to the rules of the adhesive technique. Clean the tooth in advance, making sure to remove all calculus, periodontal tissue and blood. If necessary, scale down a very broad tooth base with careful preparation and rounding with diamond grinding tools.

3. Shorten a strip of **GrandTEC** to the necessary length and set aside within reach protected from light.
4. Insert the rubber dam and dry all included tooth surfaces.
5. Etch all of these tooth surfaces including the approximal areas for the duration stated by the manufacturer.
6. Rinse afterward and dry these surfaces.
7. Apply a bond to the prepared surfaces and polymerize them according to the manufacturer's instructions.
8. Wet the bonded tooth surfaces with a flowable composite.
9. Apply the prepared **GrandTEC**-glass fiber to the tooth surfaces previously wetted with composite using hand instruments or a transparent celluloid strip.
10. Polymerize the strips in sections. Refer to the table at the end of these instructions for the times. Hold the light guide emission window as close as possible to the glass fiber strips for the duration of the curing.
11. Apply a suitable amount of flowable composite to the glass strand between the abutment teeth and on the replacement tooth prepared for the space.
12. Position this tooth in the right location, check and correct the fit and carefully polymerize the composite on all sides. Refer to the manufacturer's instructions for the times.
13. Now carefully cover the entire glass strand, junctions and approximal contacts with a flowable composite and polymerize according to the manufacturer's instructions.
14. Remove the rubber dam and check the occlusion. Finish and polish the restoration.

4. Splinting teeth with replacement of an extracted tooth, indirect method (labside-use) using a deep-drawing tray:

Please note: you will need to ensure that the material being used to fabricate the deep-drawing template will not form a bond with the dental composite. Polyethylene is well suited for this purpose.

1. Make an impression of the teeth that must be splinted. Subsequently create a plaster cast.
2. Mark the boundaries of the intended splinting on the model, while taking the patient's occlusion into consideration. Create the desired shape of the splint from wax on the model. Keep in mind that the **GrandTEC**-glass strand must be completely covered with composite.
3. Prepare the tooth to be inserted: shape the tooth to a suitable size, create recesses on the adhesive side and grind the corresponding lingual and/or palatal surfaces. Subsequently attach the prepared tooth to the model.
4. Fabricate a deep-drawing tray with a thickness of 1.5 mm. Remove the tray from the model after deep-drawing, clean it and especially remove all residual wax from the replacement tooth. The outline of the splint previously modelled from wax should be clearly visible.
5. Test the deep-drawing tray with the inserted tooth in the patient's mouth. Correct the tray and/or replacement tooth, if necessary.
6. Insert the rubber dam and dry all included tooth surfaces. Etch these tooth surfaces including the approximal areas for the duration specified by the manufacturer. Subsequently rinse and dry these surfaces.
7. Shorten a strip of **GrandTEC**-glass fiber to the required length and set aside protected from light.
8. Apply a bond to the prepared surfaces and polymerize according to the manufacturer's instructions. Wet the bonded tooth surfaces with a flowable composite.
9. Apply a small amount of a flowable composite into the cavity created by the wax model in the deep-drawing tray.
10. Take the shortened **GrandTEC**-strip and likewise position it in the cavity in the tray.
11. Take the prepared replacement tooth and carefully wet the adhesive side with flowable composite and insert it into the deep-drawing tray. Position the completed deep-drawing tray in the patient's mouth immediately afterward.
12. Now polymerize the splint in sections through the deep-drawing portion. Observe the specified polymerization times. Hold the light emission window close to the section being polymerized. This prevents undesired, premature polymerization of other segments of the **GrandTEC**-glass strand and improves the polymerization in the desired segment. Subsequently remove the deep-drawing tray.
13. Cover the splinting in the mouth with flowable composite and polymerize according to the manufacturer's instructions. Fill the mesial and distal areas of the replacement tooth with flowable composite now, if this has not already been done, and polymerize.
14. Remove the rubber dam. Finish and polish the splinting and subsequently check the occlusion.

5. Reinforcing a temporary bridge with GrandTEC:

Expansive temporary treatments are significantly more resistant to fracture when supported with **GrandTEC**-glass fiber. The following procedure has been proven to incorporate such reinforcement:

1. Take the usual situation impression before preparation.
2. Measure the distance between the outside abutments of your preparation.
3. Shorten a strip of **GrandTEC**-glass fiber to the required length and set aside within reach.
4. Apply a thin layer of material for fabricating temporary crowns and bridges in the previously taken situation impression.
5. Position the shortened glass strand on the first layer and immediately fill the impression tray with your material for fabricating temporary crowns and bridges as usual.
6. Remove the impression tray from the patient's mouth after the amount of time specified by the manufacturer and remove the temporary from the impression.
7. Deviating from the other procedure, now polymerize the region in sections where the glass fiber strips have been inserted from the occlusal and basal sides, each for 30 seconds with your polymerization light. If available, place the temporary in a "light box" and polymerize it for at least 90 seconds.

Recommended polymerization times:

Type of light	LED light with 5 watts	Halogen light with 1100 mW/cm ²	Halogen light with 550 mW/cm ²	Xenon light with 250 mW/cm ²
Polymerization time	20 sec	30 sec	80 sec	240 sec

For further informations and pictograms please visit our webpage: www.vocoamerica.com

PRODUCT ORDERING INFORMATION:

GrandTEC 5 glass fiber stands 5.5 cm each

REF 1168

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. VOCO recognizes its responsibility to replace products if proven to be defective. VOCO does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith. Descriptions and data constitute no warranty of attributes and are not binding.

CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

For questions or comments, please call 1-888-658-2584.

Keep this material out of reach of children. For dental use only.

An explanation of the symbols used in labeling can be found at www.voco.dental/symbols

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VOCO GmbH