
Carefully read instructions prior to use

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Product description:

Structur CAD is a blank for the production of highly aesthetic temporary work and long-term temporaries. The maximum clinical period of wear is 3 years. Structur CAD contains 27% by weight inorganic fillers in a polymer matrix. Thanks to the composite technology, the material displays good abrasion resistance as well as high fracture resistance and is thus also ideal for long-term temporaries. Structur CAD is fluorescent. The restorations can be inserted with a temporary cement or an adhesive technique. Structur CAD is available as a block (15.5 x 19 x 39 mm) and a disc (Ø 98.4 x 20 mm).

Shades:

A1, A2, A3

Indications:

- Temporary bridges spanning up to two pontics
- Temporary abutment crowns
- Temporary crowns

Contraindications:

Structur CAD contains methacrylates. Structur CAD should not be used for patients with a known hypersensitivity (allergy) to these constituents.

Patient target group:

Structur CAD is suitable for application on all patients without any age or gender restrictions.

Application:

Structur CAD should only be applied by a professionally trained dental practitioner.

Shade selection:

Use the VITA® shade system to determine the shade against the cleaned but unprepared tooth prior to anaesthesia and preferably in daylight conditions.

Core and cavity preparation:

As a general rule, the core/cavity preparation should be performed in accordance with the conventional rules for all-ceramic restorations. That means rounding off interior corners and edges and using a shoulder preparation with rounded interior edges or a chamfer preparation. Observe the minimum thicknesses given at the start of this IFU.

Note: Do not create any undercuts.

Minimum wall thicknesses (see Table):

Note: Pontic design: height ≥ width

Grinding/milling process:

Structur CAD is suitable for both dry and wet processing. Select the corresponding block/disc size for the designed restoration and the grinding/milling parameters for the Structur CAD. When doing so, pay attention to the software settings of the respective CAD/CAM systems. "Diamond-coated tools" and composite parameters are recommended for the CAM processing. If the grinding/milling parameters are not already available in the CAD/CAM systems' software settings, they will need to be added before you proceed. Please contact the CAD/CAM system provider for assistance. The instructions for use from the manufacturers must be observed and followed.

Check the finished restoration for defects such as cracks or chips. If the restoration displays any defects, it must be rejected.

Finishing/polishing:

Once the CAM process is complete, detach the restoration from the lugs using a carbide tipped finishing bur or suitable cutting wheel. Sand the corresponding areas with fine-toothed carbide tipped finishing burs or fine diamond-coated grinders. Pay attention to any contact points.

Note: No post-processing of the shoulder to the Ti base on abutment crowns. Polish restoration chairside with conventional composite polishers. A multiphase polishing system is advantageous. A polishing paste can be used labside in combination with goat's hair brushes and cotton/leather buffers. Work using a low speed so as to avoid excessive abrasion. The instructions for use from the manufacturers must be observed and followed.

Cementation:

Preparing the restoration

For an optimal bond, roughen the surface of the restoration to be cemented using aluminium oxide (50-100 µm) at 1-2 bar or a carbide tipped finishing bur. Use a suction device to remove the dust produced. Remove abrasive material residues carefully with an ultrasonic bath (70% ethanol) or steam cleaner. Then dry the restoration with air. Final cleaning with medical alcohol is possible. The instructions for use must be observed and followed.

Temporary cementing for crowns and bridges

Cement Structur CAD with a temporary cement (e.g., Provicol QM). Note: Use a temporary, eugenol-free cement if the permanent restoration is to be subsequently cemented with an adhesive.

Adhesive cementing for crowns and bridge

For a wearing period of > 30 days, the restoration can be cemented using a resin-based cement material (e.g., Bifix QM). The instructions for use must be observed and followed. Note: - Resin-based cements can make it more difficult to remove the temporary restoration. - Regular check-ups and follow-up examinations are required for long-term use.

Extraoral cementing of the abutment crown on the titanium base

- Observe the manufacturer's specifications for the preconditioning of the Ti base (e.g., sandblasting, cleaning).
- Precondition the adhesive surface of the Ti base with a suitable bonding agent.
- Prepare the adhesive surface of the crown with a suitable bonding agent (e.g., Ceramic Bond).

When using Ceramic Bond:

- Apply Ceramic Bond to the adhesive surface with the Single Tim or a disposable brush, allow it to act for 60 seconds and dry carefully with air.
- Use a suitable opaque resin cement to stick the restoration and Ti base together. Apply the resin cement to the abutment crown and the Ti base alike. The instructions for use from the manufacturers must be observed and followed.

Customization:

For a highly aesthetic temporary restoration, the restorations can be customized, characterised or repaired at any time using a composite/ORMOCER®. Roughen the restoration surface via grinding or sandblasting (Al₂O₃ 50-100 µm, 1-2 bar). Remove abrasive material residues/dust carefully with an ultrasonic bath (70% ethanol) or steam cleaner. Then dry the restoration with air. Apply a suitable adhesive system (e.g., Futurabond U) in accordance with the instructions for use. Using GrandioSO, Flow or Heavy Flow, for example, in combination with a stain or glaze, you can customize the restorations quickly and simply with purely light-curing techniques. The instructions for use from the manufacturers must be observed and followed.

Information, precautionary measures:

Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application.

Storage:

Store at between 39°F and 73°F (4°C and 23°C). If the material is stored in the refrigerator, bring it to room temperature before use. Do not store the product in a place where it is exposed to direct artificial light or sunlight. Do not use the product after its expiration date.

Disposal:

Dispose of the product according to local regulations.

PRODUCT ORDERING INFORMATION:

Blocks 5 x No. 40L	Disc 1 x Ø 98.4 mm x 20 mm
A1 REF 6076	A1 REF 6071
A2 REF 6077	A2 REF 6077
A3 REF 6078	A3 REF 6073

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. VOCO recognizes its responsibility to replace products if proven to be defective. VOCO does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith. Descriptions and data constitute no warranty of attributes and are not binding.

CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

For questions or comments, please call 1-888-658-2584.

Keep this material out of reach of children. For dental use only.

An explanation of the symbols used in labeling can be found at www.voco.dental/symbols

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Fijación:

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For an optimal bond, roughen the luting surface of the restoration using aluminium oxide (50-100 µm) at 1-2 bar or a carbide tipped finishing bur. Use a suction device to remove the dust produced. Remove abrasive material residues carefully with an ultrasonic bath (70% ethanol) or steam cleaner. Then dry the restoration with air. Final cleaning with medical alcohol is possible. The instructions for use must be observed and followed.

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ORMOCER® is a registered trademark of Fraunhofer Gesellschaft zur Förderung der angewandten Forschung e.V., München, Germany.
VITA® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co.
VOCO GmbH Phone +49 (4721) 719-0