# V-Print® Try-In

Carefully read instructions prior to use

**Product description:**

**V-Print Try-In** is a light-curing resin for the generative manufacturing of try-ins for full and partial dentures using CAD/CAM technology.

**Color:**  
– Beige

**Indication:**  
– Try-ins for full and partial dentures  
– Transfer and grinding templates  
– Correction and occlusal impressions

**Contraindications:**

**V-Print Try-In** contains (meth)acrylate and phosphine oxide. **V-Print Try-In** should not be used on patients with known hypersensitivities (allergies) to these ingredients.

**Patient target group:**

**V-Print Try-In** is suitable for use in all patients, without any age or gender restrictions.

**Performance features:**

The product's performance features satisfy the requirements of its intended use and the relevant product standards.

**User:**

**V-Print Try-In** should only be used by a professionally trained dental practitioner.

**Hardware and software requirements**

CAD software <sup>1</sup> dental scanner	Software for the planning and design of removable dentures. The software and dental scanner must satisfy local and current medical device specifications and allow for issuance of the patient-specific design as an STL data set.  For example: – 3Shape Design Software: 3Shape Dental System® – 3Shape Dental Scanner: TRIOS, E1, E2, E3, D500, D700, D800, D900, D750, D850, D900L, D1000, D2000
CAM-software	Software for preparation of the print order. The part will not be modified during this process. Structures that facilitate the 3D printing are simply created.  For example: – Autodesk Netfabb version 2020 or later for SolFlex 3D printing.

<sup>1</sup>The designation **Software as Medical Device SaMD** includes standalone (autonomous) software that is a medical device (MD) and not part of one.

Manufacturing equipment	For example: VOCO SolFlex 170 (FlexVat) VOCO SolFlex 350 (FlexVat & PowerVat) VOCO SolFlex 650 (PowerVat)
Post-curing devices	For example: Otoflash G171

See also: accompanying list of resources or [www.voco.dental/3dprintingpartners](http://www.voco.dental/3dprintingpartners)

All manuals and/or operating instructions for the respective programs, and for device, materials and/or parts manufacturers, which are required for the manufacturing process, must be observed.

Clarify ahead of time whether the programs, devices and/or objects that you intend to use have been designed and approved for the corresponding applications.

**CAUTION:** Non-authorized changes to the process equipment, parameters, or software could result in the **V-Print Try-In** end object not satisfying specifications.

**Use:**

**Preparation:**

For an indication-appropriate CAD construction, the following design conditions must be observed:

- The minimum wall thickness of the denture base must be at least 2.0 mm
- The internal edges on the outside should be rounded off

Prepare a print job using slicing software. In addition to the material-dependent construction specifications in these Instructions for use, please also observe the dependencies of the positioning, support type, and fit found in our other documents, for your construction. The pertinent documents can be downloaded from the **VOCO** website.

**V-Print Try-In** has been conceived for a high-precision application. It is thus recommended that a small layer thickness be selected when generating the print data set.

**Processing:**

**Note:** Use separate material containers and cleaning baths for each printing material, in order to prevent cross contamination.

**Note:** Do not shake the material before starting the printing process.

The materials container should be filled immediately before the start of the printing process. It is important to ensure that the material is free of bubbles to the extent possible, and filled to the fill level mark.

Start the print job observing the parameters that you previously selected.

Once the printing process has ended, a dripping time of approximately 10 minutes is recommended. Next, carefully detach the printed objects from the build platform.

In the following steps, the printed objects will need to be cleaned, dried and post-exposed, in order to guarantee the required product characteristics. A detailed explanation of the steps outlined above can be found under **Post-processing**.

After use, **V-Print Try-In** can be returned to its original or a similar container (HDPE, not light-transmissive, air-tight).

**V-Print Try-In** can temporarily be stored in the platform under exclusion of light and dust protected. The manufacturer's specifications on the materials container – storage of remaining material – provides information on whether the printing material can be stored in the materials container that you used.

In all cases, including storage in the materials container, make sure that the remaining material is free of contaminants and polymerized residue before further use. Thus, when transferring the printing material, use a stainless steel sieve or clean the material using the 3D printer, as the case may be.

**Recommendation:** Once your work is completed, transfer the remaining material from the materials container into the original container. This allows for the materials container to be inspected and facilitates optimal storage of the printing material.

**Post-processing:**

**Cleaning**

For cleaning purposes, we recommend the use of isopropanol (purity ≥ 98%) as a cleaning solution, in a cleaning device. An unheated ultrasonic bath or an unheated stirring bath may be used as a cleaning device.

The printed objects must be cleaned in two, or optionally in three steps. Position the unclean printed objects inside the cleaning bath so that any openings point downwards.

Use tweezers or suitable submersible baskets to fill the baths. Please ensure that the printed objects do not come into contact with one another during cleaning.

	Ultrasonic bath	Stirring bath
Pre-cleaning (optional)	Carefully pre-clean the printed objects by submerging them several times in a beaker with isopropanol.	
Preliminary cleaning*	3 minutes – may be used multiple times	3 minutes – may be used multiple times
Final cleaning	2 minutes – fresh cleaning bath	2 minutes – fresh cleaning bath

\*Note: The bath's cleaning efficacy decreases with increased use. Resin residue on the surface could indicate that the cleaning efficacy of the bath has diminished, or that the parts came into contact with one another. When the cleaning efficacy decreases, the respective bath must be replaced.

Next, the printed objects must be dried carefully using compressed air. If there is any resin residue on the printed object after the final cleaning, or if residue escapes from the undercuts when drying, the printed object can be briefly immersed once again in the final cleaning bath. Next, repeat the drying process.

**Preparation for post-exposure:**

Obstructive support structures can be removed before the post-exposure process by using a rotary instrument as close to the printed object as possible, carefully and without exerting pressure. Use a suction device. Carefully remove any remaining plastic dust using compressed air. Then, rinse the printed objects with fresh isopropanol for a few seconds. Carefully dry the printed objects once again with compressed air.

**Post-exposure:**

Conduct the post-exposure a minimum of 15 minutes after the most recent contact with isopropanol. A protective gas atmosphere is not required. It is important to ensure that the printed objects do not overlap or contact each other, as post-exposure would be negatively affected by the shadows that are cast.

Post-exposure can be conducted using the following devices:

Post-exposure device	Programme	
Xenon photoflash unit Otoflash G171	2x 2000 flashes	After 2000 flashes, observe a cooling phase of at least 2 minutes with open lid. Next, turn over and light-cure with another 2000 flashes.
UV light box LC-3D Print Box <sup>2</sup>	30 minutes	Position the objects uniformly on the support disc. Avoid the formation of shadows.

<sup>2</sup>Or a device identical in construction. If needed, please consult your device manufacturer.

**Finishing:**

In general, please work with a low contact pressure and reduced speed. This guarantees consistent results and fewer processing marks.

In order to sand the support stubs use a fine-toothed carbide bur, for example. The bur can also be used for additional finishing of special structures.

In order to achieve precise sanding, e.g. between the support stub and printed object, it is recommended that the surface be sanded in the corresponding area with sandpaper, if necessary of different grain sizes. A similar result can also be obtained using coarser or finer silicone polishers.

In order to obtain a high-gloss finish, polish the object first with a pumice stone. Then, thoroughly remove any pumice dust under running water using a brush. Finally, polish the object to a high-gloss using a buffer and high-gloss polishing paste, without applying excessive pressure.

**Customization:**

In order to achieve a more aesthetic look, **V-Print Try-In** can be customized or individualized with wax, at any time. The corresponding areas will not require polishing. A covering layer with pink wax to simulate gums in the anterior region, for example, can be applied with no extra preparation.

The tooth shape can also potentially be altered using veneering material. Sand or grind the corresponding area. Carefully remove any dust residue. Then dry the restoration with air. Apply a suitable adhesive system together with veneering material in accordance with the Instructions for use.

The instructions for use of the respective customization systems must be observed.

**Preparation of the impression:**

Before filling in the impression material, it is recommended to blast the inner sides and edges of the tray with aluminum oxide (50 - 125 µm, 2 bar) and to precondition it with a suitable bonding agent. Comply with the manufacturer's instructions for use.

**Final cleaning:**

Clean the object thoroughly. First, remove any coarse residue with the steam jet cleaner. The final cleaning can be conducted by briefly placing the object in an unheated ultrasonic water bath. In order to remove oily or fatty contaminants, a surfactant solution may be used in place of water.

**Disinfection:**

Objects manufactured from **V-Print Try-In** may be disinfected using alcohol- or aldehyde-based disinfectants (e.g. ethanol (≥ 70%), Cavex ImpreSafe by Cavex). Observe the manufacturer's Instructions for use.

**Warnings, precautionary measures:**

- Only use **V-Print Try-In** intraorally in a fully cured state. Pay attention to the finishing process.
- Contact between uncured **V-Print Try-In** and the skin/mucous membranes and eyes can cause mild irritation and should be avoided. The wearing of protective clothing is recommended. Furthermore, it is important to ensure that no vapors and/or dusts are inhaled. The wearing of a suitable mask and/or the use of suction devices is recommended. Further information on handling can be found in the safety data sheet.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for the intended purpose.

**Storage:**

Storage at **59 °F - 82 °F (15 °C - 28 °C)**. Reseal bottle immediately after use. The material will cure if exposed to light. Do not use after the expiry date.

**Disposal:**

Dispose of the product in accordance with local regulations.

**PRODUCT ORDERING INFORMATION:**

Bottle 1,000 g beige

REF 6049

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use.

VOCO recognizes its responsibility to replace products if proven to be defective. VOCO does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith. Descriptions and data constitute no warranty of attributes and are not binding.

**CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.**

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

For questions or comments, please call 1-888-658-2584.

**Keep this material out of reach of children.**

**For dental use only.**

An explanation of the symbols used in labeling can be found at [www.voco.dental/symbols](http://www.voco.dental/symbols)

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VC 60 006049 US 0121 99 © by VOCO
VOCO GmbH Phone +49 (4721) 719-0