VOCO

V-Print® c&b temp

USA | EN Instructions for use

Carefully read instructions prior to use

Product description:

V-Print c&b temp is a light-curing resin for the generative production of highly aesthetic temporaries / long-term temporaries in the CAD/CAM technology. The maximum clinical wearing time is 12 month. **V-Print c&b temp** is a thixotropic material with an anorganic filler content of 26 % w/w.

Thanks to the Viscosity Change Technology, multiple swirling is sufficient to fill the tray with **V-Print c&b temp** efficiently.

V-Print c&b temp is fluorescent. The restorations can be inserted with a temporary cement or an adhesive technique.

Shades:

A1, A2, A3, A3.5, B1, BL

Indication:

– Temporary crowns and bridges

Contraindications:

V-Print c&b temp contains (meth)acrylates and phosphine oxide. **V-Print c&b temp** should not be used for patients with a known hypersensitivity (allergy) to these constituents.

Patient target group:

V-Print c&b temp is suitable for application on all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of its intended use and the relevant product standards.

Application

V-Print c&b temp should only be applied by a professionally trained dental practitioner.

Shade selection:

Use the VITA® shade system to determine the shade against the cleaned but unprepared tooth prior to anesthesia and preferably in daylight conditions.

Core and cavity preparation:

As a general rule, the core/cavity preparation should be performed in accordance with the conventional rules for all-ceramic restorations. That means rounding off interior corners and edges and using a shoulder preparation with rounded interior edges or a chamfer preparation.

Hardware and software requirements

CAD software¹ Dental scanner	Software for the planning and design of crowns and bridges. The software and dental scanner must satisfy local and current medical device specifications and allow for issuance of the patient-specific design as an STL data set. For example: - 3Shape Dental System Version 2017 or later - 3Shape Dental Scanners: TRIOS, E1, E2, E3, D500, D700, D800, D900, D750, D850, D900L, D1000, D2000
CAM-software	Software for preparation of the print order. The part will not be modified during this process. Structures that facilitate the 3D printing are simply created. For example: - Autodesk Netfabb version 2020 or later for SolFlex 3D printing.

¹The designation **Software as Medical Device SaMD** includes standalone (autonomous) software that is a medical device (MD) and not part of one.

Manufacturing equipment	For example: VOCO SolFlex 170 (PowerVat) VOCO SolFlex 350 (PowerVat) VOCO SolFlex 650 (PowerVat) VOCO SolFlex 170 HD
Post-curing devices	For example: Otoflash G171

See also: www.voco.dental/v-print-c-b-temp_ressource_us and accompanying list of resources or www.voco.dental/3dprintingpartners

All manuals and/or operating instructions for the respective programs, and for device, materials and/or parts manufacturers, which are required for the manufacturing process, must be observed.

Clarify ahead of time whether the programs, devices and/or objects that you intend to use have been designed and approved for the corresponding applications.

CAUTION: Non-authorized changes to the process equipment, parameters, or software could result in the **V-Print c&b temp** end object not satisfying specifications.

Use:

Preparation:

For a CAD construction according to the indication the following information must be observed:

Minimum wall thickness:

– occlusal 1.5 mm

– circular 1.0 mm

Bridges with a maximum of one pontic can be produced.

Connector cross-sectional area

– Anterior region ≥ 12 mm²

– Posterior Region ≥ 15 mm²

– Note: Design of the cross-sectional areas: height ≥ width

Prepare a print job using slicing software. In addition to the material-dependent construction specifications in these instructions for use, please also observe the dependencies of the positioning, support type, and fit found in our other documents, for your construction. The pertinent documents can be downloaded from the VOCO website.

V-Print c&b temp has been conceived for a high-precision application. It is thus recommended that a small layer thickness is selected when generating the print data set.

Processing:

Note: Use separate material containers and cleaning baths for each printing material, in order to prevent cross contamination.

Multiple swirling of the bottle improves the flow properties of the material and should be done immediately before printing. It is important to ensure that the material is filled into the material tray as free of bubbles as possible while observing the filling level.

Start the print job observing the parameters that you previously selected.

Once the printing process has ended, a dripping time of approximately 10 minutes is recommended. Next, carefully detach the printed objects from the build platform. In the following steps, the printed objects will need to be cleaned, dried and postexposed, in order to guarantee the required product characteristics. A detailed explanation of the steps outlined above can be found under **Post-processing**.

Make sure that the remaining material is free of contaminants and polymerized residue before further use. Thus, when transferring the printing material, use a stainless steel sieve or clean the material using the 3D printer, as the case may be.

Recommendation: Once your work is completed, transfer the remaining material from the materials container into the original container. This allows for the materials container to be inspected and facilitates optimal storage of the printing material.

Post-processing:

Cleaning

Remove unpolymerized resin residues on the print objects using an isopropanol (purity ≥ 98 %) soaked brush.

Next, the printed objects must be dried carefully using compressed air. If there is any resin residue on the printed object after the final cleaning, or if residue escapes from the undercuts when drying, the printed object should again be cleaned with an isopropanol (purity ≥ 98 %) soaked brush.

Preparation for post-exposure:

Obstructive support structures should be removed before the post-exposure process by using a rotary instrument as close to the printed object as possible, carefully and without exerting pressure. Use a suction device. Remove remaining resin dust carefully with compressed air and, if necessary, with a brush and isopropanol (purity ≥ 98 %). Carefully dry the printed objects once again with compressed air.

Post-exposure:

Conduct the post-exposure a minimum of 15 minutes after the most recent contact with isopropanol. A protective gas atmosphere is not required. It is important to ensure that the printed objects do not overlap or contact each other, as post-exposure would be negatively affected by the shadows that are cast.

Post-exposure can be conducted using the following devices:

Post-exposure device	Program	
Xenon photoflash unit Otoflash G171	2x 2000 flashes	After 2000 flashes, observe a cooling phase of at least 2 minutes with open lid. Next, turn over and light-cure with another 2000 flashes.

See also: www.voco.dental/v-print-c-b-temp_ressource_us and accompanying list of resources

Further print objects should only be post-exposed after the device has cooled down for 10 minutes in order to avoid discoloration of the material.

Finishing/polishing:

Use, for example, a fine-toothed carbide cutter to grind down the support attachments. This can also be used for subsequent elaboration of special structures. In order to achieve a form-fitting grinding e.g. between the support attachment and printed object, it is recommended to treat the surface in the corresponding area with sandpaper, if necessary with different grains. An appropriate result can also be achieved with coarser or finer silicone polishers. Polish restoration chairside with conventional composite polishers. A multiphase polishing system is advantageous.

A polishing paste can be used labside in combination with goat's hair brushes and cotton/leather buffers.

Use only light pressure and wiping motions to control the removal. The instructions for use from the manufacturers must be observed and followed.

Customization:

For a highly aesthetic temporary restoration, the restorations can be customized, characterized or repaired at any time using a composite/ORMOCER®. Roughen the restoration surface via grinding or sandblasting (Al₂O₃ 50-100 µm, 1-2 bar). Remove abrasive material residues/dust carefully with an ultrasonic bath (70% ethanol) or steam cleaner. Then dry the restoration with air. Apply a suitable adhesive system (e.g., **Futurabond U**) in accordance with the instructions for use. Using **GrandioSO**, **Flow** or **Heavy Flow**, for example, you can customize the restorations quickly and simply with purely light-curing techniques.

The instructions for use from the manufacturers must be observed and followed.

Luting:

Preparing the restoration

For an optimal bond, roughen the luting surface of the restoration using aluminum oxide (50-100 µm) at 1-2 bar or roughen with a fine HM milling device.

Use a suction device to remove the dust produced. Remove abrasive material residues carefully with an ultrasonic bath (70% ethanol) or steam cleaner. Then dry the restoration with air. Final cleaning with medical alcohol is possible. The instructions for use must be observed and followed.

Temporary luting for crowns and bridges

Lute **V-Print c&b temp** with a temporary cement (e.g., **Provicol QM**). **Note:** Use a temporary, eugenol-free cement if the permanent restoration is to be subsequently luted with an adhesive.

Adhesive luting for crowns and bridges

For a wearing period of > 30 days, the restoration can be luted using a composite- based adhesive luting material (e.g., Bifix QM). The instructions for use must be observed and followed.

Note: – Adhesive luting causes time-consuming removal of the temporary restoration.
– Regular check-ups and follow-up examinations are required for long-term use.

Warnings, precautionary measures:

– Only use **V-Print c&b temp** intraorally in a fully cured state. Pay attention to the finishing process.

– Contact between uncured **V-Print c&b temp** and the skin/mucous membranes and eyes can cause mild irritation and should be avoided. The wearing of protective clothing is recommended. Furthermore, it is important to ensure that no vapors and/or dusts are inhaled. The wearing of a suitable mask and/ or the use of suction devices is recommended. Further information on handling can be found in the safety data sheet.

– Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for the intended purpose.

Storage:

Storage at **59 °F - 82 °F (15 °C - 28 °C)**. Reseal bottle immediately after use. The material will cure if exposed to light. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's condition and a serious risk to public health that arise or could have arisen in association with the use of **V-Print c&b temp** must be reported to VOCO GmbH and the responsible authority.

PRODUCT ORDERING INFORMATION:

Bottle 500 g A1	REF 6897
Bottle 500 g A2	REF 6898
Bottle 500 g A3	REF 6899
Bottle 500 g A3.5	REF 6900
Bottle 500 g B1	REF 6901
Bottle 500 g BL	REF 6902

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use.

VOCO recognizes its responsibility to replace products if proven to be defective. VOCO does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith. Descriptions and data constitute no warranty of attributes and are not binding.

CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

For questions or comments, please call 1-888-658-2584.

Keep this material out of reach of children.

For dental use only.

An explanation of the symbols used in labeling can be found at www.voco.dental/symbols

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VC 60 006897 US 0723 99 © by VOCO