

The tip:  
designed  
for **EASY**  
application into  
the sulcus!

The cord in the capsule!

## VOCO Retraction Paste

ASTRINGENT PASTE FOR THE EFFECTIVE TEMPORARY RETRACTION  
AND MOISTURE MANAGEMENT OF THE GINGIVAL SULCUS

# Retraction Paste

## ASTRINGENT PASTE FOR THE EFFECTIVE TEMPORARY RETRACTION AND MOISTURE MANAGEMENT OF THE GINGIVAL SULCUS

Perfect impressions, moisture control and gentle soft tissue management are essential for exact crown margins, perfect marginal adaptation, and long-lasting, esthetic restorations.

### The essential element of taking simple, effective, precision impressions!

VOCO Retraction Paste is an aluminum chloride-containing astringent paste that effectively widens and dries the gingival sulcus, creating the ideal conditions for accurate and successful digital or conventional impressions.

VOCO Retraction Paste's special unit dose capsule was engineered with two goals in mind: first, to make application easy and precise and second, to avoid injury to the gingiva. Therefore, VOCO designed a capsule with an extra-long, fine tip made of soft plastic instead of metal. This ensures that the opening of the tip does not have any sharp edges. The design of the capsule and the special rheology of the material allows for the reduction of the force required to extrude the paste. Each unit dose capsule is packaged in a single-dose ISO blister pack, reducing hygiene concerns. At 0.3 g per cap, there is enough material for up to three sulci.

VOCO Retraction Paste's two-phase consistency makes application very simple. First, the paste easily flows under low pressure for precise application. After application, the material's viscosity increases so as to remain stable in the sulcus, stop bleeding, displace moisture and ultimately lead to effective temporary widening.

Compared to packing cord alone, Retraction Paste is not only easier, it is also up to 50% faster. After the short application time of 1–2 minutes, the high-contrast VOCO Retraction Paste can be rinsed off easily. The paste can be used as an alternative or in combination with retraction cord. It is up to you whether you take a classic impression with impression material or a digital impression with an intraoral scanner.

### Indications

For temporary retraction of the marginal gingiva to provide a dry gingival sulcus when the periodontium is healthy prior to:

- analog or digital impressions
- cementation of temporary and permanent restorations
- preparation of class II and V restorations

### Advantages

- **Aluminum chloride-based** – for a clean and dry sulcus
- **Unit dose with extra-long, thin and flexible tip** – for hygienic, simple and precise application into the sulcus
- **2-Stage Viscosity** – initially flowable for easy, low-force application, then stabilizes for optimal retraction
- **Turquoise blue color** – contrasts well with the gingiva
- **Easily rinseable** – for quick clean-up
- **Up to 50% faster** – compared to just cord

### Presentation

REF 1012 Caps 25 × 0.3 g

REF 1013 Caps 100 × 0.3 g



VOCO America Inc.  
1285 Rosemont Drive  
Indian Land · SC 29707

www.vocoamerica.com  
infousa@voco.com

Toll-free phone: 1-888-658-2584  
Fax: 1-888-849-3989