# Admira Bond

Carefully read the instructions prior to use

## Product description:

**Admira Bond** is a single-component, light-curing bonding agent that produces a durable bond between filling materials and hard dental tissues (dentin/enamel). **Admira Bond** is based on complex, 3-dimensionally polymerized inorganic/organic co-polymers (ORMOCER's) that have been developed specifically for use with the total-etch technique.

When used with **Admira**, an ORMOCER-based, light-curing filling material, **Admira Bond** establishes an exceptionally durable bond. **Admira Bond** can, however, be universally used under light-curing composite filling materials.

A special feature of **Admira Bond** is that it requires only a single application within the course of a total-etch procedure. **Admira Bond** combines quick and simple application with effective, durable bonding strength.

## Indications:

**Admira Bond** is intended for use as an adhesive bond between hard dental tissue (dentin/enamel) and light-curing restorative materials. **Admira Bond** is suitable for direct restorations of all classes of fillings with light-curing restorative materials.

## Contraindications:

Do not use in patients with a known allergy to methacrylates, BIS-GMA, BHT, HEMA, acetone or organic acids.

## Application:

The adhesive bond is created within the framework of the total-etch technique.

### 1. Drying

Since light-curing restorative materials are sensitive to moisture, the tooth surface must be thoroughly dried. Do not allow the cavity to become contaminated with saliva or blood. The use of a rubber dam is highly recommended.

### 2. Preparation

Prepare the cavity in accordance with standard adhesive restoration techniques. Areas close to the pulp should be protected with a suitable liner over a calcium hydroxide base (e.g. **Calcicur**, **Calcimol**).

### 3. Etching

Etch the entire cavity with a suitable etching gel (e.g. **Vococid**), starting from the enamel margin (etching enamel for approximately 20 - 30 sec; etch dentin for no more than a maximum of 15 sec). Remove etching gel by suction, rinse for approximately 20 sec. and then remove excess moisture with an air jet. Do not overdry the dentin surface. The dentin surface should remain slightly moist. The etched enamel should have a chalk-white/matt appearance. If not, repeat the etching procedure.

### 4. Dispensing Admira Bond

When tilting the bottle for dispensing, make sure to keep the vent hole pointed upwards. Insert **Micro Tim** or other suitable applicator into the dispenser and let soak. Alternatively, dispense **Admira Bond** onto a mixing pallet and then thoroughly wet the applicator tip.

### 5. Applying Admira Bond

Apply a thin, even layer of **Admira Bond** on all prepared enamel and dentin surfaces using **Micro Tim** or a disposable brush applicator. Allow to stand for 30 sec. Then, gently disperse the **Admira Bond** with a faint air jet and polymerize with a halogen light for 20 sec. Do not remove or contaminate the inhibition layer because this layer is important for adhesion between the dental surface and the composite. **Micro Tim** fine brushes enable precise application of a thin layer of **Admira Bond**.

### 6. Filling

Carefully apply **Admira** or other selected composite filling material to the cavity walls. Build the filling up layer by layer, polymerizing each layer, and then finish. Follow manufacturer instructions for each restorative composite.

## Precautions:

- Phenolic substances (e.g. eugenol-containing temporary materials) inhibit polymerization and must be kept out of contact with **Admira Bond**.
- Avoid contact with the oral mucosa.
- Avoid contact with the eyes. Should contact occur, flush the eye thoroughly with water and consult an ophthalmologist.

## Storage:

Tightly close **Admira Bond** immediately after each use.

Storage at temperatures 39 - 73°F (4 - 23°C).

Storage at 41 - 47°F (5 - 8°C) prolongs the shelf-life of this product.

Do not expose **Admira Bond** to direct sunlight or to the chairlight.

## Product Ordering Information:

### Bottles:

- 2 x 4ml Admira Bond REF 2432
- 1 x 4ml Admira Bond REF 2479

### Micro Tim:

- Pack 100 Micro Tim REF 2249
- 5 Brushholder for Micro Tim REF 2338

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. VOCO recognizes its responsibility to replace products if proven to be defective. VOCO does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

Descriptions and data constitute no warranty of attributes and are not binding.

**CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.**

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

For questions or comments, please call 1-888-658-2584, 604-274-1975 or 604-377-9569.

**Keep this material out of reach of children.  
For dental use only.**

An explanation of the symbols used in labeling can be found at [www.voco.dental/symbols](http://www.voco.dental/symbols)

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VOCO GmbH