Carefully read instructions prior to use

**Product description:**

**Admira Fusion x-tra** is a light-curing, radiopaque nano-hybrid ORMOCER® restorative material. It is characterized by an increased depth of cure while at the same time showing reduced shrinkage stress, which means that it can be cured in layers of 4 mm thickness.

**Admira Fusion x-tra** contains 84 % w/w inorganic fillers and is used with a dentin-enamel bond. It is available in syringes and Caps for direct application.

**Shade:**

Universal (U)

**Indications:**

- Class I and II posterior restorations
- As a base in class I and II cavities
- Class V restorations
- Locking, splinting of loose anteriors
- Repairing veneers, small enamel defects and temporary C&B materials
- Extended fissure sealing
- Restoration of deciduous teeth
- Core build-up

**Contraindications:**

**Admira Fusion x-tra** contains methacrylates and BHT. **Admira Fusion x-tra** should not be used in cases of known hypersensitivities (allergies) to any of these ingredients.

**Preparation:**

Clean the teeth to be treated. If required mark contact points. Allow the material to reach room temperature before use.

**Cavity preparation:**

In principle, cavity preparation should be carried out according to the rules of adhesive filling therapy and should be minimally invasive to conserve healthy tooth structure. Clean and dry the cavity afterwards. Non-carious cervical lesions do not have to be prepared; thorough cleaning is sufficient here.

**Creating a dry working field:**

Ensure that the work area is sufficiently dry. Use of a rubber dam is recommended.

**Matrix placement:**

Use of matrices is recommended for cavities with an approximal portion. Translucent matrices that are wedged in the approximal area are advantageous. Minimal separation facilitates the shaping of the approximal contact and insertion of the matrix.

**Pulp protection:**

A suitable pulp-protective liner should be applied in close proximity to the pulp.

**Bond material:**

**Admira Fusion x-tra** is used in the adhesive technique with a dentin-enamel bond. All light-curing bond materials can be used. Follow the respective instructions for use with regard to preparation (etching technique) and application.

**Application of Admira Fusion x-tra:**

Apply **Admira Fusion x-tra** in layers that are a maximum of 4 mm thick, adapt with a suitable instrument and light-cure afterwards.

**Light-curing:**

Conventional polymerization devices are suitable for light-curing the material. When using an LED-/halogen light with a minimum light output of **800 mW/cm<sup>2</sup>** the curing time is **20 s**, with a minimum light output of **500 mW/cm<sup>2</sup>** the curing time is **40 s**. Hold the light emission tip of the device as close as possible to the surface of the filling. Otherwise, the curing depth may be reduced. Incomplete curing may lead to discoloration and discomfort.

**Finishing:**

The restoration can be finished and polished immediately after removing the matrices (e.g., fine or extra-fine diamond bur, polisher). The tooth should be fluoridated as a final step.

**Information/precautionary measures:**

- Phenolic substances, especially preparations containing eugenol or thymol, interfere with the curing of filling composites. Avoid the use of zinc oxide eugenol cements or other eugenol-containing materials in combination with filling composites.
- Each **Admira Fusion x-tra** Caps unit should be used for one patient only and should be applied exerting uniform but not excessive pressure.
- Fillings exposed to occlusal forces should be checked at least once a year for early detection of any changes.

**Storage:**

Store at temperatures between 39 °F - 73 °F (4 °C - 23 °C). To avoid exposure to light and potentially resulting polymerization, syringes should be closed immediately after dispensing.

Do not use after expiry date.

**PRODUCT ORDERING INFORMATION:**

Syringe 3gm universal  
Caps 15 x 0.2gm universal

REF 2810  
REF 2811

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. VOCO recognizes its responsibility to replace products if proven to be defective. VOCO does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

Descriptions and data constitute no warranty of attributes and are not binding.

**CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.**

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

For questions or comments, please call 1-888-658-2584.

**Keep this material out of reach of children.  
For dental use only.**

An explanation of the symbols used in labeling can be found at [www.voco.dental/symbols](http://www.voco.dental/symbols)

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VOCO GmbH