Amaris Gingiva

Carefully read instructions prior to use

Amaris Gingiva is a light-cured radiopaque restorative which contains 80% fillers in a methacrylate matrix and which can be cured with blue light. **Amaris Gingiva** is available in syringes. **Amaris Gingiva** is used with a dentin/enamel bond.

Amaris Gingiva is used to treat cervical defects and exposed, discolored or hypersensitive necks of teeth, especially in the visible anterior area. The gingival-like shade and the **Opaquer** system allow the correct reproduction of the gingiva also in its shade. The system comprises three different **Opaquers** and a **gingiva-shaded composite restorative in a basic color**.

The light-cured **Opaquers** are used to cover discolored cervical areas and to individually adapt the restoration to the surrounding gingiva.

The light-cured **Amaris Gingiva** allows a simplified and fast restorative technique.

Indications:

- Class V restorations (cervical caries, root erosion, V-shaped defects) especially in cases of retreating gingiva e.g. as a result of periodontic disease.
- Covering of exposed discolored hypersensitive necks of teeth, especially in the visible anterior area.
- Highly aesthetic correction of malpositions of teeth, for direct veneers, facings and corrections of the red-white aesthetics.

Contraindications:

Amaris Gingiva should not be used if the prescribed application technique cannot be used or if there are known allergies to its ingredients (BIS-GMA, BHT, UDMA, HEDMA, TEGDMA).

Method of use:

Preparation:

- Clean the teeth to be treated.
- Ensure a sufficiently dry workfield. The use of rubber dam is recommended.

Cavity preparation:

Generally, a minimally invasive preparation (adhesive restorative technique) should be carried out to conserve tooth substance.

It is not necessary to prepare caries-free cervical defects. Here, cleaning with a cleaning paste or pumice and a rubber cup or rotating brushes is sufficient.

The use of a retraction cord can be helpful.

Cleaning:

- Remove all residues from the cavity with the water jet.
- Any contamination with blood or saliva after cleaning has to be avoided.
- Remove water residues with a faint airjet, but do not overdry dentin. The aim is to achieve a moist, but not wet dentin surface (moist bonding).

Color design:

Dark or strongly yellowish necks of teeth impair the color effect of the restoration. Therefore, these cervical areas should be covered with the enclosed **Opaquer** system before applying the filling.

The set contains a shade guide which serves as an orientation to show the color effect of the **Opaquers** combined with the restorative. The shade guide is made from the following components:

- White - for a light gingival color
- Light - for a medium gingival color
- Dark - for a dark gingival color
- Nature - the basic gingival color

The color can be individually adapted step by step by mixing the **Opaquers**.

Bond material:

Amaris Gingiva is applied with a dentin/enamel bond in the adhesive technique. All light-curing bond materials can be used. Observe the instructions for use of the bond chosen.

Light-curing:

For light-cured the material market-available polymerization lights are suited. The output should not be lower than 500 mW/cm².

Application of the Opaquer:

Apply some **Opaquer** onto the tooth surfaces to be covered and disperse **thinly**, then light-cure for **40 sec**. If necessary, apply a second coat of **Opaquer** to cover the tooth surface, and light-cure for **40 sec**. In order to ensure the function of the non-dripping syringe do not pull back the plunger during or after use.

Application of Amaris Gingiva:

Apply fillings of more than 2 mm thickness in layers. Light-cure each layer with blue light for **40 sec**.

Use transparent strips and light-wedges. Make sure to obtain good marginal adaptation. Hold the end of the light-rod of the polymerization device as close to the filling as possible. If this distance is more than 5 mm, curing depth may be reduced and limited to the light-cone. Incomplete curing may lead to discoloration and pulpitis-like complaints or restoration failure. Contouring and finishing of the filling can be conducted immediately after removal of the shaping device, using cooling (e.g. fine or extrafine diamond polishers, polishing disks). Finally, the entire filling margin, preferable the entire tooth should be fluoridated.

Information/precautionary measures:

- No known side effects. However, sensitivity in hypersensitive patients cannot be excluded.
- Phenolic substances, especially eugenol- and thymol-containing preparations, interfere with the curing of filling composites. Do not use zinc oxide eugenol cements or other eugenol-containing materials in combination with **Amaris Gingiva**.
- The enclosed application cannulae are for single use and have to be disposed of after use. Cover the syringe with its cap after use.

Storage:

Store at room temperature 39°F - 73°F (4°C - 23°C). Close syringes immediately after dispensing. Do not use **Amaris Gingiva** after its expiry date.

PRODUCT ORDERING INFORMATION:

Amaris Gingiva kit:

1-4gm syringe (nature), 3-1.2gm Opaque (white, light, dark), accessories REF 1970

Amaris Gingiva refill: 1-4gm syringe (nature) REF 1971

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use.

VOCO recognizes its responsibility to replace products if proven to be defective. VOCO does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

Descriptions and data constitute no warranty of attributes and are not binding.

CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.

No person is authorized to provide any information which deviates from the information provided the instructions for use.

For questions or comments, please call 1-888-658-2584.

Keep this material out of reach of children.

For dental use only.

An explanation of the symbols used in labeling can be found at www.voco.com/us/symbols

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VOCO GmbH