# VOCO

## Clip Flow

USA | EN

### Instructions for use

Carefully read instructions prior to use

#### Product description:

**Clip Flow** is a flowable, light-curing material for temporary fillings, sealings and treatments. Thanks to its elastic consistency, the material is easily removable, even in case of undercuts. Therefore, post-treatment of the cavity is not necessary.

#### Indications:

- [A] Temporary inlay and onlay treatments of the cavity
- [B] Sealing of openings for implant screws
- [C] Relining material for temporary crowns and bridges
- [D] Block-out material for retentive areas in the dental arch, e.g. before taking impressions
- [E] Covering of the gingival margin
- [F] Fixing of resin matrix during filling placement
- [G] Temporary filling of cavities

#### Contraindications:

**Clip Flow** contains methacrylates and BHT. **Clip Flow** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.

#### Preparation:

If necessary, carry out preparation, lining and impression. Adhesion to dental resins or linings based on (meth)acrylate may occur. In such cases apply a suitable insulating agent (e.g. glycerine, vaseline). Avoid the use of lining materials containing eugenol since these could impede the setting of **Clip Flow**. Rinse and dry the cavity before application of **Clip Flow**.

#### Application:

Apply **Clip Flow** directly from the NDT® syringe using the respective cannula (type 46). In order to ensure the functioning of the non-dripping syringe do not pull back the plunger during or after use. If necessary, spread **Clip Flow** with an adequate instrument (e.g. probe). Light-cure with a polymerization light (minimum 500 mW/cm<sup>2</sup>) to a film thickness up to 1 mm for 10 s, (e.g. in case of covering the gingival margin), a thickness up to 5 mm for 20 s and more than 5 mm for 40 s. Hold the light emission tip of the device as close as possible to the surface of the filling. Otherwise, the curing depth may be reduced.

#### Finishing:

Excess material can be removed with a scalpel or round bur. Silicone points are suitable for finishing.

#### Removal of the temporary filling:

Insert an appropriate instrument (e.g. probe or curette) sideways into the temporary filling and extract **Clip Flow** from the cavity.

#### Warnings/precautions:

- In case of eye contact, rinse eyes immediately with plenty of water and seek medical advice.
- The enclosed application cannulae are for single use only and must be discarded after application. Close the syringes after use with the cap provided.
- Phenolic substances, especially preparations containing eugenol or thymol, inhibit the curing of **Clip Flow**. Therefore, do not use zinc oxide eugenol cements or other eugenol-containing substances in combination with **Clip Flow**.

#### Storage:

Store at temperatures of 39 °F to 73 °F (4 °C to 23 °C). To avoid exposure to light and potentially resulting polymerization, syringes should be closed immediately after dispensing.

Do not use after expiry date.

#### Product Ordering Information:

Syringe 2 x 1.8 gm, accessories REF 1282  
Application cannulae type 46, 100 pcs. REF 2148

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use.

VOCO recognizes its responsibility to replace products if proven to be defective. VOCO does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

Descriptions and data constitute no warranty of attributes and are not binding.

**CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.**

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

For questions or comments, please call 1-888-658-2584.

**Keep this material out of reach of children.**

**For dental use only.**

An explanation of the symbols used in labeling can be found at [www.voco.dental/symbols](http://www.voco.dental/symbols)

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VOCO GmbH Phone +49 (4721) 719-0