Carefully read instructions prior to use

Product description:

Fit Test C & B is an addition-curing silicone for checking the precision of fit of crowns, bridges, inlays/onlays, secondary parts and attachments, as well as for temporary blocking-out during final cementing of dentures. For the dentist **Fit Test C & B** is a diagnostic help for the optimum fitting of the finished restorative piece. With its blue color it allows to clearly visualize faulty areas on metal and full ceramic restorative pieces.

Fit Test C & B convinces with its simple handling and specifically high tear resistance.

Indications:

Control of precise fit, marginal tightness and marginal gap of

- full or partial crowns
- bridges
- inlays/onlays
- cast post cores/post crowns
- metal bases
- implant-supported crowns and bridges
- attachments and other secondary parts of dentures

Contraindications:

Do not use **Fit Test C & B** on patients with known hypersensitivities against platinum catalyst. Do not use **Fit Test C & B** on existing permanently soft silicone-based base materials (e.g. relining) since a chemical bond to these materials may occur.

The following materials are able to delay or prevent the setting reaction of **Fit Test C & B**: condensation-curing silicones, polysulfide impression materials, eugenol-containing materials, not fully cured synthetic materials as well as latex gloves. When using **Fit Test C & B** these materials may not be used.

Application:

Clean and dry restoration, preparation and abutment as usual.

Due to technical reasons slight differences in the filling level between the syringe chambers can occur prior to the first use. Therefore remove syringe cap and squeeze the material out until the material is delivered evenly from both chambers. Match the nooks on the syringe and the mixing tip and fix the mixing tip by a 90° turn clockwise. The material is mixed automatically in the correct ratio when dispensing and can be applied directly. All QM syringes come with a detachable piston. The piston can simply be pulled from the QM syringe and re-placed, without causing any leakage of material. The plugs in the syringe chambers will stay in place. This prevents any cross-contamination of the material in the syringe.

Fit Test C & B is applied directly from the mixing tip or in case of inaccessible areas with an instrument to the cleaned restorative piece or cavity. Cover all surfaces with **Fit Test C & B** and avoid air bubbles. Place the restorative piece in situ and fix (do **not** let the patient bite!).

The working time from the beginning of the mixing is at least 1 min; the intra-oral setting time is 2 min maximum – monitor setting by checking excess material. After setting remove carefully the restorative piece with the attached silicone film. Small one-surface inlays can also be removed with sticky wax.

Evaluation:

Even silicone layers show a good precision of fit. Is the silicone layer uneven, the restorative piece or the preparation have to be corrected where contact areas shine through, using standard techniques.

Precautionary measures:

Avoid contamination (e.g. of textiles) with non-cured silicone since it may be difficult to remove. The above mentioned working times apply for a temperature of approximately 73 °F (23 °C) at a relative humidity of approx. 50%. Lower temperatures prolong these times, higher temperatures shorten them.

The warranted product characteristics can only be achieved when using the supplied original mixing tips, which are also available as accessories.

Storage:

Store at temperatures 39 °F - 73 °F (4 °C - 23 °C).

PRODUCT ORDERING INFORMATION:

QuickMix syringe 5ml, mixing tips type 9

REF 2095

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use.

VOCO recognizes its responsibility to replace products if proven to be defective. VOCO does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

Descriptions and data constitute no warranty of attributes and are not binding.

CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

For questions or comments, please call 1-888-658-2584.

Keep this material out of reach of children.

For dental use only.

An explanation of the symbols used in labeling can be found at www.voco.dental/symbols

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VOCO GmbH