Product description:

Futurabond NR is a light curing self-etch-adhesive reinforced with nano-particles to build a durable bond between tooth substance and light-curing filling materials without marginal leakage. This all-in-one adhesive consists of 2 components. A separate etching step is not necessary: **Futurabond NR** shows the same adhesive properties as total-etch bonds.

Futurabond NR tolerates residual moisture.

Futurabond NR can be used universally with all light-curing composites, compomers and ORMOCER®s.

Indications:

Futurabond NR is suited for direct restorations with light-curing filling materials for all cavity classes.

Contraindications:

Futurabond NR contains organic acids, BisGMA, HEMA, TMPTMA, Campherchinon, amines (DABE), BHT, fluorides and ethanol.

Futurabond NR should not be used in case of known hypersensitivities (allergies) to any of these ingredients.

Precautionary measures:

Phenolic substances (e. g. eugenol-containing temporary materials) disturb the polymerization and should be avoided when using **Futurabond NR**.

Avoid contact with the oral mucosa. In case of contact of **Futurabond NR** with the eyes, rinse immediately with plenty of water and consult a physician.

Do not fold blister before activation. Blisters with obvious damage such as creases or kinks should not be used.

Do not damage either the blister foil or tab on either the top or the bottom side with fingernails.

Do not use any type of instrument for activation. Only the thumb and index finger should be used for activation with the thumb placed exactly on the mark and the index finger placed exactly on the opposite side of the package.

Activation and subsequent puncturing of the foil with a **Single Tim** should not take place within close proximity of the eyes.

Never simultaneously press the both reservoirs or press the liquid back and forth between the reservoirs.

Methods of use:

1. Drying

Due to the subsequent application of light-curing filling materials, drying of the work field is required. Avoid contamination with blood or saliva.

2. Preparation

Prepare the cavity according to the principles of adhesive restoration technique. Rinse cavity thoroughly with water. Remove excess moisture with faint air jet. Do not thoroughly dry dentin. The dentinal surface should remain slightly moist.

If the tooth has not been prepared (e.g. cervical defect), the surface must be thoroughly cleaned and the enamel margin must be roughened with a finishing diamond.

Areas close to the pulp should be protected by a suitable liner based on calcium hydroxide.

3. Mixing

Futurabond NR is light curing and must be immediately applied after mixing.

When using the bottle:

Mix one drop of Liquid **A** and one drop of Liquid **B** on a mixing palette with a **Micro Tim** for approximately 5 sec. to produce a self-etching adhesive.

When using the SingleDose:

Detach a *SingleDose* by folding along the perforation and pull the printed side upwards. To activate, hold the *SingleDose* blister between thumb and forefinger and press on the blister where indicated "press here" in blue. This injects Liquid **B** into Liquid **A** in the dispensing bay in a 1:1 ratio. Pierce the foil of the dispensing bay with the **Single Tim**, enlarge the hole with a circular movement and create a homogeneous mixture of **A** and **B** by repeated stirring.

4. Application of Futurabond NR

Apply a moderately thin layer of adhesive to the enamel/ dentin and massage into the tooth substance for 20 sec.

5. Drying

Dry the adhesive layer with an air jet for a minimum of 5 sec.

6. Polymerization

Polymerize with a blue light (halogen- or LED-light) for 10 sec.

Information on storage and application:

Tightly close **Futurabond NR** after use. Store in an upright position in order to guarantee flow-back of the liquid. Store at temperatures 39 °F - 73 °F (4 °C - 23 °C). Allow **Futurabond NR** to reach room temperature before use. Do not expose **Futurabond NR** to direct sunlight or chairlight. Do not use after expiry date.

Product ordering Information:

Bottle

- Set each 4ml liquid A/ B, accessories REF 1148

SingleDose

- Kit 50 *SingleDose*, accessories REF 1152

- Pack 200 *SingleDose*, accessories REF 1153

Single Tim

- Pack 100 **Single Tim** REF 2247

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. VOCO recognizes its responsibility to replace products if proven to be defective. VOCO does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

Descriptions and data constitute no warranty of attributes and are not binding.

CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.

No person is authorized to provide any information which deviates from the information provided the instructions for use.

For questions or comments, please call 1-888-658-2584.

Keep this material out of reach of children.

For dental use only.

An explanation of the symbols used in labeling can be found at www.voco.dental/symbols

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VOCO GmbH Phone +49 (4721) 719-0