# VOCO


## Grandio® disc

USA | EN


Instructions for use

Carefully read instructions prior to use

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

VITA® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co.

#### Product description:

**Grandio disc** is a nano-hybrid composite CAD/CAM disc for the fabrication of permanent single-tooth restorations. **Grandio disc** contains 86 % w/w inorganic fillers in a polymer matrix. The restorations produced must be placed with an adhesive composite luting agent in combination with a dentin-enamel bonding agent.

**Grandio disc** is available in one size (ø 98.4 x 15 mm).

#### Shades (two degrees of translucency):

High Translucent (HT): A1, A2, A3, A3.5

Low Translucent (LT): A1, A2, A3, A3.5, B1, C2, BL

Multicolor: A1, A2, A3, A3.5, B1, C2

#### Indications:

- Crowns, inlays, onlays, veneers
- Implant-supported crowns

#### Contraindications:

**Grandio disc** contains methacrylates. **Grandio disc** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.

#### Shade selection:

Use the VITA® shade system (**GrandioSO**) to determine the shade against the cleaned but unprepared tooth prior to anaesthesia and preferably in daylight conditions.

#### Core and cavity preparation:

As a general rule, the core/cavity preparation should be performed in accordance with the conventional rules for all-ceramic restorations.

That means rounding off interior corners and edges and using a shoulder preparation with rounded interior edges or a chamfer preparation.

The temporary restoration and any contamination should be mechanically removed (e.g. using pumice and a polishing brush) and the area subsequently rinsed with water.

Observe the minimum thicknesses stated below.

**Note:** Do not create any undercuts.

#### Crowns, inlays and onlays (see Fig. 1 - 3):

- Preparation angle of 4° - 6°
- Minimum occlusal wall thickness 1.5 mm
- Minimum buccal wall thickness 0.8 mm
- Minimum wall thickness below load-bearing cusp 1.5 mm
- Minimum cervical wall thickness 0.8 mm

#### Veneers (see Fig. 4):

- Minimum cervical wall thickness 0.4 mm
- Minimum labial wall thickness 0.5 mm
- Minimum incisal wall thickness 0.5 mm

#### Cement space:

70 µm (± 10 µm)

#### Grinding/milling process:

Select the correct disc size for the planned restoration as well as the milling parameters for **Grandio disc** or for composites. When doing so, pay attention to the software settings of the respective CAD/CAM systems. **Processing requires diamond instruments.** First of all, the settings for disc sizes and the milling parameters need to be created, if they are not already included in the software settings of the CAD/ CAM systems. Please contact the CAD/CAM system provider for assistance. The operating instructions from the manufacturers must be observed and followed.

Please note the software settings of the respective CAD/CAM systems and the direction of the multicolor restoration.

Check the processed restoration for defects such as fissures or material chipping. If the restoration displays any defects, it must be rejected.

**Note: Grandio disc is an already completely cured composite and must not be fired.**

#### Polishing:

The restoration can be finished and polished with standard composite polishers (e.g. medium, fine or extra fine diamond grinding tools, polishers), this can be done either intraorally or extraorally.

#### Adhesive luting:

##### Preparing the restoration:

For an optimal bond, blast the luting surface of the restoration with aluminium oxide (25 µm - 50 µm) at 1.5 bar - 2 bar. Use a suction device to remove the dust produced.

**Note:** Do not etch with hydrofluoric acid or phosphoric acid. Clean the pretreated restoration using an ultrasonic bath or steam cleaner and then dry with oil-free air. Final cleaning with medical alcohol is possible. Treat the luting surface with a silane coupling agent such as **Ceramic Bond**. Apply **Ceramic Bond** to the inside of the work and allow it to work for 60 seconds. Then dry for 5 seconds with oil-free air. The instructions for use must be observed and followed.

##### Adhesive luting of the restoration:

The processed restoration must be luted with a composite-based, adhesive luting agent. The instructions for use must be observed and followed. **VOCO** recommends the dual-curing luting system **Bifix QM** in combination with **Ceramic Bond** for an optimal bond to the restoration and a suitable dentin-enamel bond (**Futurabond DC/Futurabond U**). Alternatively, it is also possible to use a purely light-cured, flowable composite/ORMOCER® (e.g., **GrandioSO Flow**) if the restoration has thin walls (≤ 2 mm).

#### Modifications:

**Grandio disc** is a nano-hybrid composite CAD/CAM restorative material which can be modified, characterised or repaired at any time using a composite/ORMOCER®. Roughen the surface of the restoration (for example with a diamond bur) and apply a suitable adhesive system (e.g., **Futurabond U**) in accordance with the instructions for use. Then apply the composite/ORMOCER® (e.g. **GrandioSO**, **GrandioSO Flow**, **FinalTouch**) in accordance with the respective instructions for use.

**Information, precautionary measures:**

No known side effects. However, the sensitization of hypersensitive patients cannot be ruled out.

**Storage:**

Store at temperatures of 39 °F to 73 °F (4 °C to 23 °C).

If the material is stored in the refrigerator, allow it to reach room temperature before use.

Do not store the product in a place where it is exposed to direct artificial light or sunlight.

Do not use after expiry date.

**PRODUCT ORDERING INFORMATION:**

1 x Ø 98.4 mm x 15 mm

A1 LT	REF 6050	A1 HT	REF 6057
A2 LT	REF 6051	A2 HT	REF 6058
A3 LT	REF 6052	A3 HT	REF 6059
A3.5 LT	REF 6053	A3.5 HT	REF 6060
B1 LT	REF 6054		
C2 LT	REF 6055		
BL LT	REF 6056		

1 x Ø 98.4 mm x 15 mm multicolor

A1	REF 6216
A2	REF 6217
A3	REF 6218
A3.5	REF 6219
B1	REF 6220
C2	REF 6221

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. VOCO recognizes its responsibility to replace products if proven to be defective. VOCO does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith. Descriptions and data constitute no warranty of attributes and are not binding.

**CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.**

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

For questions or comments, please call 1-888-658-2584.

**Keep this material out of reach of children.**

**For dental use only.**





An explanation of the symbols used in labeling can be found at [www.voco.dental/symbols](http://www.voco.dental/symbols)

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VC 60 006051 US 0723 99 © by VOCO
ORMOCER® is a registered trademark of Fraunhofer Gesellschaft zur Förderung der angewandten Forschung e.V., München, Germany.
**VOCO GmbH** Phone +49 (4721) 719-0