Carefully read instructions prior to use

Product description:

Merlon Plus QM is a radiopaque, self-curing, fluoride-releasing, paste-to-paste resin-modified glass ionomer luting cement available in the practical QuickMix syringe for quick and simple direct application.

Merlon Plus QM offers optional removal of excess material following light-curing.

Indications:

Cementation of

- metal-based inlays, onlays, crowns and bridges
- high-strength (zirconium dioxide-based) all-ceramic crowns and bridges
- high-strength all-ceramic inlays
- metal, ceramic and fiber posts

Contraindications:

Merlon Plus QM contains polyacrylic acid, peroxide, methacrylates and BHT. **Merlon Plus QM** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.

Method of Use:

1. Preparation:

Ensure that the surfaces to be cemented are clean. Thoroughly clean and dry enamel and dentin, and the inside surfaces of the restoration. Before cementation, clean the cores (e.g., pumice, water), rinse thoroughly and dry, ideally with a foam pellet or gentle blasts of air (take care not to overdry!). Prepare and treat the restoration as per the manufacturer's instructions.

Remove the protective cap from the QuickMix syringe. Due to technical reasons, there may be slight differences in the filling levels of the syringe chambers prior to first use. For this reason, when using a new syringe, dispense paste until the material exits both openings evenly. When reusing syringes, check the openings for any blockages and unblock as necessary. Attach a type 10 mixing tip in such a way that the notches on the mixing tip and the syringe are aligned. Turn the mixing tip 90° clockwise to lock the tip into place. To ensure even mixing of the two pastes, dispense and discard a little material before commencing the application.

Curing behavior:

Working time:	2 min (73 °F (23 °C), from start of mixing)
Removal of excess material:	a) in gel state (approx. 2 min after placement in mouth) or b) following light-curing (approx. 5 - 10 s per segment)
Setting time:	4.5 min (after placement in mouth)

2. Placing the restoration:

Apply material to the preparation and/or the prepared restoration. Then place the restoration, applying gentle pressure so as to expel excess material. At 73 °F (23 °C), the working time is 2 min from the start of mixing. Remove the used mixing tip and seal the syringe with the original twist cap until the next time it is required.

3. Removal of excess material:

Excess cement is best removed in the gel phase, either following light-curing (5 - 10 s per segment) or following self-curing (approx. 2 min after placement), using a suitable instrument (e.g., scaler). Excess material in hard-to-reach areas should be removed first. Clean contact points and interproximal areas with dental floss. No material residue should remain in the sulcus, as this can result in irritation of the gingiva.

4. Curing:

Curing is complete 4.5 min after placement. The restoration must be held in position during curing. Cured excess material can be removed with a finishing diamond or flexible polishing wheels. Finish the approximal areas with finishing/polishing strips.

Information, precautionary measures:

In case of contact with the eyes, rinse immediately with plenty of clean water and seek medical advice.

Ensure sufficient residual dentin thickness to prevent pulp irritation. If necessary, protect the pulp or build up the core. Do not use any phenolic substances (e.g., temporary cements containing eugenol), as these can impair the curing of **Merlon Plus QM**. Only use eugenol-free temporary cements for temporary cementation.

Refit the original twist cap to prevent the **Merlon Plus QM** pastes drying out during storage.

The warranted product characteristics can only be achieved when using the supplied original mixing tips, which are also available as accessories.

Storage:

Store **Merlon Plus QM** at temperatures of 39 °F to 73 °F (4 °C to 23 °C) with the original twist cap fitted. Do not use after expiry date.

Allow **Merlon Plus QM** to reach room temperature before use!

PRODUCT ORDERING INFORMATION:

QuickMix syringe 8.5 g, accessories	REF 1727
QuickMix syringe 3 x 8.5 g, accessories	REF 1728
Mixing tips type 10, 50 pcs.	REF 2189

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. VOCO recognizes its responsibility to replace products if proven to be defective. VOCO does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to

determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

Descriptions and data constitute no warranty of attributes and are not binding.

CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

For questions or comments, please call 1-888-658-2584.

Keep this material out of reach of children.

For dental use only.

An explanation of the symbols used in labeling can be found at www.voco.dental/symbols

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VOCO GmbH