# VOCO Quick Up

USA | EN

## Instructions for use

Carefully read instructions prior to use

### Product description:

The kit includes **Quick Up**, a gingiva-colored, self-curing luting composite for attachments. This material is used for bonding of attachments such as ball, locator® and telescopic anchors as well as other anchors in composite-based dentures.

The material can also be used for reattaching secondary elements in a denture such as bar retainers. The kit also includes **Quick Up Adhesive** for permanent and firm bonding, **Quick Up LC**, a gingiva-colored, light-curing correction material for subsequent corrections, and **Fit Test C & B** for checking the accuracy of fit and blocking out.

### Fields of application:

Fixing attachments such as ball, Locator® and telescopic anchors, as well as other anchors in partial and full dentures both chairside and in the laboratory. The kit can also be used as a repair material for reattaching secondary elements (e.g. bar retainers) in dentures.

### Please note:

Attachments and other secondary components must have mechanical retention or must have been roughened mechanically, e.g. by sandblasting, before being bonded, if there is no such retention.

### Contraindications:

Do not use on patients with known allergies to the constituents of the products.

### Application:

#### 1. Preparing the denture:

Place all attachments onto the abutments and check their fit carefully. The position of the attachments can be transferred onto the denture with the help of articulating paper, a suitable felt-tip pen, occlusal spray, toothpaste or the like. Make oval recesses in the denture of sufficient depth to receive the attachments. The attachments must not touch the denture base. Approximately 2 mm space must be left for **Quick Up**. Undercuts in the denture base increase the hold of the attachments. Optionally, lingual openings can be made below the teeth to allow excess material to flow out which enables a visual check.

#### 2. Checking accuracy of fit of the attachments in the denture:

The enclosed **Fit Test C & B** can be used to check whether the openings in the denture base provide enough space to receive the attachments. The instructions for use for **Fit Test C & B** should also be followed when applying the product. Fill the openings for receiving the attachments with the **Fit Test C & B** and place the denture in the patient's mouth. Allow the patient to bite with moderate pressure and wait 2 minutes for the **Fit Test C & B** to set. Remove the denture and check the results: an even silicone layer around the attachment indicates that there is sufficient space in the denture base to receive the attachment. If there is an uneven silicone layer or if the attachment is in contact with the denture base, the openings in the denture must be enlarged and the test repeated if necessary. Remove the attachments and remove **Fit Test C & B** completely from the denture using suitable instruments (e.g. scaler or curette).

#### 3. Insulating and blocking out:

Using Vaseline for example insulate all areas that are not to be permanently bonded, i.e. surfaces of dentures, abutments/implants, as well as the neighboring teeth. To prevent an accidental and irreversible adhesion of the denture to implants or abutments in the mouth, undercuts and the cervical free gingiva around the implants must be blocked out or covered with the **Fit Test C & B** or another suitable, stable material such as silicone or wax. To do this apply **Fit Test C & B** around the undercuts and the cervical gingiva of abutments or other male parts and wait 2 minutes for it to set. Use the enclosed mixing tips type 11 to which the intraoral tip type 1 is attached.

#### 4. Applying the Quick Up adhesive:


Dispense a few drops of **Quick Up Adhesive** into the mixing tray and apply them to the parts of the denture to be bonded using the enclosed **Single Tim** application brush. Then dry the **Quick Up Adhesive** layer with an air syringe or dry for 30 seconds in the air. The bottle must be resealed immediately after use and stored upright.

#### 5. Preparing the Quick Up cartridge:

Twist the protective cap on the double-chamber syringe of the **Quick Up** to 90° and then remove it. For technical reasons there may be slight differences in the fill levels of both cartridges chambers. When opening a new cartridge, you should therefore continue to press out the material until an equal amount has been squeezed out of both outlets. Then attach the enclosed mixing tip type 10 and rotate to 90° to lock it in place. The material is mixed in the mixing tip as it is dispensed and can be applied immediately. After use leave the mixing tip on the cartridge or replace the protective cap.

#### 6. Applying Quick Up:

Quickly fill with **Quick Up** a maximum of two thirds of each of the prepared openings in the denture base described in points 1 and 2. You have an extraoral working time of 30 seconds. Place the denture in the patient's mouth and allow the patient to bite with moderate pressure. Check the fitting of the denture and occlusion. You have 30 seconds to carry out this check in the mouth. After placing the denture wait 2.5 minutes before removing the denture from the mouth (a total time of 3.5 minutes from the start of mixing). **Quick Up**'s curing time will be extended if adjustments are being made to a model because the prevailing temperature will be lower than that of the oral cavity. The setting time can be controlled by spreading the **Quick Up** on a mixing pad. Higher temperatures accelerate curing. Use suitable rotary instruments if excess material has to be removed, ensuring that the attachments do not come into contact with the rotary instruments. If insufficient material has been applied or if there are areas requiring correction, use the **Quick Up LC** which is enclosed in the kit (see point 7).

#### 7. Making corrections with the Quick Up LC:

The **Quick Up LC** contained in the kit is used for correcting shortfalls and surface defects. To make these corrections clean the surfaces and apply the **Quick Up Adhesive**. Then dry the **Quick Up Adhesive** layer with an air syringe or dry in the air for 30 seconds. The bottle must be resealed immediately after use and stored upright. To apply the **Quick Up LC** lock the enclosed application tip type 44 in place on the black syringe by twisting it in a clockwise direction. Replace the protective cap after use. The plunger must not be retracted during or after use, otherwise the syringe may drip or run. The material must be applied in layers and cured if corrections exceeding 1 mm in thickness are necessary. Allow at least 10 seconds curing time per layer. Conventional polymerization devices are suitable for light-curing **Quick Up LC**. The light output of the LED halogen lamp should be at least 500 mW/cm<sup>2</sup>.

#### 8. Finishing and final inspection:

The material used for blocking out must be removed from the abutments after bonding the attachments. The surfaces are finished using the usual instruments (carbide bur) and polished with a flannel buffer. Then place the denture in the patient's mouth and check its fit and the occlusion.

### Warning/precautionary measures:

Phenolic substances (e.g. eugenol) impair polymerization and therefore must not come in contact with the materials. **Quick Up Adhesive** contains hydroxyethyl-methacrylate and acetone. **Quick Up** contains benzoyl peroxide, amines, BHT and methacrylates. **Quick Up LC** contains methacrylates and amines. Do not use on patients with known allergies to the constituents of the products. Avoid skin contact. Any areas of the skin that have come into contact with the products should be washed thoroughly with soap and water. In case of eye contact, rinse well with plenty of water and consult an ophthalmologist. The warranted product characteristics can only be achieved when using the supplied original mixing tips, which are also available as accessories.

### Storage instructions and application method:

Store at 39°F - 73°F (4°C - 23°C). If the material is stored in the refrigerator, allow it to assume ambient temperature before use. Do not use the products after their expiry date. Always keep the **Quick Up Adhesive** bottle tightly sealed and store upright. The black syringe used with the **Quick Up LC** should be sealed after use. Do not expose the materials to direct sunlight or the operating light. The enclosed application and mixing tips, as well as the brushes, are intended for single use only and must be disposed of after each use.

### Product ordering information:

Kit 7.5gm QuickMix syringe <b>Quick Up</b> , 2gm <b>Quick Up LC</b> , 4ml <b>Quick Up adhesive</b> , 5ml <b>Fit Test C &amp; B</b> , accessories	REF 1625
7.5gm QuickMix syringe <b>Quick Up</b> , accessories	REF 1626
2gm syringe <b>Quick Up LC</b> , accessories	REF 1627
4ml bottle <b>Quick Up</b> adhesive	REF 1628

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. VOCO recognizes its responsibility to replace products if proven to be defective. VOCO does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

Descriptions and data constitute no warranty of attributes and are not binding.

**CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.**

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

For questions or comments, please call 1-888-658-2584.

**Keep this material out of reach of children.**

**For dental use only.**

US Patent No. 8,641,416 B2

An explanation of the symbols used in labeling can be found at [www.voco.dental/symbols](http://www.voco.dental/symbols)

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**VOCO GmbH**