# VOCO

## Ufi Gel® hard C

USA | EN

### Instructions for use

Carefully read instructions prior to use

#### Product description:

**Ufi Gel hard C** is a cold-curing hard relining material for permanently relining dentures. It is simple and quick to use for direct as well as indirect relinings. An adhesive for obtaining a durable bond is included in the set.

#### Indications:

- hard permanent total or partial relinings for restoring partial and complete dentures
- lengthening denture margins

#### Contraindications:

- **Ufi Gel hard C** contains benzoyl peroxide (catalyst), BIS-GMA, HEDMA, UDMA and acetone (adhesive). Do not use in cases of hypersensitivity (e. g. allergies) to any of these ingredients.
- If the patient is uncertain about allergies, a patch test should be performed by extruding a small amount of **Ufi Gel hard C**, allowing it to set, applying it to the patient's arm and covering it piece of gauze or with a bandaid for 24 hours. DO NOT apply uncured **Ufi Gel hard C** to the skin.

#### Methods of use:

##### 1. Preparation of the denture

Check occlusion and carry out corrections if necessary. Thoroughly clean denture with a brush and dry. Roughen all areas to be relined incl. buccal and labial margins with a suitable bur. Then clean and dry. Insulate all areas not to be relined with e.g. vaseline (e.g. artificial teeth, labial and buccal area beyond the roughened surfaces) to avoid unwanted contact with **Ufi Gel hard C**.

**IMPORTANT: Partial dentures with clasps, anchors and attachments must be blocked out with thin flowing silicone or wax with low melting point. Otherwise, the denture cannot be removed once the material has completely polymerized.**

##### 2. Application of the adhesive

Exchange the transportation cap for brush cap. Apply the adhesive to all surfaces to be relined and let air dry (approx. 30 sec.).

**CAUTION: Immediately close the bottle tightly after use and store in an upright position.**

**IMPORTANT: For a safe bond use only Ufi Gel hard C Adhesive.**

##### 3. Preparation of the cartridge

Put the **Ufi Gel hard C** cartridge into a VOCO dispenser (type 2) or any dispenser with a suitable plunger (1:1). Slight differences in the filling level of the chambers may occur. Remove the transport cap from the cartridge and extrude material until both chambers deliver uniformly if this happens. Then fix and lock the mixing tip with a 90° turn. The material is extruded and automatically mixed in the correct ratio by pressing the dispenser lever. After use, leave the mixing tip on the cartridge or use the transport cap.

##### 4. Application of Ufi Gel hard C

Apply **Ufi Gel hard C** evenly onto the prepared margin and/or base of the denture, avoiding excess.

**CAUTION: Do not apply an excessive amount of Ufi Gel hard C to the surface. Otherwise, the patient may swallow Ufi Gel hard C upon application to the denture. Remove excess material with a suitable instrument.**

Reinsert denture and have the patient exert slight pressure in occlusal position for 1 minute, then carry out functional, chewing and swallowing movements for 2 minutes.

**WARNING: Make sure that no material flows into the throat when relining, especially at the transition from hard to soft palate (A-line).**

**Tell the patient to alert the dentist if, at any time, they feel relining material moving towards their throat.**

##### 4.1 Partial and total denture with undercuts

**IMPORTANT: Strictly observe the specified timing**

Remove denture after 3.5 minutes after the begin of application and remove excess immediately with scissors or a scalpel. Reinsert into the mouth for a final occlusion and check/let cure for 2 minutes. Instead of reinserting the denture, curing can also be completed in warm water, i.e. in a pressure polymerization unit at approx. 104 °F (40 °C).

**CAUTION: Do not cure material by contact with air since oxygen causes an uncured inhibition layer on the surface.**

##### 4.2 Total dentures without undercuts

**IMPORTANT: Strictly observe the specified timing**

Excess material can be removed intra or extra orally before final polymerization.

Intraorally: remove excess material after 3 minutes (beginning of application) with a suitable instrument. Let the denture cure for a further 2 minutes in the mouth until **Ufi Gel hard C** is completely cured.

Extraorally: see section 4.1

#### 5. Finishing

The whole relined surface (also basal) must be treated by finishing (tungsten carbide bur, silicone polisher) and polishing (polishing paste, polishing disc).

beginning of application	0.5 min	apply
insert denture into the mouth	1 min	exert pressure
	2 min	functional movements
remove denture if necessary	1.5 min	remove excess
insert denture into the mouth	2 min	curing
remove denture	2.5 min	finishing, polishing
	<b>9.5 min</b>	

- extraorally     intraorally     intra- or extraorally

#### Precautions:

- Avoid contact of adhesive to skin. If contact occurs, thoroughly wash skin with soap and water.
- In case of contact with the eyes, thoroughly rinse with water and consult an ophthalmologist.
- The stated time periods refer to a room temperature of 75°F (≈ 25°C) and a temperature in the mouth of 95°F (≈ 35°C). Polymerization will be slower at lower temperatures, higher temperatures will accelerate polymerization.
- The warranted product characteristics can only be achieved when using the supplied original mixing tips, which are also available as accessories.

#### Storage:

Keep containers tightly closed and store between 39°F (≈ 4°C) - 73°F (≈ 23°C). Do not use after the expiration date. Shelf life of the components are imprinted on packaging.

Store tightly closed adhesive bottle in an upright position.

#### Cleaning indications:

Dentures relined with **Ufi Gel hard C** can be cleaned with the usual cleansers and procedures. This applies to domestic as well as professional cleaning.

#### Product ordering information:

Set cartridge 80 g, adhesive bottle 10 ml, accessories	REF 2215
Set + dispenser cartridge 80 g, adhesive bottle 10 ml, accessories	REF 2218
Cartridge 80 g, accessories	REF 2216
Adhesive bottle 10 ml	REF 2217
Dispenser – type 2	REF 2158
Mixing tips type 8, 50 pcs.	REF 2187

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. VOCO recognizes its responsibility to replace products if proven to be defective. VOCO does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

Descriptions and data constitute no warranty of attributes and are not binding.

**CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.**

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

For questions or comments, please call 1-888-658-2584.

**Keep this material out of reach of children.**

**For dental use only.**

An explanation of the symbols used in labeling can be found at [www.voco.dental/symbols](http://www.voco.dental/symbols)

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VOCO GmbH    Phone +49 (4721) 719-0