VisCalor® bulk

in accordance with DIN EN ISO 4049

Carefully read instructions prior to use

Product description:

VisCalor bulk is a light-curing, radiopaque nano-hybrid restorative material. It is characterized by an increased depth of cure while at the same time showing reduced shrinkage stress, which means that it can be cured in layers of 4 mm thickness.

VisCalor bulk contains 83 % w/w inorganic fillers and is used with a dentin-enamel bond.

VisCalor bulk is pre-warmed before application. The increase in temperature leads to reduced viscosity and thus effortless insertion into the cavity as well as very good adaptation of the composite to the cavity walls. When filling of the cavity is complete, viscosity increases again as the composite cools down. In this state the material is very easy to sculpt.

Shades:

Universal (U), A1, A2, A3

Indications:

- Class I and II posterior restorations
- Base in class I and II cavities
- Class V restorations
- Locking, splinting of loose teeth
- Repairing veneers, enamel defects and temporary C&B-materials
- Extended fissure sealing
- Restoration of deciduous teeth
- Core build-up

Contraindications:

VisCalor bulk contains methacrylates and BHT. **VisCalor bulk** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.

Preparation:

Clean the teeth to be treated. If required mark occlusal contact points. Bring a preheating device for composites (e.g. **Caps Warmer**) to 155 °F (68 °C) in accordance with the operating instructions. Insert caps and allow them to warm up (approx. 3 min.). Leave in the preheating device for max. one day.

Alternatively, a **VisCalor Dispenser** can be used. Please refer to its accompanying instructions for use.

Shade selection:

Clean the teeth prior to shade selection. The shade is selected in comparison with the tooth while it is still moist.

Cavity preparation:

In principle, cavity preparation should be carried out according to the rules of adhesive filling therapy and should be minimally invasive to conserve healthy tooth substance. Clean and dry the cavity before treatment. Non-carious cervical lesions do not have to be prepared; thorough cleaning is sufficient here.

Creating a dry working field:

Ensure that the work area is sufficiently dry. Use of a rubber dam is recommended.

Matrix placement:

Use of matrices is recommended for cavities with an approximal portion. Translucent matrices that are wedged in the approximal area are advantageous.

Minimal separation facilitates the shaping of the approximal contact and insertion of the matrix.

Pulp protection:

A suitable pulp-protective liner that can be covered with a stable cement, if necessary, should be applied in close proximity to the pulp.

Bond material:

VisCalor bulk is used in the adhesive technique with a dentin-enamel bond. Any light-curing bond material can be used. Follow the respective instructions for use with regard to preparation (etching technique) and application.

Application of VisCalor bulk:

After warming it in the **Caps Warmer**, insert the Cap into the Applicator and apply the material within 20 s. After 20 s, **VisCalor bulk** still has a temperature of 142 °F (61 °C). When using the **VisCalor Dispenser**, select programme 1 (149 °F / 65 °C), which has a warm-up time of 30 seconds. The material is then kept warm for 2 min 30 sec. Please also refer to the instructions for use of the **VisCalor Dispenser**. Insert **VisCalor bulk** directly into the prepared cavity, starting at the lowest point, and fill the cavity from bottom to top while keeping the tip of the Caps submerged. Ensure that the material is applied in a slow and steady flow to prevent air bubbles from forming. Apply **VisCalor bulk** in layers that are a maximum of 4 mm thick, adapt with a suitable instrument and light-cure afterwards

Light-curing:

Conventional polymerization devices with a wavelength in the range of 400 - 500 nm are suitable for light-curing the material.

Depending on the light output and shade the curing time for increments of 2 - 4 mm is:

LED-/halogen light	Universal (U)	A1, A2, A3
≥ 1000 mW/cm ²	10 s	20 s
≥ 500 mW/cm ²	20 s	40 s

The curing time for increments of **up to 2 mm** is 10 s for all shades when using an LED-/halogen light with a minimum light output of **500 mW/cm²**. Hold the light emission tip of the device as close as possible to the surface of the filling. Otherwise, the curing depth may be reduced. Incomplete curing may lead to discoloration and discomfort.

Finishing:

The restoration can be finished and polished immediately after removing the matrices (e.g., fine or extra-fine diamond bur, polisher). The tooth should be fluoridated as a final step.

Information/precautionary measures:

- Phenolic substances, especially preparations containing eugenol or thymol, interfere with the curing of filling composites. Avoid the use of zinc oxide eugenol cements or other eugenol-containing materials in combination with filling composites.
- Each **VisCalor bulk** Caps unit should be used for one patient only and should be applied exerting uniform but not excessive pressure. Caps should only be heated once (**VisCalor Dispenser**).
- Fillings exposed to occlusal forces should be checked at least once a year for early detection of any changes.
- Extract **VisCalor bulk** after warming in a preheating device, as high extrusion forces occur in the unwarmed state. The device must be allowed to cool down completely (approx. 2-3 minutes) between two uses with the **VisCalor Dispenser** so as to avoid overheating.

Storage:

Store at temperatures of 39 °F to 73 °F (4 °C to 23 °C). Do not use after expiry date.

PRODUCT ORDERING INFORMATION:

Set Dispenser 80 x 0.25gm (16 x universal, 16 x A1, 32 x A2, 16 x A3, REF 6062
 VisCalor Dispenser

Set Caps Warmer 80 x 0.25gm (16 x universal, 16 x A1, 32 x A2, 16 x A3, REF 6063
 Caps Warmer

16 x 0.25gm in the following shades

universal (U)	REF 6065
A1	REF 6066
A2	REF 6067
A3	REF 6068

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. VOCO recognizes its responsibility to replace products if proven to be defective. VOCO does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

Descriptions and data constitute no warranty of attributes and are not binding.

CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

For questions or comments, please call 1-888-658-2584.

**Keep this material out of reach of children.
 For dental use only.**

An explanation of the symbols used in labeling can be found at www.voco.dental/symbols

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VOCO GmbH Phone +49 (4721) 719-0