

Rebilda Post – Revision of an endodontic filling

VOCO GmbH, Knowledge Communication Department

Anton-Flettner-Str. 1-3
27472 Cuxhaven, Germany

Tel.: +49 (0)4721-719-1111
Fax: +49 (0)4721-719-109

info@voco.de
www.voco.com



Conventional ceramic or metal posts cannot be drilled out, but can only be trepanned with a special, tube-like, active root post (e.g. Komet RepairPost®), which means further loss of tooth substance in an already weakened tooth. The professional journals qualify the use of glass fibre posts as a fundamental advantage, because they can easily be drilled out in the event of a needed revision and facilitate a gentler procedure.

Determining the orientation and length of the post in question by taking an x-ray is necessary for the revision of a glass fibre post. The high radiopacity of the Rebilda Post is truly indispensable here.

- 1.) After opening the coronal entrance, the exposed post is first centrally pre-drilled to half the length of the post using a ca. 1 cm long, small diameter pilot drill. The correct direction of this course should be radiographically verified.
- 2.) The post is now hollowed out with preferably a reamer (Torpan type, 15,000 rpm) which has approximately half the diameter of the post.
- 3.) The corresponding Rebilda Post drill is used for final drilling. Radiographic verification is also necessary for this step. After rinsing and disinfection, a new post can be adhesively luted.

This technique is described in the literature^[1,2] and also recommended by other manufacturers.

Conclusion: Post endodontic treatments with glass fibre posts have an inherently low fracture predisposition due to their dentine-like loading behaviour. In the event of a necessary revision, however, the radiopaque Rebilda Post additionally facilitates the preservation of the most tooth substance possible.

[1] R. S. Schwartz, J. W. J. Robbins, *J. Endod.* **2004**, *30*, 289-301.

[2] W. G. de Rijk, *Am. J. Dent.* **2000**, *13*(Spec No), 19B-21B.