VOCO

Admira® Fusion x-base

EN Instructions for use
MD EU Medical device
 In accordance with DIN EN ISO 4049

**Product description:**

Admira Fusion x-base is a flowable, light-curing, radiopaque nano-hybrid ORMOCEP® base restorative material. It is characterised by an increased depth of cure while at the same time showing reduced shrinkage stress, which means that it can be cured in layers of 4 mm thickness. When used as a base/liner **Admira Fusion x-base** has to be overlaid occlusally with a covering layer, of at least 2 mm, of a universal or posterior composite/ORMOCER®.

Admira Fusion x-base contains 72 % w/w inorganic fillers and is used with a dentine-enamel bond.

Shade:

Universal (U)

Indications:

- Base in class I and II cavities
- Cavity lining under direct restorative materials in class I and II cavities
- Small, non-occlusal stress-bearing class I restorations according to minimally invasive filling therapy
- Class III and V restorations
- Extended fissure sealing
- Undercut blockout
- Repair of small enamel defects
- Repair of small defects in aesthetic indirect restorations
- Repair of temporary C&B-materials
- Core build-up

Contraindications:

Admira Fusion x-base contains methacrylate and BHT. **Admira Fusion x-base** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

Admira Fusion x-base is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

Admira Fusion x-base should only be used by a professionally trained dental practitioner.

Use:**Preparation:**

Clean the teeth being treated. If applicable, mark occlusal contact points. Bring material to room temperature before use.

Cavity preparation:

In principle, the cavities should be prepared according to the rules for adhesive restorative treatment; preparation should be minimally invasive in order to protect healthy dental hard tissue. For anterior restorations, angle off the enamel margins. Then clean and dry the cavity. Caries-free lesions in the cervical area do not need to be prepared; thorough cleaning is sufficient in this case.

Creating a dry working field:

Ensure that the work area is sufficiently dry. Use of a rubber dam is recommended.

Matrix placement:

Take appropriate measures for areas where the applied **Admira Fusion x-base** can flow away. Translucent matrices that are wedged in the approximal area are advantageous. Minimal separation facilitates the shaping of the approximal contact and placement of the matrix.

Pulp protection:

A suitable pulp-protective liner that can be covered with a stable cement, if necessary, should be applied in close proximity to the pulp.

Bonding materials:

Admira Fusion x-base is applied with a dentine enamel bond, using the adhesive technique. Any light-curing bonding materials can be used. Observe the relevant instructions for use during preparation (etching technique) and processing.

Application:

Lock the enclosed application tip on the syringe with a clockwise turn and directly apply **Admira Fusion x-base**. Do not pull the plunger back on the syringe during or after use to guarantee the function of the non-dripping syringe.

Begin applying the material directly into the prepared cavity at the deepest point and fill the cavity from bottom to top. Keep the tip of the metal cannula immersed all times. To prevent the formation of air bubbles ensure a slow, even flow. If necessary, air bubbles or excess material can be removed before polymerisation.

When used as a base/liner the distance to the occlusal margin of the cavity must be at least 2 mm. **Admira Fusion x-base** can be polymerised in increments of up to 4 mm.

The attached application tip is intended for single use only. Remove it after use and seal the syringe again tightly with its original cap. The use of a protective barrier is recommended to protect the syringe against contamination with bodily fluids, dirt from hands or oral tissue.

Light-curing:

Conventional polymerisation devices are suitable for light-curing the material. When using an LED-/halogen light with a minimum light output of 800 mW/cm² the curing time is 20 s, for a minimum light output of 500 mW/cm² the curing time is 40 s.

Bring the light aperture as close as possible to the light-curing unit as otherwise the curing will be less effective. Inadequate curing can cause discolouration and complaints.

Final occlusal layer:

When used as a base/liner **Admira Fusion x-base** has to be overlaid with a covering layer, of at least 2 mm, of a universal or posterior composite/ORMOCER®. For working and finishing directions follow the manufacturer's instructions for use.

Finishing:

The restoration can be finished and polished immediately after removing the matrices (e.g., fine or extra-fine diamond burs, polisher). The tooth should be fluoridated as a final step.

Warnings, precautionary measures:

- Phenolic substances, especially products containing eugenol and thymol, lead to curing failures in restorative composites. Consequently, the use of zinc oxide/eugenol cements or other materials containing eugenol should be avoided in combination with restorative composites.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):

Barium aluminium borosilicate glass, ORMOCEP® resin, silicon dioxide, fumed silica, initiators, stabilisers, pigments

Storage instructions and application method:

Storage at 4°C-23°C. Seal the syringes immediately after removing the material in order to avoid exposure to light and subsequent polymerisation. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Admira Fusion x-base** must be reported to VOCO GmbH and the responsible authority.

Note:

The Summary of Safety and Clinical Performance of **Admira Fusion x-base** can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at www.voco.dental.

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ORMOCER® is a registered trademark of Fraunhofer Gesellschaft zur Förderung der angewandten Forschung e.V., München, Germany.