



Firma / company	Drucker Modell / printer model	V-Print							
		model beige	splint clear	SG clear	cast blue	Try-In beige	dentbase pink	tray blue	model fast
VOCO	SolFlex 650	✓	✓	✓	✓	✓	✓	✓	✓
VOCO	SolFlex 650 SMP	✓	✓	✓	✓	✓	✓	✓	✓
VOCO	SolFlex 350	✓	✓	✓	✓	✓	✓	✓	✓
VOCO	SolFlex 350 SMP	✓	✓	✓	✓	✓	✓	✓	✓
VOCO	SolFlex 170	✓	✓	✓	✓	✓	✓	✓	✓
VOCO	SolFlex 170 SMP	✓	✓	✓	✓	✓	✓	✓	✓

W2P	SolFlex 650	✓	✓	✓	✓	✓	✓	✓	✓
W2P	SolFlex 650 PLUS	✓	✓	✓	✓	✓	✓	✓	✓
W2P	SolFlex 350	✓	✓	✓	✓	✓	✓	✓	✓
W2P	SolFlex 350 PLUS	✓	✓	✓	✓	✓	✓	✓	✓
W2P	SolFlex 170	✓	✓	✓	✓	✓	✓	✓	✓
W2P	SolFlex 170 Plus	✓	✓	✓	✓	✓	✓	✓	✓

ASIGA	ASIGA MAX UV	✓	✓	✓	✓	✓	✓	✓	✓
-------	--------------	---	---	---	---	---	---	---	---

Rapid Shape	Rapid Shape D20 II	✓	✓	✓
Rapid Shape	Rapid Shape D30 II	✓	✓	✓
Rapid Shape	Rapid Shape D40 II	✓	✓	✓
Rapid Shape	Rapid Shape D20+	✓	✓	✓
Rapid Shape	Rapid Shape D30+	✓	✓	✓

Straumann	Straumann P20	✓	✓	✓
Straumann	Straumann P30	✓	✓	✓
Straumann	Straumann P40	✓	✓	✓
Straumann	Straumann P20+	✓	✓	✓
Straumann	Straumann P30+	✓	✓	✓

Dentalwings	Dentalwings D20 II	✓	✓	✓
Dentalwings	Dentalwings D30 II	✓	✓	✓
Dentalwings	Dentalwings D40 II	✓	✓	✓
Dentalwings	Dentalwings D20+	✓	✓	✓
Dentalwings	Dentalwings D30+	✓	✓	✓

Alle Druckpartner finden Sie unter /  
 All print partners can be found here  
[www.voco.dental/3dprintingpartners](http://www.voco.dental/3dprintingpartners)

✓ validiert / validation done 24.06.2020, mit kontinuierlichem Update/ updated continuously

Bitte stellen Sie sicher das Ihre CAM / Slicing-Software in einer aktuellen Version vorliegt. /  
 Please make sure that your CAM / Slicing software is available in a current version.