

VOCO Digital




Valid from 01.04. - 30.09.2019

Grandio® blocs · Grandio® disc

Nano-ceramic hybrid CAD / CAM material

- For the fabrication of definitive crowns, inlays, onlays, veneers and implant-supported crowns
- Highest filler content (86 % w/w)
- Resembles natural teeth perfectly
- Excellent physical values for flexural strength and abrasion
- No firing required
- Can be processed using any conventional milling unit
- Can be polished and repaired optimally
- Ideal for milling even in cases with thin margins
- Based on the nano-hybrid technology



NEW



Grandio blocs

5 × No. 12

5 × No. 14L

(REF see table below)

(REF see table below)

103,00 €

110,50 €

Set blocks 2 × No. 12 (A2 LT, A3 HT), 3 × No. 14L (A3 LT, A3 HT, A3.5 LT), Bifix QM QuickMix syringe 10 g universal, Futurabond U *SingleDose* 5 pcs., Ceramic Bond bottle 5 ml, Dimanto set, accessories

REF 6000 ~~266,00 €~~

save **25%**

199,50 €



Grandio disc

1 × 15 mm (ø 98.4 mm)
LT (low translucent)

1 × 15 mm (ø 98.4 mm)
HT (high translucent)

A1 REF 6050

A1 REF 6057

A2 REF 6051

A2 REF 6058

A3 REF 6052

A3 REF 6059

A3.5 REF 6053

A3.5 REF 6060

398,00 €

398,00 €

No. 12	A1	A2	A3	A3.5	B1	C2	BL	Mixed	No. 14L	A1	A2	A3	A3.5	B1	C2	BL	Mixed
LT	6003	6004	6005	6006	6007	6008	6009	6033	LT	6018	6019	6020	6021	6022	6023	6024	6034
HT	6012	6013	6014	6015	-	-	-		HT	6027	6028	6029	6030	-	-	-	

Mixed: A1 LT, B1 LT, C2 LT, BL LT, A1 HT

SolFlex 170 · 170 SMP

Additive manufacturing for digital dentistry

- DLP-Technology with a long-lasting UV-LED light source enables printing of clear transparent items such as splints
- Patented Flex-Vat technology – increasing printing speed while minimizing support structures, saving time and material
- High vat volume allows production over night
- Processing of widely-used STL format

SMP

- Maximized building speed with high process reliability due to SMP technology
- Increase in printing speed by up to 40 % with SMP technology



SolFlex 170
3D printer

REF 9100

12.900,00 €

SolFlex 170
SMP – 3D printer (incl. sensor technology)

REF 9101

14.900,00 €

V-Print

Material for 3D printing

- **V-Print model beige**
for the entire spectrum of models in dentistry
- **V-Print SG**
for the generative production of dental surgical guides
- **V-Print splint**
for splints used in treatments, auxiliary parts and functional parts for use in dental diagnostics, as well as bleaching trays (home bleaching)

**UV-LED
385 nm**



V-Print model
Bottle 1000 g beige

REF 6042

299,00 €

V-Print SG
Bottle 1000 g clear

REF 6043

432,00 €

V-Print splint
Bottle 1000 g clear

REF 6044

432,00 €

*compared with individual purchase

YOUR ORDER

Please contact your dental depot or your VOCO dental consultant to order the products presented here. Offers are only valid at participating retail partners. Special prices offered by your depot are always taken into consideration. All prices are RRP and do not include VAT.

THIS BROCHURE WAS HANDED OUT BY: