


Vococid®

ETCHANT FOR ACID-ETCH TECHNIQUE





Vococid®

ETCHANT FOR ACID-ETCH TECHNIQUE

Indications

Vococid is intended to etch only the part of the enamel or the dentine (in total-etch technique), which will subsequently be covered by a composite material:

- for composite fillings
- for adhesive luting of Maryland bridges, inlays, brackets and temporary splinting
- for sealing as caries prophylaxis
- for linings to create a permanent bond (according to manufacturer's instructions)

Advantages

- VOCO offers proven quality in suitable consistencies for application in every treatment situation
 - gel contains 35 % phosphoric acid whenever precision is essential for application: stable, non-drip, precisely accurate to apply with fine cannula
 - fluid contains 36.1 % phosphoric acid and can also be used for treating deep fissures and larger surfaces
- Stockage
 - the 60 ml refill syringe containing the proven Vococid etching gel offers a great price-performance ratio
 - the 2 ml syringes are easy and quick to fill and use
- Blue colouring
 - simplifies application and permits it to be checked visually
 - excess material can be identified and removed
- Achieves an optimal etching pattern



Presentation

REF 1063	Syringe 5 ml, accessories
REF 1066	Bottle 2 × 3 ml
REF 1224	Syringe 5×2 mI, accessories
REF 1256	Empty syringe 10×2 ml
REF 1257	Syringe 60 ml, empty syringe 2×2 ml, accessories
REF 2310	Application cannulae type 42, 30 pcs.
REF 2132	Application cannulae type 48, 50 pcs.

Cover picture: Dr. Marcelo Balsamo, São Paulo / Brazil

VOCO GmbH Anton-Flettner-Straße 1-3 27472 Cuxhaven Germany

Freecall: 00 800 44 444 555 Fax: +49 (0) 4721-719-140

info@voco.com www.voco.com

aße 1-3	Available from:
44 444 555 21-719-140	

